



MASSAGE HEAD ATTACHMENTS



ROUND

Use on both large and small muscle groups like the glutes, hamstrings and quads.



FLAT

Designed for more dense muscle groups such as the shoulder blades and IT bands.



BULLET

Pinpoint smaller muscle groups including the hands and feet.



SPINAL

Suitable for trigger points and the lower back.


CORENINE MICRO MASSAGER



INTERCHANGEABLE MASSAGE HEAD ATTACHMENTS



AIR VENTS



LED CHARGE INDICATOR

- Charge is greater than 40%
- Charge is less than 40%
- Charge is less than 10%

CHARGING PORT

POWER/SPEED (ON/OFF)

CHARGE INDICATOR

INSTRUCTIONS

CHARGE the CORENINE MICRO MASSAGER for one hour before using the first time or after a prolonged period of no use.

POWER/SPEED: Press and hold the power button to turn on and off. When the massager is on, press the power button to change the speed.

PLEASE NOTE: The massager will automatically turn off in 15 minutes.

RECOMMENDED USAGE: One minute per body region.

MEDICAL WARNINGS

Despite its benefits, massage isn't meant to replace regular medical care.

If you have a medical condition, do not use this product without consulting your physician. This product is not intended for the self-treatment of medical conditions that should be managed by a qualified health care provider.

- **NEVER** use without consulting your physician if you are pregnant, if you have been recently injured or have had surgery, or if you have a malignant tumor.
- If you feel any abnormal discomfort while using this product, stop using and consult your physician.
- If bruising occurs, do not use on bruised areas until completely healed.
- Use only on dry skin surfaces or over clothing.
- Use caution when using on soft tissue that connects, supports, or surrounds the organs of the body.
- **NEVER** use on the head or on sensitive or swollen areas of the body.
- **NEVER** use on bony areas such as the tibia, the top of the feet, or the back of the hands.

SAFETY PRECAUTIONS



Connect the charging cord to a properly grounded power outlet. Grounding provides a path of least resistance for electrical current to reduce the risk of electric shock if the device malfunctions or breaks down

- **NEVER** charge this device with a damaged power cord. If the charging cord is damaged, contact your authorized retailer for a replacement.
- To avoid electric shock, **NEVER** use wet hands to connect or disconnect the charging cord from a power outlet.
- Keep fingers, hair, and other body parts away from the percussion axle.
- The massager should **NEVER** be used anywhere near jewelry worn, belts, or other loose clothing items.
- **NEVER** place objects or liquids into the ventilation holes of this device.
- If liquid is accidentally spilled onto the device, immediately turn it off the power and unplug the power cord.
- **NEVER** place the charging cord or the massager into water or near damp or moist areas such as swimming pools, sinks, or bathrooms.
- Keep the charging cord and massager away from heated surfaces, fireplaces, direct sunlight, or excessively hot, humid, or dusty environments.
- To avoid electric shock, **NEVER** operate under a blanket or pillow.
- Do not use any attachments that are not provided.
- Do not insert, rotate or remove the massage head attachments when the power is ON.
- Do not change or modify this device in any way.

This product is only recommended for healthy adults; it is not suitable for children.