# Neck & Shoulder Shiatsu Massager

## Warnings

- **Read all instructions thoroughly and retain for future reference.**
- This appliance can be used by children aged 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and maintenance shall not be made by children without supervision.
- This appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance using adapter IVP1200-2000 supplied.
- If the power lead or AC/DC adapter are damaged, it must be replaced by the manufacturer or its service agent or similarly qualified person in order to avoid a hazard. Contact Customer Care to arrange a replacement.
- Do NOT use this appliance if it is not working properly, or if it has been dropped or damaged.
- Do not leave the device unattended. Turn off before you leave the room.
- Do NOT immerse in water.

## Know your unit

- |                  |             |
|------------------|-------------|
| 1 Massage pillow | 5 Heat      |
| 2 Controls       | 6 Speed     |
| 3 Arm loops      | 7 Direction |
| 4 Power adapter  | 8 On/Off    |

## How to use

- Connect the power adapter and plug into the mains supply.
- Place on the area you wish to massage and press the On/Off button.  
**Note:** The massage will last for 15 minutes before automatically turning off.
- Pressing the On/Off button again will turn the unit off.

### DIRECTION:

- Pressing the Direction button will change the direction of rotating massage heads.

### SPEED:

- To select between the 3 speed levels (low, medium & high) press the SPEED button.

- **Note:** To increase intensity apply pressure through the Arm Loops.

### HEAT:

- Pressing the Heat button will turn the heat On/Off.

## After Use

- You may experience a little muscle stiffness after first use. This will ease and lessen with regular use.

## Cleaning & Maintenance

- Unplug the device and clean the unit using a soft damp cloth, allow to air dry

## Important Information

There are no user serviceable parts inside the appliance.

## Disposal

 The Waste of Electrical and Electronic Equipment (WEEE) Directive (2012/19/EU) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill. This product is classed as Electrical or Electronic Equipment so please ensure that at the end of the life of this product it is disposed of in the correct manner in accordance with local authority requirements. It must NOT be disposed of with household waste.

## Specifications

Model: BT-MANS  
Input: 12V  2000mAh

 Read these instructions carefully before use and then retain for future use.

 This product is for indoor use only.

 This unit complies with current EU product legislation.

 Adapter IVP1200-2000  
100-240V 50/60Hz  
0.5A(UK) 0.7A(EU)

 The appliance must only be supplied at the extra low safety voltage indicated on the rating plate.

 Alternating current

 Direct current

 Centre pin positive

**UK  
CA**

## Guarantee

To register your two year guarantee online, visit: [www.bodi-tek.co.uk/warranty](http://www.bodi-tek.co.uk/warranty)

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©The Dezac Group Ltd 2021, PO Box 17, Cheltenham Spa, England GL53 7ET.
© The Dezac Group Ltd 2021, PO Box 17, Cheltenham Spa, Inglaterra GL53 7ET.