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-If the temperature is above 37.8°C, a double «Beep» signal will sound 2 to 10 times to indicate that your baby may have fever.  
-Remove the pacifier thermometer from your baby's mouth and read the temperature displayed on the LCD screen.  
-Stop the thermometer by pressing the On/Off button again.  
-To protect the battery, we recommend that you switch off the thermometer once you have taken the temperature measurement. However, the thermometer will turn itself off automatically if the temperature measurements are taken within 10 minutes.  
-The «Filoo» pacifier thermometer is designed to work for 200 hours. The «» symbol will flash if the battery's low. When this happens, simply replace your pacifier with a new «Filoo» pacifier thermometer, ensuring you comply with the regulations concerning waste disposal.  
-If for any reason you need to interrupt the measurement, start the procedure described above.  
NOTES:  
-Oral temperature measurements are usually 0.5 to 1°C lower than the actual body temperature.  
-If the body temperature is also 42.0°C, the LCD screen will display «4H». If it is below 32.0°C, the LCD screen will display «4L».

MAINTENANCE AND CLEANING  
-You can sterilize the thermometer by immersing the silicon latex teeth ONLY in warm water.  
Please note:  
-You must NEVER completely immerse the pacifier thermometer in boiling water, as you risk damaging the electronic components.  
-Do not use abrasive detergents or chemical solvents.  
ADVICE REGARDING THE DISPOSAL OF WASTE  
-All the packaging is composed of materials that pose no hazard for the environment and can be disposed of in your local sorting centre to be used as secondary raw material. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.  
-When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the battery and discard it in a collection point that it can be recycled.  
Please note:  
-Under no circumstances should spent batteries be disposed of with household waste!!!

LIMITED WARRANTY  
LANAFORM® guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, unless in the circumstances described below.  
The LANAFORM® guarantees does not cover damage caused as a result of normal wear or tear. In addition, the guarantee covering this LANAFORM® product does not cover damage caused by abuse or inappropriate or incorrect use, accidents, the use of unsuitable accessories, changes made to the product or any other circumstances, of whatever sort, that is outside LANAFORM®'s control.  
LANAFORM® may not be held liable for any type of circumstantial, indirect or specific damage.  
-All implicit warranties relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.  
On receipt, LANAFORM® will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM®'s authorised service centres. Any attempt to maintain this product by a person other than LANAFORM®'s Service Centre will render this guarantee void.

EC 0029 DECLARATION OF CONFORMITY  
Medical device  
Category Ia.  
We, the company LANAFORM,  
Declare that the medical device «Filoo» LA090110 complies with Council Directive 93/42/EEC concerning medical devices and that it meets the requirements of the CE mark.  
APRAGAZ, 156 chaussée de Wivroville, 1120 Brussels (Belgium), notified body no. 0029, has checked this product (according to Annex IV).

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Symbol	Definition	Identification	Information provided by the manufacturer
	Manufacturer's details		LANAFORM SA B-4141 SPRIMONT
	Batch number		001/54131012
	BF device		
	Please note: see instructions for use.		Please note: consult the enclosed documents. Please read these instructions before using the device. Please also follow these instructions carefully when using the device.
	Notified body		001/54131012
	Electrical and electronic waste (WEEE)		

**PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR «PACIFIER THERMOMETER», ESPECIALLY THIS BASIC SAFETY INFORMATION:**  
- Do not place the thermometer in direct sunlight, or expose it to high temperatures, humidity or dirt.  
- Do not use this device in a room in which aerosols (sprays) are used or in a room in which oxygen is administered.  
- Do not immerse the device in water.  
- Do not drop the device and do not use it if it is damaged.  
- This device must be used at an ambient temperature (between 10°C and 40°C).  
- Do not expose the thermometer to extreme shocks.  
- Do not expose this thermometer to electronics of 30°C or below -20°C.  
- Do not use the device if the relative humidity is above 90% or below 20%.  
- Do not use the device in the vicinity of large electromagnetic fields, such as cellphones or portable telephones.  
- It is recommended to place the thermometer in the room for around 10 minutes before use.  
- Use of this thermometer is not a substitute for consulting your doctor.  
- If there is a problem with your thermometer, please contact your reseller.  
- Do not try to repair the thermometer yourself.  
- This unit is not intended for use by people (including children) with reduced physical, sensory or mental capacities, or persons with lack of experience or knowledge, unless a person responsible for their safety can supervise or instruct them beforehand on how to use the device. Children should be supervised to make sure they do not play with the device.  
- Only use the thermometer for its intended purpose, as described in the manual.  
- Before each use, examine the pacifier, checking for any cracks, tears, dirt or other signs of use. Do not use the pacifier thermometer if there is the slightest sign of damage.  
- Do not try to pull the pacifier or twist it, or take it out of the device by force, as you might damage the temperature sensor and electronic components, and consequently the sensitivity of the thermometer.  
- Do not try to pull the pacifier or twist it, or take it out of the device by force, as you might damage the temperature sensor and electronic components, and consequently the sensitivity of the thermometer.  
- Do not let children play with the device. It is not a toy. The thermometer is a sensitive electronic device which could be damaged if it is not used correctly in the presence of an adult.  
- Do not attach a ribbon or string to the pacifier thermometer as your child could be strangled by it.  
- Do not use the pacifier thermometer for temperature measurements taken orally.  
- Do not use the pacifier thermometer like a normal thermometer.  
- The battery is not replaceable. Do not try to remove it.

FEATURES  
- LCD screen to display the temperature.  
- Last measurement in memory.  
- Sound alert to announce that the maximum temperature has been reached.  
- Double sound alert when the temperature exceeds 37.8°C.  
- Automatic function for 10 minutes.  
- Protective cap to transport the product hygienically

DESCRIPTION OF THE THERMOMETER:  
1- On/Off button  
2- LCD screen  
3- Pacifier  
4- Sensor  
5- Protective cap

INSTRUCTIONS FOR USE  
The following table is a guide showing average human body temperatures according to age.  
**Normal temperatures according to age.**

Normal temperatures according to age	Temp. °C
0 to 2 years	36.4 to 38.0
3 to 10 years	36.1 to 37.8
11 to 65 years	35.9 to 37.6
>65 years	35.8 to 37.5

NOTES:  
-The temperature of the human body varies throughout the day. It can also be influenced by a number of external factors: age, gender, type and thickness of skin, etc.

TECHNICAL SPECIFICATIONS

Normal conditions of use	Temp. °C
Display resolution	0.1°C
Operating temperature	10 to 40°C
Storage temperature	-20°C to 60°C
Humidity level	20% <H< 90%
Weight	+/- 14 grams

**Measurement range**  
32°C to 42.9°C  
Precision  
+0.1°C 35-38°C  
+0.2°C 38-42.0°C  
Automatic stop  
10 mins

UTILISATION  
-Appuyez sur le bouton On/Off pour allumer le thermomètre. Les symboles «188.8», «°C», «M» et «» will appear for approximately 2 seconds on the LCD screen. This means that the unit is working correctly.  
-The last temperature measurement will appear once on the screen.  
-Place the pacifier in your baby's mouth and let them suck on the pacifier for about 1 minutes.  
-Make sure the pacifier is correctly placed in your baby's mouth and that the sensor is touching their tongue.  
-The «L» symbol will flash until the maximum temperature has been reached.  
-A sound alert will announce that the temperature has been reached.  
-If the temperature is normal (below 37.8°C), the «Beep» signal will be taken 10 times.

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