



iHealth HS6 CORE Body Analysis

1 Device Installation

1. Remove the battery protection located in the batteries compartment. SET is displayed on the scale. Note - If SET is not displayed on the scale, open batteries compartment, remove and reinstall one battery then press the SET button.
2. Launch iHealth MyVitals > MENU > Set up my devices > Select New Device > Scales > Core (HS6)
3. Tap Scan QR code then scan the QR code printed on the back of the scale to register a new user in the scale.
4. Enter your current weight, then tap Done.
5. Configure the Wi-Fi connection
6. Enter your Wi-Fi password

2 Connecting to SNUG

1. Download Snug from the Google Play or Apple App Store and create your Snug account
2. Access "My Devices" through the app menu in the top left-hand corner
3. Select "Add a data source", choose "Human API", then press "Okay"
4. Choose "iHealth", add your iHealth account email and password (created in Step 1), then "Sign In"
5. Once Human API has processed, you can close this window manually to return to Snug
6. "iHealth" will appear in your list of Data Sources
7. Allow up to 30 minutes for data to flow into Snug from iHealth
8. For any connectivity issues contact support@snughealth.com.au

