

MODEL # 4832

**TAYLOR**<sup>®</sup>

## Mechanical Dial Bathroom Scale

Capacity: 330lbs

### INSTRUCTIONS FOR USE & WARRANTY INFORMATION

#### CUSTOMER SERVICE INFORMATION PLEASE READ!

For questions regarding this product, please contact us at:

**TAYLOR** division of Lifetime Brands, Inc.  
2220 Entrada Del Sol, Suite A | Las Cruces, New Mexico 88001 USA  
**1-866-843-3905** | [www.taylorusa.com](http://www.taylorusa.com)  
Email: [taylorusa@lifetimebrands.com](mailto:taylorusa@lifetimebrands.com)

We value your thoughts on our product and services. Please post an on-line review to help us continue to meet your needs and expectations.

**Please contact us before returning to retailer.**

#### OPERATING INSTRUCTIONS

- Place your scale on a hard, flat surface, do not use on uneven tile, flooring or carpeting. Try to keep your scale in the same location or position as much as possible.
- Before using scale, allow it to reach room temperature; then put scale through several weighing cycles consisting of:
  - *Step on scale and allow dial or pointer to come to rest.*
  - *Step off scale and allow dial or pointer to come to rest.*
- If scale is relocated, readjustment to zero may be required. Repeat above procedure. When weighing, place feet on platform so weight is evenly distributed.
- After weighing cycles are completed, adjust dial or pointer to zero, if necessary (adjustment knob is located under edge of platform or, on some models, under edge of dial), and repeat weighing cycles until zero is stabilized.

#### CARE & MAINTENANCE

An occasional wiping with a soft cloth and mild soap will keep your scale attractive. Do not use abrasive cleaners. Your scale does not require oiling or lubrication of any kind.

# Contact Information

## Distributed By:

### **Taylor, a division of Lifetime Brands Inc.**

2220 Entrada Del Sol, Suite A  
Las Cruces, New Mexico 88001 USA

**1-866-843-3905**

*www.taylorusa.com*

*Email: taylorusa@lifetimebrands.com*

**MADE IN CHINA.**

### ***Five (5) Year Limited Warranty***

This scale is warranted against defects in materials of workmanship for five (5) years for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. Should this scale require review (or replacement at our option) while under warranty, please pack the item in the original packaging and return it prepaid, along with store receipt showing date of purchase and a note explaining reason for return to:

#### **Taylor division of Lifetime Brands Inc.**

2220 Entrada Del Sol, Suite A  
Las Cruces, New Mexico 88001 USA

**1-866-843-3905**

*www.taylorusa.com*

*Email: taylorusa@lifetimebrands.com*

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information please contact us through:

*www.taylorusa.com*

If assistance is required, do not return to retailer. For information call 1 (866) 843-3905 from 8:00 am to 5:00 pm, Mountain Standard Time, Monday through Friday. To assist us in serving you, please have the model number and date of purchase available.

**Not legal for trade.**

**Made to our exact specifications in China.**

© 2020 Taylor, a division of Lifetime Brands Inc. All rights reserved.

Ref 4832

12.20 WC

