



LUMIA



## LUMIA RF GEL

THE LUMIA RF GEL IS UNIQUELY FORMULATED TO PREPARE YOUR SKIN FOR THE FULL BENEFITS OF THE EXPERTO AND EXPERTO + TREATMENT DEVICES. IT ENSURES THE DELIVERY THE RF ENERGY INTO THE CORRECT LAYERS OF THE DERMIS AND PROTECTS THE SKIN FROM DAMAGE.

APPLYING A THIN LAYER OF THE GEL TO CLEAN AND DRY SKIN WILL SUPERCHARGE YOUR DEVICE AND OPTIMIZE THE RESULTS OF YOUR TREATMENT

# TREATMENT TIPS

## STEP 1

Cleanse the area. The area should be cleared of all creams and ointments.

---

## STEP 2

Lumia RF Gel must be applied before treatment

---

## STEP 3

Spread a thin layer over the treatment area.

---

## STEP 4

RF mode: slow linear motions over the treatment area.

---

## STEP 5

Gently clean your skin.

# LUMIA

— TECHNOLOGY MEETS BEAUTY —

# RÉGEL



[WWW.THELUMIA.COM](http://WWW.THELUMIA.COM)