

LEG WRAP: QUICK START GUIDE


⚠ ATTENTION: CONSULT ACCOMPANYING DOCUMENTS BEFORE USE.

1

Place Gel Pack in freezer for 12 hours prior to use.



2



Attach cold Gel Pack inside of wrap using velcro.

3



Place supplied cloth over leg. Squid also works over clothing.

4

Wrap Squid around the leg.*



⚠ ATTENTION: CONSULT ACCOMPANYING DOCUMENTS BEFORE USE.


5

Fasten the wrap around the leg snugly using the velcro panels.



6

Attach connector tube to pump and wrap.



7



Begin 15 Minute Treatment.

*

TREATMENT TIP:

The Squid Leg wrap can be positioned in a variety of ways to better treat affected areas:

- QUADRICEP
- HAMSTRING
- KNEE
- CALF
- SHIN