

LUCE

DERMAL DEFINER HOW TO USE

With minimal time investment, you can achieve noticeable improvements in the overall tone, texture, and vitality of your skin, making the LUCE Dermal Definer a valuable addition to your daily skincare regimen.

SUPPLIES NEEDED:

- + LUCE Dermal Definer Multipurpose tool
- + Gentle Face wash
- + Moisturizer or serum for roller

HOW TO USE THE BRUSH PART

Prepare Your Face:

Wet your face with lukewarm water to open pores.

Apply Cleanser and Use the Brush:

Apply a suitable cleanser to your brush or directly to your face. Turn on the brush and gently use circular motions to cleanse your entire face.

Customize Your Approach:

Adjust the pressure and speed as needed for comfort.

Focus on problem areas that require extra attention.

Rinse and Clean: Rinse your face thoroughly with lukewarm water. After use, clean the brush by rinsing and drying the bristles.

Recommended use: Twice a day. Mornings and before bed.

HOW TO USE THE ROLLER PART

1 PREP

Start with a clean face: Ensure your face is free from makeup, dirt, and impurities.

Apply a serum or moisturizer: This helps the roller glide smoothly and enhances product absorption.

2 ROLLING TECHNIQUE

Begin at the neck: Start rolling from the base of your neck, moving upward toward the jawline. Use gentle, upward strokes. Move to the jawline and cheeks. Work on the forehead: Roll horizontally across your forehead, moving from the center toward the hairline.

Under-eye area: Use the smaller end of the roller to roll gently under the eyes, moving from the inner corner outward.

3 CLEAN AND STORE

After use, rinse the roller with warm water and mild soap. Pat the roller dry with a clean towel or let it air dry.

Combining this multi purpose tool in your skincare routine offers deep cleansing, gentle exfoliation, and enhanced product absorption, promoting a healthier, more radiant complexion. The gold roller massaging effect can help relax facial muscles, potentially reducing the appearance of fine lines and wrinkles

