

## Using TePe Compact Tuft™

TePe Compact  $Tuft^{\mathsf{TM}}$  is a dome-shaped, single tufted brush with tightly packed short bristles that gives a firm feel. This versatile brush works well for implants, attachments for overdentures, fixed braces, along the gum line, emerging or partially erupted teeth, misaligned teeth, and removing surface stain. It is important to use and care for your TePe Compact  $Tuft^{\mathsf{TM}}$  properly—please read below for instructions and helpful tips.



## Simple Steps to Using Your TePe Compact Tuft™:

Due to its versatility, instructions for using TePe Compact Tuft™ are unique to your situation.

- 1. When using along the gumline, hold the brush tip at a 45-degree angle towards the gum line, moving the tip gently back and forth.
- 2. When using on the top of the tooth, hold the brush tip at a 90-degree angle towards the biting surface of the tooth. Move the brush tip gently with small, circular movements. This will clean the grooves and help prevent cavities. This works for emerging teeth or partially erupted teeth.
- 3. When cleaning around orthodontic appliances, implants, or attachments for overdentures use minimal pressure, moving in a back and forth motion.

Scan the QR code to watch a video about
TePe Compact Tuft™



## Tips for using your TePe Compact Tuft™

- It is best to complete this in a systematic approach to ensure you don't miss any areas.
- Replace brush once bristles start to fray or become worn.









