

## Using TePe® Implant Orthodontic Brush

TePe® Implant Orthodontic Brush has two rows of soft bristles that create an extra narrow brush head to facilitate cleaning around orthodontic appliances, as well as the narrow area where implants meet the gum line. It is important to use and care for your TePe® Implant Orthodontic Brush properly—please read below for instructions and helpful tips.



## 4 Simple Steps to Using Your TePe® Implant Orthodontic Brush:

- 1. Use daily to disrupt bacteria.
- 2. When using to clean the narrow area where the implants meet the gum line, hold at a 45-degree angle, moving the bristles in a small back and forth motion along the gums.
- 3. When using this brush for orthodontics, place the bristles under the brackets and wire and move back and forth. Repeat this process, placing the bristles above the brackets and wire.
- 4. After use, rinse and air dry in an upright position.

Scan the QR code to watch a video about TePe® Implant Orthodontic Brush



## Tips for using your TePe® Implant Orthodontic Brush

- (i) It is best to complete this in a systematic approach to ensure you don't miss any areas.
- If using along a tooth or implant, focus on the gumline, because this is where harmful bacteria resides.
- Angle the neck in hot water to adjust to your liking.
- Replace brush every 3-4 months or once bristles start to fray or become worn.







