



# Using TePe Mini Flosser™

TePe Mini Flosser™ is a preloaded floss holder with a unique bite plane to gently insert the floss between teeth. The floss is both strong and slim, fitting into even the narrowest of spaces. It is important to use and care for your TePe Mini Flosser™ properly—please read below for instructions and helpful tips.



## 4 Simple Steps to Using Your TePe Mini Flosser™:

1. There are two ways to insert the device between the teeth
  - Using the unique bite plan, you can simply bite down on the device controlling the speed and pressure with your jaw movement.
  - The other option is to move the string of floss in a horizontal back and forth motion until you get through the contact area.
2. Once your floss is inserted, move the device up and down against each tooth surface, holding it tightly against the tooth.
3. Remove the TePe Mini Flosser™ device from between the teeth.
4. TePe Mini Flosser™ is intended for single use, to be used throughout the entire mouth; dispose after use

Scan the QR code to watch a video about TePe Mini Flosser™



## Tips for using your TePe Mini Flosser™

- ① It is best to complete this in a systematic approach to ensure you don't miss any areas.
- ① Looking in a mirror can help confirm you are using the correct technique.

