



Using TePe Universal Care™

TePe Universal Care™ has a slim brush head design with a unique angled neck to enable easy access for hard-to-reach areas, specifically the inside surfaces of teeth or implants. It is important to use and care for your TePe Universal Care™ properly—please read below for instructions and helpful tips.



5 Simple Steps to Using Your TePe Universal Care™:

1. Use daily to disrupt bacteria in those hard-to-reach areas.
2. Evaluate the area in your mouth where you will be using the brush.
3. Using small back and forth motions is preferred around single implants, tongue side of lower molars to prevent gag reflex, the back of the last molar, tilted teeth, and around cemented retainers.
4. A sweeping motion is ideal for under bridges or implant-retained dentures.
5. After use, rinse and air dry in an upright position.

Scan the QR code to watch a video about TePe Universal Care™



Tips for using your TePe Universal Care™

- It is best to complete this in a systematic approach to ensure you don't miss any areas.
- If using along a tooth or implant, focus on the gumline, because this is where harmful bacteria resides.
- Angle the neck in hot water to adjust to your liking.
- Replace brush every 3-4 months or once bristles start to fray or become worn.

