

5

MAINTENANCE

1. Wipe off with a dry cloth after each use. Do not wash off, Mask is not water resistant.
2. Turn off before wiping down this mask.
3. Always unplug after use and do not disassemble the product.

STORAGE

1. Store this product in a cool and dry place avoiding humidity and direct sun light or store it in the box.
2. Do not submerge the product in water or liquid
3. Keep out of reach of children.

TROUBLESHOOTING


Trouble	Action for trouble
The light is not very bright	Charge the mask
Charging issues	Change the plug/ usb
The light flashes	Charge the mask
Others	Contact Customer care

6

SPECIFICATIONS


Product Name	7 Color LED Beauty Mask
Product Size	195*165*125mm
Product Weight	76g
Materials	ABS
Output Voltage	5V / 1A
Accessories	Mask, Manual, Cable Line
Power Supply	Rechargeable Battery
Charging Time	2 Hours

PRODUCT APPEARANCE



Verdermis

7 Color LED Beauty Mask



1

User Manual

Please Read This Manual Carefully Before Using

4

Verdermis

7 Color LED Beauty Mask



TO MAXIMIZE RESULTS SCAN HERE
 For more information about the device please visit
WWW.VERDERMIS.COM
 Dist. by Y&E Ent., LLC
 Atlanta, GA 30338 USA
 71-75 Shelton Street, London, UK
 Hauptstrasse 22 Neugersdorf 02727 Germany

PRODUCT BENEFITS

- ※ The mask uses different types of light sources including, RED, BLUE, GREEN, YELLOW, PURPLE, CLEAR BLUE, and BRIGHT colors.
- ※ Light therapy may assist in improving overall skin condition and appearance.

FUNCTIONS

Red Light (650-730nm)	May contribute to overall skin appearance.
Blue Light (525-550nm)	Associated with skin nourishment.
Green Light (430-450nm)	Calming the skin.
Yellow Light (565-570nm)	May assist with toning of the skin.
Purple Light (380-440nm)	Balancing the skin.
Clear Blue (160-195nm)	Balancing, calming, and soothing the skin.
Bright light	May contribute to overall skin conditions.

3


DIRECTIONS

1. Cleanse your face.
 2. Apply your favorite skincare products (preferably Verdermis).
 3. Press ON/OFF button for 3 seconds to Turn On.
 4. Click ON/OFF button to select desired LED light.
 5. Place Mask on face for 10 minutes.
 6. To Turn Off press the ON/OFF button for 3 seconds.
- ※ Use as needed (auto timer is set for 10 minutes per use).
 - ※ Do NOT use while charging.

NOTES

1. This product is only recommended for use by an adult.
2. Mask is designed to distribute lights evenly on the face with no side effects.
3. The linear area is only the light guide and not the real source of light.
4. The circuit inside this product should not be pressed or folded and it should be placed flat to protect the product after daily use.
5. In the case of sensitivity stop using this mask.
6. Using the Mask after applying skincare (preferably Verdermis) shall improve results.
7. Use as needed. Each application for up to 10 minutes.
8. In the rare case of any discomfort stop using this mask and contact customer care.
9. While charging, red light next to charging port indicates that the mask is being charged. When the red light is turned off, mask is fully charged.

折叠方法



反面456

正面123