

IN THE BOX

Accessories included in your package:

- A Control Device
- B USB Charging Cable
- C Waist Belt
- E Operations Manual
- D Extension Belt
- F Quick Start Guide



USING THE CONTROL DEVICE

- 1 BUZZ / SOUND**
+ Press to turn device sounds on/off.

- 2 TIME ADJUSTMENT**
+ Device defaults to 20 minutes.
+ Press the Time button to increase time by 10 minute increments. (10 - 60 minutes)

- 3 TIME INDICATOR LIGHTS**
+ The light closest to the **Time button** = 10 minutes.
+ The light closest to the **Buzz button** = 60 minutes.
- NOTE:** The device will also beep once for 10 minutes, twice for 20 minutes, and so on.



- 4 EMS TONE CONTROLS**
+ Press (+) to increase intensity. (up to level 20)
+ Press (-) to decrease intensity.

- 5 LIPO LIGHT TRIM CONTROLS**
+ Press TRIM to turn Red LED Lipo-Lights on / off.



- 6 USB PORT**
+ To charge battery, plug provided USB charging cable into a USB port on your computer or certified wall charger.

- 7 ON/OFF BUTTON**
+ Press to turn ON (will beep once).
+ Press and hold to turn OFF (will beep twice).

- 8 BATTERY LIFE LIGHTS**
+ If fully charged, the 8 battery indicator lights will be illuminated. If only 1 light is illuminated, the device needs charging.



QUICK START GUIDE

STEP 1

CHARGE THE BATTERY

- + Push the ON/OFF button located on bottom of control device.
- + 8 lights = charged. If device needs charging only 1 light will illuminate.
- + Use included USB cord to charge device. Lights will flash while charging.
- + When all lights are illuminated, device is fully charged



STEP 2

WRAP BELT AROUND WAIST

- + Use a damp rag or cloth to *slightly* wet the silicone portion of the belt (non-fabric area). **DO NOT** immerse in, or soak belt in water.
- + Wrap belt around your waist in direct contact with skin and secure. Do not use over clothing or other materials.
- + Use the included extension belt if required.



STEP 3

ATTACH CONTROL DEVICE TO WAIST BELT

- + Match the 4 magnetic buttons on the back of the control device to the 4 buttons on the belt.
- + Match the 2 gold buttons to the 2 gold buttons on the belt
- + Match the 2 silver buttons to the 2 silver buttons on the belt.



STEP 4

TURN ON DEVICE

- + Press ON/OFF button
- + Press TRIM button to turn on the Lipo-Lights



STEP 5

SELECT TREATMENT TIME

- + Time defaults to 20 minutes
- + Press TIME button located on top of control device
- + Increase by 10 minutes up to 60 minutes



Refer to Using the Control Device #3 "Time Indicator Lights"

STEP 6

INCREASE / DECREASE INTENSITY

- + Press (+) to increase intensity up to level 20 (will beep twice when increased to highest and lowest level)
- + Press (-) to decrease intensity
- + There is a 2 minute "Warm Up" program with lower pulse frequency.



BE CAREFUL OF INTENSITY LEVEL AS PULSE FREQUENCY WILL SIGNIFICANTLY INCREASE AFTER 2 MINUTES

STEP 7

REMOVE AND STORE

- + Control Device will automatically shut off when time is up (will beep twice when turning off)
- + Press and Hold the ON/OFF Button to end treatment at any time.
- + Remove Belt wrap, remove Control Device from belt and store separately.



PLEASE REFER TO THE OPERATIONS MANUAL FOR COMPLETE INSTRUCTIONS BEFORE USE.