

DO NOT use the Omnilux Contour™ LED system to treat any other conditions apart from those listed in the indications for use. The Omnilux Contour™ LED system has not been tested for any other conditions than those listed and the risk is unknown.

DO NOT use the Omnilux Contour™ LED system if you are pregnant, planning on becoming pregnant or breastfeeding. The Omnilux Contour™ LED system has NOT been tested on pregnant or breast-feeding women and therefore the risk to the fetus, neonate or pregnant women is unknown.

DO NOT use the Omnilux Contour™ LED system if you suffer from lupus erythematosus, photosensitive eczema or albinism. If you use the Omnilux Contour™ LED system to treat lupus erythematosus, photosensitive eczema or albinism you may cause a severe skin reaction.

DO NOT use the Omnilux Contour™ LED system if you suffer from any photosensitive disorder (sensitization to light). If you use the Omnilux Contour™ LED system and you suffer from a photosensitive disorder you may cause a severe skin reaction.

DO NOT use the Omnilux Contour™ LED system if you are taking any medication that can cause photosensitivity. If you use the Omnilux Contour™ LED system and you are taking any medication that can causes photosensitivity you may cause a severe skin reaction.

DO NOT use the Omnilux Contour™ LED system if you suffer from light induced headaches.

DO NOT use the Omnilux Contour™ LED system if you suffer from any genetic conditions of the eye.

Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs, and diuretics. If you are unsure about any medication you may be taking consult your healthcare provider.