Arch Refresh

FOOT MASSAGER

Before Use

Thank you for choosing the Arch Refresh Foot Massager. Please read the user manual carefully before use in order to understand the correct operations. Please keep the user manual for reference. *Note: product updates and design specifications are subject to change without notice.*

Technical Specifications

Product Name: Arch Refresh

SKU Number: FM160

Output: 24 Vdc, 1.5 A, 36 W

Location of Controls

Features

- Ergonomic design to let you enjoy a comfortable and professional massage.
- High quality material and durable design. Improve circulation and reduce stress with this acupressure massage designed to improve overall health from the comfort of your own home or office.
- 4 in 1 Massage Design Rolling massage, acupressure, air pressure massage and thermotherapy in one massager.
- Safe Dual Protection Design. Protection from overheating and electric current overloading

Plus

- Full foot massage: toes, insteps, soles and heels.
- Choice of three levels of adjustable intensity.
- 15-minute automatic timing.
- Washable inner massage fabric cover

To Operate

1. **ON/OFF:** Press the On/Off button, then the massager starts working. The massager begins the default mode: Vibration 2, Mode 1, Air Pressure 1, heat function on, and the indicator light of corresponding function is on.
2. **MODE:** The Mode function defaults as MODE 1. Press the Mode button to change to MODE 2. Press the Mode button again for MODE 3. Different modes correspond to different indicators. Each mode adopts different rhythm, inflating and deflating movements to simulate the real effect of a manual massage. When the button is pressed again, it returns to MODE 1 and cycles.
3. **VIBRATE:** The massager has 4 vibration modes: Low (indicators: about 2 flashes in 1 second) the slowest vibratory rhythm; Medium (indicators: about 1 second flash) medium vibratory rhythm, High (indicators: long light) the fastest vibration rhythm; and Pure Air Bag mode: (the indicator light is not on) the vibration and roller stop, and only the air bag is used for aeration – air release massage cycle.
4. **AIR PRESSURE:** Circularly control the pressure intensity from level 1 to level 2 to level 3 (the massage strength increases in turn), and the corresponding indicator light is on.
5. **HEAT:** The thermotherapy function is automatically on and the light will be red. Press the HEAT button to turn it off.
6. The massager will automatically turn off after 15 minutes.

Warning

1. Only use the adapter with voltage DC24V.
2. Unplug the power cord from the outlet after use.
3. Do not open the massager or change any parts by yourself.
4. Do not use the massager more than 30 consecutive minutes.
5. Do not use and/or store the massager in damp or high temperature environment.
6. If the power cord gets damaged, it must be replaced by original manufacturer, or authorized distributors.
7. Do not stand on or place objects on top of the massager.
8. If you feel light headed, dizzy or nauseated, please stop using the product immediately.
9. If you smell a strong electrical odor unplug the power cord immediately.
10. The massager has overheat protection.
11. Please consult your physician if you are pregnant, have skin diseases, heart diseases, skin redness, brittle bones, hypertension, tumors, pacemaker or other similar implanted devices, or osteoporosis before using the product.
12. Keep the massager away from flammable, explosive or chemical solutions.
13. Please use a soft cotton cloth and neutral soapy water to clean the massager.

Trouble Shooting

PROBLEM	SOLUTION
The product does not work.	<ol style="list-style-type: none">1. Check the connection between the power cord and the outlet.2. Check the Power button to see if the product has timed out.
The massage function does not work.	<ol style="list-style-type: none">1. The motor is over heated. Please turn off the product and allow it to cool off for 30 minutes.2. Check the Power button to see if the product has timed out.
The massage intensity is too strong.	<ol style="list-style-type: none">1. Please make sure the product is stationary.2. It is best to adapt to the massage intensity level. Start from the lower intensity level and then move up to the higher levels.

The product has overheating protection. It will automatically stop working when it reaches too high of a temperature. Please unplug the power cord and let it cool down for at least 30 minutes.

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