


geniani

PLEASE READ AND SAVE
THESE INSTRUCTIONS

INSTRUCTION MANUAL

ROCK CORDLESS MASSAGER




**Inquiries
Help
Issues**

1 (404) 225-3955

Monday to Friday 9 AM to 5 PM EST

info@geniani.com


QUICK START GUIDE



Make sure that the massager is fully charged before use. When fully charged, the battery's color is constantly blue.




Screw in the head.
Use one of the bundled heads.



Press the Power / Modes button for 3 seconds to turn on/off. The massager will start working in Auto Mode. Press for 1 second to change the mode to Powerful. Press again for 1 second if you want to change the mode to Uniform Speed.

QUICK START GUIDE



a.



b.



c.

- Increase the vibration intensity
- Press the Power / Modes button for 3 seconds to turn on/off
- Decrease the vibration intensity



WARNING

For your safety and to avoid any damage to property, please read this manual carefully and follow it as indicated.

The non-compliance with the following notes can cause serious injuries or death and/or severe property damage.

Do not use the Geniani Cordless Massager in any other way, except the one described in the manual. This product is intended for personal use at home. It is a non-professional appliance designed to provide a soothing massage only.

The commercial use will cause premature breakdown of the product. Before using it, make sure the adapters and the product are compatible with the power supply.

SAFETY INSTRUCTIONS



Do not use it while near water or in damp environments. Do not use the massager with wet hands.



Do not use the product while driving a vehicle.



Keep the product out of the reach of children or pets.



Do not attempt to repair, disassemble, or modify the appliance as it has no user-serviceable parts inside.

SAFETY INSTRUCTIONS




Do not use the product on a heat-generating material. This may cause the batteries to overheat and cause product malfunction.



Do not disassemble the product and/or batteries. This may cause product malfunction, electric shock, or fire.



Use only the charger provided by Geniani. An unspecified charger or adapter may cause overheating, explosion, electric shock, or fire.



The recommended temperature for using the product is from 5°C (40°F) to 35°C (95°F). Use common sense when using or storing the Geniani Cordless Massager. Leaving it in overly hot places, such as inside a parked car on a hot day or next to a window with direct sunlight, may cause batteries malfunction or overheat.

SAFETY INSTRUCTIONS



Do not use the product for more than 15 minutes at a time. This is the usage time, recommended by doctors.



Do not insert any object into the charger input. This may cause damage to the product or electric shock.

If the Geniani Cordless Massager stops working due to overload or overheat, turn it off immediately. Let it cool down for 5 minutes and then reset the massager by plugging it into the charger for 3 seconds. The charger must be plugged into an electrical outlet.

1. Do not use it when folded.
2. Do not drop, use excessive force, or put heavy objects on the Geniani Cordless Massager to avoid malfunction and/or damages.
3. Do not take apart or modify the parts. This may cause electric shock and the manufacturer's warranty will no longer cover the product.

SAFETY INSTRUCTIONS

4. We do not recommend using it for personal sexual needs. This is a powerful massager that may cause severe bodily injuries.
5. Do not keep or use the product in areas where temperature changes drastically or where humidity is high. This may cause damage to electrical circuits and batteries.
6. If the product is not in use for an extended period of time, charge the unit once a month. If the batteries remain discharged for too long, this may degrade the product's life cycle and performance.
7. Do not keep the product near a high magnetic field or where it can be affected by a magnetic field, for it can cause batteries to discharge.
8. Do not use the Geniani Cordless Massager with an extension lead or switched timer.
9. Do not use your massager if it has malfunctions or has been damaged in any way. Have the appliance inspected by a qualified person.
10. Never cover the massager when it is in operation.

SAFETY INSTRUCTIONS

11. It is recommended that individuals with pacemakers should consult their doctor before using this product.
12. If you have any concerns regarding your health or suffering from an illness that may give you cause for concern, consult your doctor before using the product.
13. Never use this product on swollen or inflamed areas of skin.
14. Do not use this product as a substitute for medical help.
15. Do not use this product before going to bed as the Geniani Cordless Massager stimulates blood circulation and may delay sleep.
16. Always unplug the unit from the power source when the appliance is not used, or before cleaning.
17. If the supply cord is damaged, it must be replaced by a qualified person.
18. The adapter must be removed from the power socket before cleaning and having other maintenance procedures, or before moving the appliance.
19. Please make sure that the power cable is kept free from knots.

CONTRAINDICATIONS

Avoid using the Geniani Cordless Massager if experiencing any of the following medical conditions:

- Sensitive skin, swollen or inflamed areas; areas with poor circulation, skin eruptions, or if feeling unexplained calf or abdominal pain
- Pain in the throat area
- Pain in genital areas
- Unconscious or sleeping state
- Frostbite
- Skin irritability
- Deep vein thrombosis
- Recent burn
- Recent surgery
- Varicose veins
- Acute conditions requiring first aid or medical attention
- Acute flare-up of inflammatory conditions, such as rheumatoid arthritis
- Osteoarthritis, gout, or fibromyalgia
- Autoimmune conditions (lupus, scleroderma, multiple sclerosis, etc.)
- Hypertension
- Other conditions that affect your blood vessels (atherosclerosis, peripheral artery disease, deep vein thrombosis or arteriosclerosis)
- Osteoporosis (bone degeneration)
- Muscular dystrophy or other muscle disorders

CONTRAINDICATIONS

Do not use it during pregnancy. Pregnant should consult with their healthcare provider before using the massager.

If you have a pre-existing medical condition, including implants or pacemakers, consult with your physician prior to use.

Do not use it if you have blood clotting or are at risk for blood clotting.

Do not use it if you have bleeding disorders, low blood platelet counts, or are currently taking blood-thinning medications such as Warfarin. Consult with a physician prior to use.

Do not use it on the body areas with blood clots, fractures, open or healing wounds, skin infections, weakened bones (such as from osteoporosis or cancer), or where there has been a recent surgery.

Avoid any direct pressure on a tumor. Cancer patients should discuss any concerns about the massage therapy with their oncologist.


The Geniani Cordless Massager does not constitute medical treatment.

MASSAGE STICKS



U-Shaped Head Stick (ABS)

For calf and arm muscles



Four-Head Stick (ABS + Silicone)

For large areas, back, waist and legs;
delicate areas as well




Scalp Stick (ABS + Silicone)

Scalp stimulator for relaxation
and increased circulation



Cushion Stick (ABS + Silicone)

Soft, cushion for general massage



Six-Head Stick (ABS + Silicone)

Deep tissue/ sports massage for large muscle
groups



Point Stick (ABS)

Acupressure and reflexology for targeted
tension points and scar tissue

OPERATING INSTRUCTIONS

The Geniani Cordless Massager is a fast, rechargeable portable percussion massager. The massager is efficiently designed to be lightweight and yet yield outstanding performance and power.

The massager is powered by a lithium-ion battery and includes an AC charger. Also, it comes with six massage sticks and a storage pouch. It is a powerful cordless massager and should be used with caution in mind. Please read all the warnings before using it.

Power Button

The power button (Fig. 1) allows turning on the unit by pressing it for 3 seconds.

1st Press: Auto Mode

2nd Press: Powerful Mode

3rd Press: Uniform Speed Mode

+/- button controls the vibration intensity. When the unit is on, there are three levels of vibration: Low, Middle, and High.




Figure 1.
The power button

OPERATING INSTRUCTIONS

Charging

For the utmost convenience of use, charge the massager before using it.

1. Plug the charger into an electrical outlet.
2. Plug the charger into the Geniani Cordless Massager.
3. The charger has an LED light that will be flashing red or blue to indicate the massager's charge state. A blue light on the charger indicates a 100% charge.
4. Remove the charger from the massager before using it.
5. Now, the Geniani Cordless Massager is ready for use. Please check the operating instructions.

Battery Life

1. The battery's continuous use is up to 120 minutes (if no load is applied).
2. The battery tests were performed in the manufacturer's facility with its own equipment. The battery life depends on usage.

OPERATING INSTRUCTIONS

3. If the product is not in use for a long time, please turn it off. Batteries last longer if kept charged versus discharged.
4. When storing, keep the power off and store the product in a cold place.

Massage Sticks

The 4 threaded attachments (cushion, point stick, four-head back stick, six-head stick) screw on clockwise. Make sure you screw them on tightly but avoid overdoing.

The cushion stick provides a softer massage with its cushion absorbing mechanization. It can be used on most areas and is recommended for sensitive areas where there is more bone to muscle area, such as joints, shoulder, knee, or ankle.

OPERATING INSTRUCTIONS



For a deep tissue massage, use the point stick or six-head stick. It is recommended for areas where muscles are particularly stiff.

For a wider area massage, pick the six-head stick. If it hurts your skin, you can use it over clothes.

Enjoy the scalp stick for a soothing scalp massage. If you experience a bounce, try to switch to a higher speed setting or adjust your hand's position on the grip.

The U-Shaped Head Stick is the best choice for calf and arm muscles.

CASES OF USE



CARE AND MAINTENANCE

- 1.** Keep the massager in a safe, cold, and dry place when not in use.
- 2.** Always unplug the massager before cleaning.
- 3.** After use, wipe the massager with a soft cloth lightly dampened with a non-abrasive cleaner. If necessary, clean with a sterilizer.
- 4.** Do not submerge the Geniani Cordless Massager in water or any other liquid.
- 5.** Keep away from all solvents and harsh detergents.
- 6.** Do not dispose of the Geniani Cordless Massager in a landfill. It contains a Li-ion battery and other electrical components that are considered as hazardous E-Waste.

SPECIFICATIONS

Battery: 7.4V Lithium-ion 2200mAh 18650 x 2 Cell

Battery Indication:

Battery Energy <20%: Red Light

Battery Energy >20%: Blue Light

When charger:

Battery Energy <80%: Red Light is Blinking

Battery Energy >80%: Blue Light is Blinking

Battery Energy 100%: Blue Light

Charging Time: 120 - min charge lasts a week

Recommended Duration of Use: 20 - 30 min per day

Variable Speed Change: 1,500 rpm - 3,700 rpm

Percussion Massager Motor: 7.4V, DC, 3700 rpm

Charger Input Voltage: 100V-240V, 50/60Hz, 0.5A Max

Charger Output Voltage: 9V, DC, 1A

Contents: Massager, 6 Heads, Storage Pouch, AC adapter

TROUBLESHOOTING

The massager does not work.

1. Make sure the massager is fully charged. The power light should be blue.
2. Reset the massager by turning the power off. Plug the charger into an outlet, then plug the charger into the product. Wait for 3 seconds and unplug both sides of the charger. Turn the product on and start massaging.
3. The Lithium-ion battery may have reached its life cycle limit.



**Inquiries
Help
Issues**

1 (404) 225-3955

Monday to Friday 9 AM to 5 PM EST

info@geniani.com

2-YEAR WARRANTY

The manufacturer's warranty is void under the following circumstances but not limited to:

- Improper or inadequate maintenance or modification.
- Accident, misuse, abuse, contamination, or other external causes.
- The use of unspecified adapter and accessories.
- Loss or damage during the transportation.
- Damages that occur as a result of the user's failure to follow the instructions.

CUSTOMER SUPPORT



**Inquiries
Help
Issues**

1 (404) 225-3955

Monday to Friday 9 AM to 5 PM EST

info@geniani.com

*We are proud of our perfect
customer service, and we promise
to treat you like family.
Feel free to contact us!*

To get more information,
please visit our website
www.geniani.com



 SCAN ME