

# 40K fat exploding instrument R.F multipolar radio frequency User manual

## Product function introduction

## Precautions

## About maintenance

### RF radio frequency:

RF radio frequency acts on dissolving fat, detoxifying lymph, tightening skin, and enhancing skin elasticity. It can directly reach the fat layer. It has the function of directional radio frequency output. It accelerates the movement of fat cells to generate heat, and removes excess in the body through sweat glands, hepato-intestinal circulation and lymphatic circulation. Fat and toxins are excreted from the body to dissolve fat; When acting on the skin, the collagen fibers in the dermal layer can be heated to 45° -65° , and the collagen fibers shrink, so that the loose skin wrinkles are tightened. At the same time, the stimulated dermal layer collagen proliferation promotes the thickness and density of the dermis. , So as to achieve the purpose of removing wrinkles, removing scars, and restoring skin elasticity and luster to achieve firmness and wrinkle removal;

### Please read the following carefully before use:

1. Clean the skin first, then use special gel or essential oil;
2. When using the instrument, there must be an aqueous medium (such as slimming essential oil, etc.), otherwise there will be a certain loss of the instrument;
3. The operating probe should not stay in one place for too long, avoid the place with bones;
4. Do not use disinfectant to clean the product, it is recommended to use wet cotton or towel;
5. The operator and the person accepting the use take off the metal decorations;
6. If the machine is used continuously for 1 hour, please pause for 10 minutes before using it;
7. Some customers may experience tinnitus when using the fat-exploding probe, which is a positive phenomenon. Normal phenomenon, because this band has exceeded the range of human sound waves.
8. Keep it out of the reach of children;
9. Unplug the power when not in use.

After each use, it is recommended to clean the surface of the instrument. When cleaning, you must turn off the power and unplug the power plug. If you encounter stubborn stains, you can wipe it with a sponge or cloth moistened with a diluted neutral detergent. Please do not use strong alkali, Wipe with alcohol, volatile oil, delustering liquid, etc.  
This instrument is not waterproof, and you cannot use it in water or wash the host and power supply with water.

Note: Due to the irregular upgrade of the product, the model of the manual may be different from the actual product. This manual is for reference only.

### 40K fat exploding instrument:

The 40KHZ fat-blasting probe utilizes the 40,000HZ sound wave emitted by the fat-blasting probe that gathers strong sound waves. After entering the human body, the fat cells of the human body can produce a strong impact, which can effectively consume calories, water and fat cells to burst and make fat The cells are reduced to achieve the effect of removing fat.

## Not suitable for the crowd

1. If there is a wound on the skin, avoid the wound or do not do it when using it;
2. It is forbidden for people with epilepsy;
3. It is forbidden for patients with heart disease (such as arrhythmia), especially those with a pacemaker;
4. It is forbidden for cancer, nausea and tumor patients;
5. It is forbidden for patients with anxious sputum or infectious diseases;
6. It is forbidden for pregnant women, breastfeeding women, and during menstruation;
7. Serious diseases, such as heart disease, gastric ulcer, severe gastric disease, duodenal ulcer, high blood pressure, abdominal surgery less than three months old, diabetes, etc.;
8. It is forbidden to do it if there is metal or silica gel in the body;

## Warranty Card

1. This product enjoys a free warranty within one year from the date of purchase. It has been used in normal operation and has not been disassembled and repaired. After the warranty personnel confirm that it is a product quality problem, the warranty service is provided, and repairs are provided under any normal use.
2. The warranty card should be filled out at the time of purchase. The warranty card has nothing to do with our company's warranty is invalid. Please bring the warranty card when repairing the product.

Customer name						
The way of buying						
Product number						
Purchase date	year	month	day-	year	month	day
Warranty period	year	month	day-	year	month	day

## Product operation instructions

1. Before use, wash away the sebum and dirt accumulated on the surface of the facial skin with warm water, so that the skin with open pores can better absorb nutrients.
2. Connect the handle to the corresponding socket, connect the power supply and turn on by default, press the Menu button to select the function, then press the Timer button to adjust the working time, and then press Level+ or Level- to adjust the appropriate working intensity and start working. When the time countdown ends, it will automatically stop working.

## Instrument parameters

Power input voltage: 100V-240V  
Power output voltage: 24V  
Whole machine power: 50W  
RF output frequency: 5M  
Number of handles: 3  
Screen: Led digital display  
Gear position: 7 gears  
3 modes of radio frequency: body, face, eyes  
Time: 5, 10, 15, 20, 25, 30 minutes  
Ultrasonic fat explosion working frequency: about 40kHz  
Three modes of ultrasonic burst fat: Mode:A/Mode:B/Mode:C



LW-111

- Thank you very much for your purchase
- Please read the manual carefully before use for normal operation. After reading it, please keep it properly for inspection at any time.
- The use of this beauty device is based on the care of healthy skin and for the purpose of beauty. It is prohibited to use or sell it for the purpose of treatment.

## Part name introduction

