


# Belmint

PRESS ESCAPE™



 [www.belmint.com](http://www.belmint.com)

 [support@belmint.com](mailto:support@belmint.com)

 1 (844) 423-5646

 [@belminproducts](https://www.facebook.com/belminproducts)

 [@belminproducts](https://www.instagram.com/belminproducts)

## FULL BACK MASSAGER

## ► FULL BACK MASSAGER

### User's Guide

 <b>Model</b>	 <b>Power Consumption</b>	 <b>Voltage</b>
BEL-FBM	30 W	DC 12V   2.5A

### Congratulations on your purchase!

Transform any seat in your home into a luxurious massage experience. Feel the tension leave your body every time you sit down. Enjoy a spa getaway any time you want.

Is the pressure at work getting to you? Just lean back and enjoy the relaxing and soothing feeling of full back massage with the gently rotating nodes replicating the movements used by professional spa reflexologists.

Do you feel a slight pinch anywhere in your back? The full back massager works on your lumbar area with each node strategically positioned to target specific points on your back similar to a shiatsu massage. The massager targets your upper back, lower back and your waist to reduce the tension and pinching sensation.

All Belmint products are designed by Belmint engineers and crafted from high quality, durable materials to ensure maximum comfort with kneading nodes that target hard-to-reach areas on your back. So, go on. Sit back, read a book or watch a movie! Let our Belmint Full Back Massager relax and soothe you.

To celebrate your comfort,  
**The Belmint Team**

## ► Sometimes peace is hard to find

### But good service shouldn't be

At Belmint, we strive to provide our customers with ultimate relaxation. If you encounter any issues with the performance of your device, our authorized service personnel are here to help. Check back cover for information.

## ► What's in the Box?

- 1 Full Back Massager
- 1 User's Guide
- 1 Power Cord
- 1 Car Adaptor

## ► Product Details

### Dimensions

24.5" x 16" x 4"

### Weight

7 lbs.

### Heat

98 °F — 118 °F

### Auto Shut-Off

15 minutes



## ► **Parts & Features**

### **Customization Options**

**Unit ON, Heat ON** - Click Power Button

**Unit ON, Heat OFF** - Click Power Button

**Unit OFF** - Click Power button for 3 seconds

**Clockwise and Anti-Clockwise Mode** - Click program button to change between clockwise and counter-clockwise motion of the massage nodes

### **ERGONOMIC DESIGN FOR COMFORT**

The unique design is engineered to ensure maximum comfort with the 12 massage nodes strategically placed against the back's pressure points and designed to mimic a real massage therapist's hands.

### **SWITCHABLE HEAT**

The switchable heat function promotes an increase in blood and oxygen flow to soothe your muscles

### **CHARACTERISTIC FUNCTIONS**

- Easy control with one button.
- 12 massage nodes ensure that the entirety of the back is well-massaged
- Lightweight and easy to carry and store.
- 15-minute auto-shutoff

## ► **Let the Bliss Begin!**

### **Using your device**

1. Connect the AC adapter to the plug
2. Connect the adapter end of the plug to the unit outlet.
3. Plug your Belmint Full Back Massager into the appropriate outlet. If the plug does not fit, please contact Belmint Customer Service (Information on the back.) DO NOT alter the plug in any way.
4. Make sure you have the massage secured firmly in place on your seat. Do not use your massager on a plush surface to avoid overheating.
5. Gently sit and place your back against the massager.
6. Turn the unit ON by pressing the POWER button.
7. Upon turning ON, the massage nodes will immediately start to rotate. The HEAT is set to ON (the massage nodes will glow RED).
8. To turn the HEAT OFF, press the POWER button again. When the HEAT is OFF, the nodes will no longer glow RED.
9. Turn POWER OFF by pressing the POWER button for 3 seconds.

**NOTE:** Please do not use the back massager for more than 15 minutes at a time



### **Before things get awesome..**

To prevent potential injury or property damage, please read this entire manual before operating your device. Keep it on file for future reference.

## ► **Cleaning Your Device**

### It is as easy as 1-2-3

To clean the outer shell of your massager, use a soft, slightly damp sponge. Do not allow water or any liquid to come in contact with the unit.

**DO NOT** use any chemicals or harsh abrasives such as cleaners, brushes, gasoline, kerosene, glass/furniture polish or paint thinner to clean.

**NEVER** immerse your back massager or any electrical appliance in water.

**NEVER** clean when the unit is on.

**NOTE:** Make sure the massager is turned off and unplugged from the adapter before cleaning.

Allow the massager to cool down before cleaning.

## ► **Storing Your device**

Unplug the cord from the power source.

Place the massager, along with the cord, in its box or in a safe, dry and cool place.

**DO NOT** wrap the power cord.

**AVOID** contact with sharp edges or pointed objects.

Your relaxation and peace of mind is very important to us. That's why your Belmint Full Back Massager is made from high quality durable materials.

## ► **Keeping The Peace** Troubleshooting Guide

CONDITION	REASON	SOLUTION
Massager won't start or turn ON	Unit may be unplugged	Check that the unit is plugged into the appropriate outlet or power source. Check that the unit is plugged into the adapter. Press the POWER button. You should feel it click
Massager will NOT turn OFF	Unit may be on wrong MODE	Check to see if unit is on HEAT OFF mode. If so, press POWER button for 3 seconds to turn OFF.
Massager will NOT turn OFF	Unit may be on wrong MODE	Check that the unit is plugged into the appropriate outlet or power source. It may have come unplugged. Check that the adapter is plugged into the unit
	Unit may have overheated	Unplug the unit from the outlet and let it cool before using again.
	Unit may automatically shut off.	Press the POWER button. The unit may have shut OFF due to the automatic timer.