

▼  
**NECK  
MASSAGER**  
▲



[www.belmint.com](http://www.belmint.com)



[support@belmint.com](mailto:support@belmint.com)



1 (844) 423-5646



@belmintproducts



@belmintproducts

# NECK MASSAGER



## ▶ User Guide

 <b>Model</b>	 <b>Voltage</b>	 <b>Power Consumption</b>
BEL-NECK	12V	24 W

### Congratulations on your purchase!

You're on your way to Spa-level experience

Your Belmint Neck Massager electronic shiatsu massager with heat makes it possible for you to get a spa-level massage anywhere, anytime.

Like all Belmint products, your neck massager is crafted with high-quality, durable materials to ensure your comfort. Specially designed to ease tension in tired and sore muscles, your massager is equipped with rolling and kneading nodes that enables you to target those hard-to-reach areas of your body. For a truly custom massage experience, your Belmint neck massager is equipped with optional heat, directional control and nodular speed control—because some days we need more kneading than others. Cleanup and storage are easy. Simply wipe with a damp cloth after use and store it in its included carry case.

You won't want another neck massager. You've got a Belmint

To celebrating your comfort,  
The Belmint Team

## ▶ Sometimes peace is hard to find But good service shouldn't be

At Belmint, we strive to provide our customers with ultimate relaxation.

If you encounter any issues with the performance of your device, our authorized service personnel are here to help. Check back cover for information.

## ▶ What's in the Box?

- 1 Neck Massager
- 1 User Manual
- 1 AC/DC Plug Adapter

## ▶ Product Details

### Dimensions

14 x 6.8 x 5.3 inches

### Weight

3.5 pound

### Auto shut-off

15 minutes

## ▶ How to Use

1. Plug-in appropriate Power adapter when using at home or in a vehicle
2. Press the button on the controller to turn on the massager.
3. Use the Control Panel to select your Desired Massage Options
4. After use, Power off the Massager and unplug the power adapter from the socket or power supply

