

**FOOT & CALF
AIR COMPRESSOR**



www.belmint.com



support@belmint.com



1 (844) 423-5646



@belmintproducts



@belmintproducts

**FOOT & CALF
AIR COMPRESSOR**



Model	Voltage	Power	Default Time
BEL-ACWH	12V	18W	20 minutes

Congratulations on your purchase!

You're on your way to better self-care

Are you feeling particularly sore in your foot and calf area after a long walk? Do you work where you are required to stand for long periods of time. Are you engaged in activities requiring you to put pressure in your legs and feet?

Then, the foot and calf air compressor is for you! With the Belmint Foot and Calf Air Compressor's lightweight and portable features, you don't have to wait to get home to feel its relaxing effects. Bring it with you wherever you need to go to enjoy its relaxing effects!

So, go on. Sit back, read a book, or watch a movie! Let your Belmint Foot and Calf Air Compressors Relive your Pain

To celebrating your comfort,

The Belmint Team

► Quick User Guide

Box Content

2 x Foot & Calf Air Compressor | Control Adapter | User Manual

How to Use

1. Wrap the air compressor bags around your Legs and Feet
2. Connect the Two Air Pressure Tubes to the controller unit
3. Plug into the wall



Power Button

1. Press to Start
2. Press again to adjust intensity between LOW - MEDIUM - HIGH
3. Press & Hold 5 Seconds to shut off

Heat Button

Choose between 3 Levels of heat or turn off

Mode Button

Choose between 6 different Modes

1. Auto
2. Leg Only
3. Foot Only
4. Pinch
5. Press
6. Relax

► Cleaning & storing your device

1. For best results, clean and store your device after each use.
2. Always unplug the device and allow it to cool completely before cleaning.
3. keep the device air tubes & Connectors free of dirt, dust, and other obstructions
4. Use a clean, dry cloth to wipe the outside of the device
5. Do not use any harsh chemicals or cleaning products