





▼  
**FOOT & CALF  
AIR COMPRESSOR**  
▲



# FOOT & CALF AIR COMPRESSOR



## ▶ User Guide

 Model	 Voltage	 Battery	 Rated time
BEL-FAC	DC 3.7V	2000 mA Li	40-60 minutes

### Congratulations on your purchase!

#### Your footsies are in for a treat!

Now, you can enjoy a luxuries massage especially targeting your foot and calves!

Are you feeling particularly sore in your foot and calf area after a long walk?

Do you work where you are required to stand for long periods of time?

Are you engaged in activities requiring you to put pressure in your legs and feet?

Then, the foot and calf air compressor is for you!

With the Belmint Foot and Calf Air Compressor's lightweight and portable features, you don't have to wait to get home to feel its relaxing effects. Bring it with you wherever you need to go to enjoy its relaxing effects!

All Belmint products are designed by Belmint engineers and crafted from high quality, durable materials to ensure maximum comfort and ease of use.

So, go on. Sit back, read a book or watch a movie! Let our Belmint Foot and Calf Air Compressors soothe you.

To celebrating your comfort,

The Belmint Team

## ▶ Sometimes peace is hard to find

But good service shouldn't be

At Belmint, we strive to provide our customers with ultimate relaxation.

If you encounter any issues with the performance of your device, our authorized service personnel are here to help.



Phone

**1-844-423-5646**



Email

**support@belmint.com**



Website

**www.belmint.com**



Instagram

**@belmintproducts**



Facebook

**@belmintproducts**

## ► What's in the Box?

- 2 Foot & Calf Air Compressor Bags
- 1 Remote Control
- 1 Charger
- 1 User Guide

## ► Product Details



**Note :** The Foot and Calf Air Compressor will automatically Shut Off in 15 minutes.

## ► Getting Started

### How to use

1. Carefully wrap the air compressor bags around your feet and calves.
2. Plug your Belmint Foot and Calf Air Compressor into the appropriate outlet. If the plug does not fit, please contact customer support.
3. Press the POWER button.
4. Adjust the air pressure by clicking the air pressure intensity button.
5. Click Mode 1 for the manual mode or Mode 2 for the automatic mode.
6. Upon turning ON, the massager nodes will immediately start to rotate. The HEAT is set to ON. (the massage nodes will glow RED).
7. Press POWER button again to turn the massager off.
8. Charging time is 2 hours. For continuous operation, make sure the battery is fully charged before use.

## ► Massage Modes

### 2 Types of Mode

**Mode 1 :** This massage mode is the manual mode with 2 alternating positions in the massage

**Mode 2 :** This massage mode is the auto mode with 4 positions of circular-motion massage.

## ► Air Pressure Intensity

### 3 Levels of Intensity

**Low :** The corresponding LED light is GREEN.

This level of intensity is perfect for those who are sensitive to pressure.

**Middle :** The corresponding LED light is YELLOW.

This level of intensity is perfect for those who are fine with more pressure in their massage.

**High :** The corresponding LED light is RED.

This level of intensity is for those who prefer the highest pressure setting.

## ► Efficiencies

Keep using **1 OR 2 TIMES PER DAY** to relax muscles & improve the circulation; The intensity level should be set where you will feel involuntary muscle movement & vibrations in the muscle.

## ► Features

### Ergonomic Design for Comfort

This unique design is engineered to ensure maximum comfort every time you use the massager.

### Rechargeable

The Foot and Calf Air Compressor's rechargeable function assures that you do not need your massager tethered into a cable as long as you make sure it's fully charged.

### Characteristic Function

Easy to use remote control

3 levels of air pressure intensity to suit your preference

Lightweight and easy to carry and store

15-minutes auto-shut off

## ► Storing your Belmint Foot & Calf Air Compressor

Unplug the cord from the power source.

Place the massager, along with the cord, in its box or in a safe, dry and cool place.

DO NOT wrap the power cord.

AVOID contact with sharp edges or pointed objects.

## ► Cleaning your Belmint Foot and Calf Air Compressor

Your relaxation and peace of mind is very important to us. That's why your Belmint Foot & Calf Air Compressor is made from high-quality durable materials.

### Cleaning your unit is as easy as 1-2-3

To clean the outer shell of your massager, use a soft, slightly damp sponge. DO NOT allow water or any liquid to come in contact with the unit.

DO NOT use any chemicals or harsh abrasives such as cleaners, brushes, gasoline, kerosene, glass/furniture polish or paint thinner to clean.

NEVER immerse your back massager or any electrical appliance in water.

NEVER clean when unit is on.

Make sure the massager is turned OFF and unplugged from the adapter before cleaning.

Allow the massager to cool down before cleaning.