User Guide: Back Stretcher (Waist Relax Mat)

3 Stretch Levels

The Back Stretcher has 3 different stretch levels, which can be manually adjusted:

- Level 1: For initial user with moderate stretching
- Level 2: For intermediate user with medium stretching
- Level 3: For advanced users with strong stretching

How to use

- 1) Choose the right level for you
- 2) Sit on floor with
- 3) Place the Back Stretcher right behind you
- 4) Slowly go from sitting position towards a lying position over the BackStretcher
- 5) Relax and stretch for 5-10 min.

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Intended Use:

The Back Stretcher is a device which is intended to use for stretching the back and waist of adult users. The Back Stretcher is intended to be used for 5-10 minutes per session. Sessions can be repeated.

Precautions:

When using the Back Stretcher, the users must monitor the effect the product has on their body to avoid risk of injury.

In case of (but not limited to) discomfort or pain, stop using the Back Stretcher and consult a physician immediately.

In case of recent surgeries, back injuries or pregnancy do not use the product before consulting a physician.

Do not use this product within 20min after meals.

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