*Innovative Ideas That Make Life Easier™*

# **Hot & Cold Therapeutic Point Massager**

**Model No. HCM01**

**Important Safety Instructions**

# Important Safety Instructions

## READ ALL INSTRUCTIONS BEFORE USE

### FEATURES

- Heats to 102° F\*
- Cools to 57° F\*
- Simultaneously Massage While Heating or Cooling
- Convenient Portable Size is Great for Travel
- Eliminates Hassled Need for Ice Packs and Heating Pads
- Uses 4 (AA) Batteries (Not Included) or USB Cable (Included)

\*At Ambient Temperature Level of 70°

### OPERATING INSTRUCTIONS

1. Press the Power Button to turn the unit ON (The Light will Turn Green)

#### HOT & COLD INSTRUCTIONS

Press the C/H Button to Switch between Hot and Cold

- Green Light will Light when not in Hot or Cold Mode
- Blue Light will Light When in Cold Mode
- Red Light will Light When in Hot Mode

Hold the Silver Circular surface to the desired area

#### MESSAGE MODE INSTRUCTIONS

1. To use the Massage Mode feature press the Vib Button.

2. Press the Vib button again, to turn the Vibration off.

2. Press the Power Button to turn the unit OFF.

## **USB POWER**

1. To operate the unit via USB power, connect the USB cable to a powered USB port.
2. Connect the power tip to the power port at the bottom of the unit. Your unit is now ready to operate.

## **DANGER/WARNING**

To reduce the risk of burns, fire, electric shock, or injury to persons:

- Do not block fan area while in use.
- Use this unit only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this unit if it is not working properly, if it has been dropped or damaged, or dropped into water. Return unit to

manufacturer for examination and repair.

- Never drop or insert any object into any opening.

## **CARE AND MAINTENANCE**

1. Use only a soft damp cloth or non-abrasive cleaner to clean.

## **BATTERY REPLACEMENT**

The battery compartment is located on the bottom side of the unit. To replace the batteries, use these instructions:

1. Before attempting to replace batteries, turn unit off.
2. Turn the unit upside-down and pull the battery compartment cover to open the battery cover from the closed to open position.
3. The compartment holds four (4) AA batteries. All four batteries are required and must be properly positioned for the unit to operate.
4. Once the batteries have been inserted, replace the battery compartment cover.

5. Your unit is now ready for you to use.

## CARE/DISPOSAL OF ALKALINE BATTERIES

1. Use only fresh batteries of the required type and size.
2. Do not mix old and new batteries, different types of batteries, or rechargeable batteries of different capacities.
3. The batteries may explode if mistreated. Do not attempt to recharge or disassemble the batteries.
4. Depleted batteries are to be removed from the item.
5. Dispose of the old batteries properly, and according to your local regulations. Ensure that batteries are kept out of reach of young children.

## 90 Day Limited Warranty

Zadro products include a 90-day limited warranty from date of purchase against defects in materials and workmanship. This warranty does not include damage to the product as a result of misuse or accident.

Should your product become defective during the warranty period, Zadro Products, Inc.™ will repair or replace (at our sole discretion) your product, at no charge. However, you are responsible for properly packing your product and any applicable shipping charges to Zadro Products, Inc.

Please contact our Customer Service at 800-468-4348 or email [customerservice@zadroinc.com](mailto:customerservice@zadroinc.com) on how to obtain an RMA (return authorization number). Please have available and/or include a copy of your receipt, model number (found above the barcode on the packaging) and the issue you are experiencing with your product.

Any product returned for preference is subject to a 20% restocking fee of the product and no shipping charges will be refunded.

Zadro Products, Inc.™ shall not be liable for any incidental or consequential damages arising out of defects in this product, whether or not caused by its negligence.

Some states do not allow the exclusion or limitation of incidental or consequential damages or allow limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights and you may also have other rights, which vary from state to state.

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© 2015 Zadro Health Solutions