

AGELOC BOOST™



Let your skin glow brighter with an innovative device system from Nu Skin. ageLOC Boost features unique, Variable Pulse Microcurrent Technology and potent ingredients to visibly brighten and awaken your skin. Now that's simply brilliant!

MEET AGELOC BOOST

Say goodbye to dull, tired-looking skin and hello to skin that glows. Introducing ageLOC Boost, a breakthrough beauty device system that brings out skin's brilliant best. We're talking visibly brighter, plumper, bouncier skin with youthful radiance and a more even-toned complexion. Effortless to use (just two minutes a day), ageLOC Boost gives you a beauty boost with every use.

WHY YOU'LL LOVE IT

- Boosts your skin's natural beauty for brighter, plumper, and bouncier skin.
- Promotes cellular turnover, encouraging a more radiant, youthful complexion.
- Provides intense brightening benefits for a vibrant, glowing appearance.
- Helps ease signs of premature aging and prolongs youthful-looking skin.
- Offers gentle exfoliation to reveal visibly smoother, luminous skin.
- Hydrates, so skin feels softer and appears fuller.

WHY YOU'LL LOVE IT EVEN MORE

- Features patent-pending Variable Pulse Microcurrent Technology.
- Delivers lasting benefits.
- Offers spa-like results at home.
- Keeps you at your selfie- and video-chat-ready best!
- Super easy and quick to use (just two minutes a day).
- Effortless add-on to most skin care routines.
- Features a bottle made with 15% post-consumer resin (PCR).

WHAT POWERS IT?

- **Brightening Blend**—this concentration of brightening ingredients promotes brighter, more radiant-looking skin.
- **Caviar Lime Fruit Extract**—encourages exfoliation and increased cellular turnover for visibly smoother, even-toned skin.
- **Plumping Blend**—hydrates and softens, while giving skin a fuller, plumper appearance.
- **ageLOC® Blend**—Nu Skin's exclusive technology that targets the sources of aging to help maintain a youthful, vibrant appearance.

WHAT YOU'LL RECEIVE

- ageLOC Boost device with magnetic charger
- Activating Serum

WHAT'S VARIABLE PULSE MICROCURRENT TECHNOLOGY?

ageLOC Boost sends intermittent, variable pulses to the skin, eight times a second. The pulses are charge-balanced, so they switch from positive to negative every 2.4 seconds. Why? When used with Activating Serum, the microcurrent optimizes and activates skin for a more youthful complexion.

AGELOC BOOST ACTIVATING SERUM—ADDING OOMPH TO YOUR SKIN CARE

Formulated to work with the device, ageLOC Boost Activating Serum will help you achieve a glowing, youthful complexion. It contains efficacious brightening ingredients, a hydrating plumping blend, a gentle exfoliator, and our proprietary ageLOC Blend to address the early signs of


AGELOC BOOST™

aging at the source. The outcome? Visibly bouncier skin with a healthy glow.


CLINICALLY PROVEN RESULTS

Nu Skin commissioned a third-party clinical study to evaluate the efficacy of ageLOC Boost (device and ageLOC Boost Activating Serum). Study participants, 25-40 years-old, used the system once a day for 12 weeks, and the study included clinical grading and self-perception evaluations.

12-WEEK CLINICAL GRADER RESULTS



12-WEEK SELF-PERCEPTION RESULTS



HOW TO USE IT

Super easy to use, ageLOC Boost transforms your skin from blah to brilliant. After cleansing and toning, treat one half of your face for one minute with ageLOC Boost and ageLOC Activating Serum, and then treat the other half for another minute. That's it—you're done! Follow up with your daily serum, moisturizer, and sunscreen. You'll immediately notice your skin feels more hydrated and softer. Please note: only use ageLOC Boost Activating Serum with the ageLOC Boost device, since it was scientifically formulated to conduct the unique Variable Pulse Microcurrent. Store ageLOC Boost Activating Serum in a cool, dry place, out of direct sunlight.

SMART DESIGN

To provide the best user experience, we created ageLOC Boost in partnership with an award-winning, world-renowned design firm. The device features an ergonomic feel that's comfortable and easy to use.

LEARN MORE ABOUT IT

When should I use ageLOC Boost?

You can use ageLOC Boost anytime. However, we recommend using it as part of your morning skin care routine, so you'll get a beauty boost that lasts throughout the day.

How are ageLOC Boost and ageLOC Galvanic Spa® different?

You can use ageLOC Boost every day, whereas you should use ageLOC Galvanic Spa two to three times a week. The main difference between the two devices is they use a different type of current. ageLOC Boost uses exclusive Variable Pulse Microcurrents to optimize skin, which help reveal a more youthful complexion. ageLOC Galvanic Spa uses direct microcurrents to deliver treatment products to visibly improve aged skin.

Can I use ageLOC Boost in the shower?

You should not use ageLOC Boost in the shower or bath because the device is not waterproof (it's water-resistant). Also, you'll want to use the system after showering because Activating Serum is a leave-on product, not one you wash off.

How does ageLOC Boost support Nu Skin's commitment to become more sustainable?

ageLOC Boost's Activating Serum bottle is made with 15% post-consumer resin (PCR), and the product and device cartons are made from virgin paper. Please check your local capabilities and recycle the bottle and cartons where possible. By purchasing ageLOC Boost, you're joining us in taking an active, engaged stride to ensure the future is bright for both the planet and each of us. Please visit www.nuskin.com/sustainability for more information and details on how to recycle.

Can I use any product with ageLOC Boost?

To enjoy the best treatment experience, you should only use ageLOC Boost Activating Serum with your device. That's because we scientifically formulated Activating Serum to optimally conduct the unique Variable Pulse Microcurrent.

AGELOC BOOST™

Why do I need to put Activating Serum on the device's treatment head, rather than directly on my skin?

Dry hands are a must when using the device, which is why we've developed a treatment protocol that lets you dispense ageLOC Boost Activating Serum onto the device head, so your hands remain dry.

WHAT'S IN IT

Water (Aqua), Glycerin, Butylene Glycol, Dimethicone, Polyglyceryl-6 Distearate, Tetrahexyldecyl Ascorbate, Niacinamide, Sodium Acetylated Hyaluronate, Sodium Hyaluronate, Sodium PCA, Citrullus Lanatus (Watermelon) Fruit Extract, Lens Esculenta (Lentil) Fruit Extract, Microcitrus Australasica Fruit Extract, Pyrus Malus (Apple) Fruit Extract, Pancratium Maritimum Extract, Sodium Lactate, Ergothioneine, Tocopheryl Acetate, Jojoba Esters, Cetyl Alcohol, Lecithin, Mannitol, Polyglyceryl-6 Polyricinoleate, Polyglyceryl-3 Beeswax, Xanthan Gum, Hydroxypropyl Methylcellulose Stearoyl Ether, Sodium Acrylates Copolymer, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Propanediol, Chlorphenesin, Ethylhexylglycerin, Phenoxyethanol, Benzoic Acid, Tetrasodium Glutamate Diacetate, Aminomethyl Propanol.