

THE FIRST STEP TO A HEALTHY LIFE

Thank you for purchasing the USJ-815 VITALITY SWING.
Please read these instructions carefully before operating this machine.
Please keep this instruction manual in a safe place, along with the purchase receipt for future reference.

Note: Due to constant upgrades, you may find slight differences between the design and functions states in this manual and those on the actual U.S. Jaclean, Inc. device.

IMPORTANT SAFETY INSTRUCTIONS

- Use the unit under its intended purposes only, as described in this manual. Do not attach any components that are not specifically made for this product.
- Never carry or drag the unit by its cord. Do not use the cord as a handle. At the end of each massage session, please remember to disconnect the power plug.
- Keep the cord away from heated surfaces at all times.
- Do not use the unit within 30 minutes before and after meals.
- Individuals undergoing the following conditions should consult a doctor before using this unit:

Skin disease	Heart disease
Recent surgery	Unhealed fractures
Pregnancy	

FEATURES

- Extra heavy-duty motor to handle heavier bodies
- 15-minute automatic shut-off timer
- Extra padded vinyl ankle cradle for support
- Easy-to-use remote control
- Adjustable speed

HEALTH BENEFITS

- Increases metabolism and promotes weight loss
- Increases oxygen to the cells
- Promotes cell repair and rejuvenation
- Relieves back pain and relaxes muscles
- Strengthens the immune system
- Eliminates fatigue and improves sleep

OPERATION

1. Place the unit on a hard, flat surface, such as the floor or a mattress. Insert the power plug into an electrical outlet.
2. Lie down on your back and place both ankles into the unit's cradles.
3. To turn on the unit, flip the power switch on one corner of the unit to the "on" position.
4. Use the round dial on the handheld controller to adjust your desired speed of massage. We recommend that you start each session with a slower swinging speed, and increase your speed as your body gets used to the movements.
5. At the end of each session, flip off the power switch and remain in a lying position for 2-3 minutes before getting up.
6. Drink a glass of water before and after each session. Immediately following a meal, wait 30 minutes before using the unit.

MAINTENANCE

1. Disconnect the power plug from the wall socket before cleaning the unit. Wipe all surfaces with a soft, damp cloth. Do not immerse the unit into any liquids for cleaning.
2. When the unit is not being used, store it in a dry location. Do not store near intense heat or an open flame.
3. Do not leave this unit under direct sunlight for prolonged periods of time.
4. Do not twist or wrap the power cord around the unit.
5. Keep the unit out of the reach of pets and young children.

SPECIFICATIONS

Model #:	USJ-815
Power Supply:	AC 120V- 60Hz
Power Consumption:	40W
Weight:	18 lbs.
Dimensions:	16 x 11 x 15 inches

USJ-815 VITALITY SWING

Consumer Limited Warranty: ONE-YEAR WARRANTY

This warranty guarantees that the product purchased is to be free from defects in material and workmanship under normal usage and local voltage for a period of one year. This warranty expires in **TWO (2) years on parts, and ONE (1) year on labor** after the date of purchase.

U.S. Jaclean, Inc. (referred to as "the Warrantor") will repair this product with new or rebuilt parts in case of defects in material or workmanship, free of charge, in the USA provided that:

- 1) The product is operated in accordance with the printed instructions supplied.
- 2) The product is used under normal conditions for which designed.
- 3) The product is not subjected to misuse, negligence, or accident.
- 4) The product receives proper care, protection, and maintenance.

This warranty is extended to the original end-user when purchased from an authorized U.S. Jaclean, Inc. distributor. A purchase receipt or other proof of date of the original purchase is required for service and parts replacement under this warranty. This warranty does not cover shipping charges for defective items if a claim is made more than 30 days after the date of purchase.

This warranty only covers failures due to defects in materials and workmanship, and does not cover normal wear or cosmetic damages. This warranty does not cover damages which occur in shipment, or failures which result from: accident, misuse (including, without limitation, failure to comply with instructions and use limitations contained in the Operating Instructions), abuse, neglect, mishandling, faulty installation, set-up adjustments, improper operation or maintenance, alteration, modification, power line surge, lightening, improper line voltage, commercial use such as hotel, office, restaurant, or other business or rental use of product, service by anyone other than from U.S. Jaclean, Inc., or damage that is attributable to "acts of God."

Continued use of this product after discovering a defect voids this warranty. This warranty shall be governed by the state of California, United States of America.

IMPORTANT: DO NOT return defective product to the store.

Contact U.S. Jaclean, Inc. for repair services.

Please retain this Warranty Certificate, together with the original purchase receipt, and present it to our technician for inspection during service.

For assistance and more information, please contact U.S. Jaclean, Inc.

VITALITY SWING

DELUXE MODEL

AEROBIC EXERCISER



SAVE THESE INSTRUCTIONS

Thank you for purchasing the Vitality Swing
Please read this manual carefully to ensure optimum performance and safety.

INSTRUCTION MANUAL