

Water Flosser EWFO03

Q1: How does a cordless water dental flosser work for teeth?

A cordless water dental flosser works by using a pulsating stream of water to effectively remove plaque and food debris from between the teeth and along the gumline. It directs a jet of water into hard-to-reach areas, providing a thorough cleaning and gum massage.

Q2: What are the advantages of using a cordless water dental flosser

over traditional flossing methods?

A cordless water dental flosser offers several advantages over traditional flossing methods. Firstly, it provides more effective plaque and debris removal from between the teeth and along the gumline. Secondly, it is gentle on the gums and can contribute to improved gum health. Additionally, it is easier to use, especially for individuals with limited dexterity or those wearing braces.

Q3: Can a cordless water dental flosser replace regular brushing?

No, a cordless water dental flosser should not replace regular brushing. It is essential to brush your teeth with a toothbrush and toothpaste twice a day to remove plaque and maintain good oral hygiene. The cordless water dental flosser is a supplementary tool that can enhance your oral care routine when used in conjunction with regular brushing.

Q4: How should I aim the jet stream of the flosser to properly clean my

teeth and gums?

To achieve effective cleaning, it is recommended to direct the jet stream of the cordless water dental flosser at a 90° angle to your gumline. This angle helps efficiently remove plaque and debris from the gum line and between teeth, promoting better oral health.

Q5: Can I use the cordless water dental flosser while it is charging?

No, it is not possible to use the cordless water dental flosser while it is charging. The product cannot be used during the charging process. To ensure proper usage and avoid any potential safety issues, it is recommended to wait until the device is fully charged before using it. The charging indicator light will turn green when the flosser is fully charged.

We hope these tips were useful. If you have any inquiries about our products or services feel free to contact us via email, call our customer service at 855.822.6999 M-F 9am-5pm CT or visit us at healthcare-manager.com