Smart Body Fat Scale

Lescape P3

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THANK YOU FOR YOUR PURCHASE!

The Scale uses bio-electrical impedance analysis (BIA) technology to help you conveniently measure body weight, body fat rate, visceral fat, body water, skeletal muscle rate, muscle rate, muscle mass, bone mass, protein, BMR, body age and more.

Scan the QR code to access the latest instructions or visit our website www.lepulsefit.com to get more information and support.


WARNING

- DO NOT use the scale if you have implanted medical devices such as pacemakers.
- DO NOT stand on the edge of the scale or jump on it.
- DO NOT overload the scale (Max 396lb/180kg/28st).
- DO NOT drop the scale or drop objects on it as this may damage the sensors.
- DO NOT immerse the scale in water or use chemical cleaning agents. Clean the scale with a slightly damp cloth.
- Any measurement obtained using this device is for reference only and should not be considered as a medical opinion.
- DO NOT use to diagnose or treat any medical condition. You should consult your primary care physician before making changes to your diet, exercise plan or physical activities.
- Always place the scale on a hard, dry and flat surface before measurement.
- Make sure your feet are dry before stepping on the scale.
- DO NOT use the device if damaged. The continuous use of a damaged unit may cause injury or improper results.
- Please check the device before each use.
- Be careful when using on the wet and slippery surface.
- For people under 16 years old, the scale can only measure body weight.

PRODUCT SPECIFICATIONS

Product size:	310x355x50mm
LED screen display:	68x109mm
Weight unit:	lb/kg/st
Weight limit:	111lb-396lb/5kg-180kg
Weight division:	0.2lb/0.1kg
Power Supply:	300mAH lithium battery
Operating temperature:	10-40°C
Operation humidity:	40%-80% RH

PRODUCT OVERVIEW

DISPLAY SCREEN

Note: The pictures shown in this manual are only for illustration. Please refer to the actual product for details.

BMI TABLE

BMI (Suitable for men and women)					
Age	Underweight	Healthy weight	Overweight	Obesity	Severe obesity
18-100	<18.5	18.5-23.9	24-27.9	28-34.9	≥35

Body Mass Index (BMI) = Weight (kg) ÷ Height² (m)


For example: BMI = 70kg ÷ (1.75m * 1.75m) = 22.86 In the range of healthy weight

Note: The BMI index is a commonly used standard to measure body weight and health status.

POWER ON

For first use

Press the reset/power button or charge the scale to cancel the shipping mode and turn on the scale.

For subsequent use

Press to wake the body composition scale. The pressure should be greater than the minimum range of the scale (5 kg).

Download and Install the App

- 1 Search Fitdays in the App Store or Google Play, or scan the QR code below to download and install the Fitdays App.


- 2 Register an App account: Enter the email address and password, select [I have read and agree Privacy Policy], and tap [NEXT].

Note: The feature of creating a strong password for iOS is not supported now.

A screenshot of the 'Email Registration' screen in the Fitdays app. The screen has a white background with a light blue border. At the top, it says 'Email Registration'. Below that, it asks 'Please enter your Email' and shows 'example@gmail.com' in a text field. Then it asks for 'Password' with a note '6-16 digits or letters or symbols' and shows 'Enter password' in a text field with a dropdown arrow. Below that, it asks for 'Re-enter password' and shows 'Enter password' in a text field with a dropdown arrow. At the bottom, there is a link 'I already have an account.' and a checkbox 'I have read and agree Privacy Policy'. A large teal button labeled 'NEXT' is at the very bottom.

3

Enter personal information, including the gender, nickname, birthday, height, and target weight. Select to enable the athlete mode if needed (for users who engage in long term physical work or exercise more than 3 times a week). Tap [Confirm] to save.



Note: Ensure that all the personal information is correct, since age, height, and gender will affect the accuracy of the measurement results.



4

Change the password.

If you forget your password, go back to the login page, enter your e-mail address and select [I Forgot my password]. Then Tap [Get Verification Code] and enter the verification code provided in the sent e-mail within 15 minutes.





If you want to change your current password, select [Account] --> [Settings] --> [Change Password], enter the old and new password, and tap [Confirm].

First Connection to the App

Note: When the screen of the scale is on, the Bluetooth of the scale is enabled. When the screen is off, its Bluetooth is immediately disabled, and the App cannot find the device.


- 1 Turn on the Bluetooth of the mobile phone and open the Fitdays App. For the Android system, you need to further enable the location.
- 2 Enter [Account] page --> tap [Device] --> tap + in the upper right corner (alternatively, enter the [Measure] page --> tap + in the upper right corner --> select [My Device] --> tap + in the upper right corner). Gently step on the scale to wake up the screen. When the App displays the searched device model, tap [Connect]. Now the device is successfully connected to the App.

Start the Measurement

- 1 Place the scale on a solid flat surface. Carpeted or uneven floors may affect the accuracy of the measurement results.
- 2 Calibration: Gently step on the scale to wake it up and wait until the scale displays "0.00", which indicates that the calibration is completed (Remember to re-calibrate the scale if it has been moved).
- 3 Connecting to the App: Turn on the Bluetooth of the mobile phone and open the App. Enter [Measure] page. Gently step on the scale to wake it up and the measure page will display "Connected" at the top, which indicates that the device is successfully connected to the App.

Note: Please take the first measurement with the App connected; and subsequent measurements can be taken offline without connecting to the App.

4 Correct posture for measurement

- ✓ Step on the scale with two bare feet. Ensure that the soles of the feet are in full contact with four metal electrodes and there is no contact between the left and right lower limbs (including thighs, calves and feet).


Your feet should not touch each other.

Your feet should be in full contact with the electrodes


- ✓ Hold the handle with both hands. Ensure that the fingers evenly cover the four electrodes on the handle with no contact between the hands.


Your fingers should be in full contact with the electrodes

Thumbs should not rest against the edge of the display screen.

- ✓ Keep your arms straight at a 30- to 45-degree angle.

Note: Please take off your shoes and socks, otherwise the scale can only measure your weight.


- 5 Keep the above measurement posture until all the data is displayed on the screen of the scale or the App prompts that the measurement is complete before getting off the scale.

Note: Do not step off the scale until the measurement is complete.

- 6 The scale will automatically shut down after displaying for about 15 seconds.

- 7 The scale is equipped with an automatic take-up reel. When you have finished measuring, please slightly lift up the handle and then gently put it back to the scale.



View Detailed Data of Body Composition

When the measurement is complete, the measured weight value, the Fat standard ratio and Muscle standard ratio of each part of the body will be displayed in [Measure] page.

Then, the screen displays the Weight body fat and the BMI changes compared with the last measurement (tap [Contrast] below to customize the historical data to be compared).


Swipe down to view the details. Tap [More] to view the measured values of all body composition items. You can also select a specific indicator for its description and analysis.

Tap  to view and share body composition analysis report.


View the Trend Chart

Enter [Charts] page to view the recent/weekly/monthly/yearly trend chart of 15 physical indicators.


View Historical Measurement Data

Enter [Charts] page --> tap [Click to view all data] --> select a specific date. All the measurement history of the selected date will be presented. Tap any measurement data to view the specific body composition value.

Delete Historical Measurement Data



Enter [Charts] page- [Click to view all data]- select the specific data to be deleted - slide left and tap [Delete]

Share and Export Data

Share measurement data details

- 1 Enter [Measure] page- tap - select [Save] or [Share].



- 2 Enter [Measure] page --> tap + in the upper right corner --> tap [Share] --> select [Save] or [Share].

3




Enter [Charts] page --> tap [Click to view all data] --> select the data --> tap  --> select [Save] or [Share]

Share the trend



1

Enter [Charts] page --> tap  --> tap [Custom management] to customize the indicator to be shared --> tap the share icon  --> select the way of sending the data.



Share the professional report

- 1 Enter [Measure] page --> tap  --> tap  --> select [Save to local]/[Print]/[Share]/[Cancel].


- 2 Enter [Charts] page --> tap [Click to view all data] --> select [User data] --> swipe down to the bottom --> tap [Report printing] --> tap  --> select [Save] or [Share].

Switch the Unit


When the device is connected to the App, choose [Account] --> [Settings] --> [Switch Units], and select your preferred unit.


Note: If the unit is changed in the App when the device is not connected, the device unit will not be synchronized until it is successfully connected to the App next time. (i.e. if the device is disconnected from the App, the device unit will remain unchanged despite any adjustments in the App). Therefore, it is recommended to set the unit in the App with the device connected.

Add, Switch or Delete the User

Add a user: Enter [Account] page --> tap [Add User] --> enter the user information --> tap [Confirm].


Switch the user: Enter [Measure] page --> tap the avatar --> select the user.

Note: Every new user needs to create his/her own user profile before measurement and ensure that the personal information is correct.

Delete the user: Enter [Account] page --> select the user --> slide left to delete.

Baby Mode

Enter [Measure] page --> confirm the user --> tap + in the upper right corner --> tap [Baby Mode] --> weigh the adult alone --> weigh the adult and the baby together --> View the results in the APP.

Offline Measurement

Prerequisites for offline measurement


- 1 The user has created correct personal information in the App and completed a measurement with the App connected, i.e. the user information is sent via the App and stored on the scale.
- 2 The scale can recognize the user, which requires the user's weight varies within ± 2 kg. If the measured weight change is more than ± 2 kg, connect to the App for measurement.

Offline measurement for multiple users

Conditions for the scale to distinguish users: The weight difference between users is greater than ± 2 kg, i.e. if the weight difference between two or more users is within ± 2 kg, the scale is not able to distinguish users. In this case, it is recommended to connect to the App for measurement.

Offline upload of data to the App


- 1 The scale can store the user data (24 sets).
- 2 The historical data will not be automatically uploaded to the App as soon as the device is connected to the App. The user needs to take a measurement with the App connected, and the historical data will be transmitted to the App together with the latest measurement data.
- 3 If there are multiple users, the scale will automatically assign the historical data to the users in the App according to the weight difference (± 2 kg) and send unidentified data to the data claim center. Tap the message icon in [Measure] page, and manually claim or delete the data.

Timestamp	Weight	Action
11:27 Mar 13, 2023	57.2 kg	claim
11:35 Mar 13, 2023	57.1 kg	claim
11:33 Mar 13, 2023	57.1 kg	claim
11:32 Mar 13, 2023	57.6 kg	claim
11:30 Mar 13, 2023	57.6 kg	claim
11:29 Mar 13, 2023	57.6 kg	claim


Other Settings

Select [Account] --> [Setting]. You can adjust the language or theme, Check the [FAQ], and leave a message for help in [Help].

Manual Recording

Enter [Measure] page --> tap + in the upper right corner --> tap [Manual recording] --> input the weight, circumference or height --> tap [Confirm].

Data Synchronization to Third-party APPS


Fitdays supports data synchronization across four Apps -- Apple Health, Google Fit, Fitbit and Samsung Health. The synced data items are as follows:

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Note: The feature is only for the primary user of the account.

Apple Health

Data synchronization

Click on "Apple Health"

Click on the headshot


Click on "Apps"

Choose "Fitdays"


Click on "Turn On All"

View data


1 In [Summary], you can view health data and edit which items are displayed.

2 In [Sharing], tap [Apps]. Select [Fitdays] to view the synced data of Fitdays.


3 In [Browse], you can view health data by health items.


Fitbit

Data synchronization

Follow steps 1-6 in the screenshot. The screen will display "Sync Successfully" when the setting is completed.

View data


The data in Fitdays will be synced to Fitbit when the setting is completed.

Note: If Fitbit cannot receive the data, check whether the Fitbit account to be synced is the same as the Fitbit account you actually view.


Google Fit

Data synchronization

Follow steps 1-6 in the screenshot. The screen will display "Paired Successfully" when the setting is completed.

View data

The measurement data in the Fitdays App will be synchronized to Google Fit when the setting is completed.

TROUBLESHOOTING

Battery is low.
Please plug in USB to charge.

The scale is overloaded.
More than 180kg is on the scale.

The sensor is incorrectly zeroed.
Place the device on a solid flat surface.

The sensor cannot be zeroed.
Replace the sensor.

Accuracy problems

Problem Description	Possible Causes	Solution
Inaccurate weight value.	The device is not placed on a flat floor or is used on a carpet.	Place the device on a solid flat surface.
	The device is not calibrated when it is used for the first time or after it is moved.	When the device is used for the first time or after it is moved, calibrate it before measurement. Place the device on a solid flat surface, gently step on it to wake it up and wait until the screen show "0.00" before measuring on the scale.
	Uneven standing position, too far to left/right/front/rear.	When measuring, keep the center of gravity stable in the center of the scale, and ensure that there is no contact between the left and right lower limbs (including thighs, calves and feet).
	The device is damaged. It does not weigh or the number jumps randomly after weighing.	Replace the device.
	Refer to the reasons of inaccurate weight value.	Refer to the solutions for inaccurate weigh value.

<p>Discontinuous measurement results / excessive deviation of multiple measurements.</p>	<p>The time of measurement, due to normal human metabolism, can cause a difference of about 1.5kg in a person's weight between morning and evening, and metabolism can affect fluctuations in water and electrical impedance</p>	<p>Take your measurements at the same time every day to minimize the deviation caused by natural laws. It is recommended to take measurements after getting up or before going to bed with less clothing.</p>
	<p>Large fluctuations in body data due to strenuous exercise or overeating.</p>	<p>Do not measure immediately after strenuous exercise or overeating.</p>
	<p>Changes of measurement posture.</p>	<p>Try to keep a correct posture for measurement. Step on the scale with two bare feet. Ensure that the soles of the feet are in full contact with four metal electrodes. Ensure that there is no contact between the left and right lower limbs (including thighs, calves and feet). Hold the handle with both hands. Ensure that the fingers evenly cover the four electrodes on the handle with no contact between the hands. Keep your arms straight at a 30-to-45-degree angle. Ensure that the arms are not too close to your body.</p>
<p>Inaccurate body composition data.</p>	<p>Incorrect measurement posture.</p>	<p>Ensure that the measurement posture is correct for each measurement. Step on the scale with two bare feet. Ensure that the soles of the feet are in full contact with four metal electrodes. Ensure that there is no contact between the left and right lower limbs (including thighs, calves and feet). Hold the handle with both hands. Hold the handle with both hands. Ensure that the fingers evenly cover the four electrodes on the handle with no contact between the hands. Keep your arms straight at a 30-to-45-degree angle. Ensure that the arms are not too close to your body.</p>

	Incorrect personal information (age/height/gender).	The calculation of the user's body composition data involves the user's age, gender and height. To ensure the accuracy of the measurement results, please enter the correct age, gender and height.
--	---	---

Measurement problems

Problem Description	Possible Causes	Solution
Weight value only, no other data	The user takes the measurement with shoes or socks on, or the product bag is not removed.	The user must measure with bare feet. Ensure that the soles of the feet are in full contact with four electrode plates/metal electrodes and there is no contact between the left and right lower limbs (including thighs, calves and feet).
	The user gets off the scale before the measurement is complete.	Keep the measurement posture until all the data is displayed on the scale screen or the App prompts that the measurement is complete before getting off the scale.
	No user data on the scale end.	Every new user needs to create a user profile in the App and enter the correct personal information. The first measurement needs to be taken with the device connected to the App and the user in [Measure] page is correct, so that the user data can be transferred to the scale and subsequent measurements can be done offline without connecting to the App.
	Incorrect personal information (height/age/gender).	Incorrect personal information usually refers to young age (<16 years old) or abnormal height, which are beyond the measurement range of the scale. In this case, the body fat will not be measured. Check whether the user's basic information is correct.
	The user is younger than 16 years old.	If the user is younger than 16 years old, only weight is measured.

The historical body composition data suddenly disappears.	The user's weight varies too much (greater than ± 2 kg) for the scale to recognize the user.	Conditions for distinguishing the user during offline measurement: The user's weight varies within ± 2 kg, otherwise the user needs to connect to the App for measurement (If the scale cannot recognize the user, the scale screen will only displays the weight value and "--" for other body composition items).
	The weight values of two or more users are so close (the difference is less than ± 2 kg) that the scale cannot distinguish.	Conditions for distinguishing multiple users on the scale during offline measurement: The weight difference between users is greater than ± 2 kg, i.e. if the weight difference between two or more users is within ± 2 kg, the scale is not able to distinguish users. In this case, it is recommended to connect to the App for measurement (If the scale cannot distinguish users, the scale screen will only display the weight value and "--" for other body composition items).

Connection problems

Problem Description	Possible Causes	Solution
Bluetooth connection failed.	Bluetooth is not enabled on the mobile phone.	For iPhones, manually turn on the Bluetooth. For Android phones, tap [Allow] when the request pops up.
	Mobile Bluetooth unauthorized.	For iPhones, tap [Bluetooth] in [Privacy and Security], and turn on the Bluetooth authorization for Fitdays.
	The GPS is not enabled on the Android phone.	Bluetooth 4.X provides the location function. Android 6.0 and iOS 10.0 and above versions can manage and authorize the Bluetooth 4.X device as a location device like GPS. When the App needs to search for a Bluetooth device, it requires the location permission. Since the body composition scale transmits data via Bluetooth, the GPS must be enabled in Android 6.0 and iOS 10.0 and above versions.

	<p>The App is disabled to use the network.</p>	<p>Please make sure you allow the App to use the mobile network or WiFi for later pairing and using. This is because a server connection is required to store the data in the cloud.</p>
	<p>Bluetooth is not enabled on the scale end.</p>	<p>When the screen of the scale is on, the Bluetooth of the scale is enabled. When the screen is off, its Bluetooth is immediately disabled, and the App cannot find the device. Please ensure that the screen is on throughout the connection.</p>
<p>Data transmission failed.</p>	<p>No online measurement is completed after the device is connected to the App.</p>	<p>The user needs to take a measurement with the App connected, and the historical data will be transmitted to the App together with the latest measurement data. (The historical data will not be uploaded until an online measurement is completed.)</p>
	<p>The device has been connected to multiple mobile phones, but the data is only transmitted to the currently connected mobile phone.</p>	<p>If multiple users use different mobile phones to connect to the same scale, the scale will upload all historical data to the currently connected mobile phone and then clear the history (For example, A and B use their own mobile phones to register the App and connect to the same scale. Assuming that there are 10 offline measurement records of A and B respectively on the scale, when the scale is connected to A's phone, all 20 offline measurement records will be uploaded to A's App, and B's phone will not receive these 20 records when it connects to the scale again).</p>

Unit problems

Problem Description	Possible Causes	Solution
The unit cannot be switched.	The user does not know how to switch the unit.	Switching the unit with the App: When the device is connected to the App, choose [Account] --> [Settings] --> [Switch Units], and select your preferred unit.
The unit on the device remains KG after setting to LB in the App. (or vice versa).	After setting the unit in the App, the user does not connect to the scale for synchronization.	If the unit is changed in the App when the device is not connected, the device unit will not be synchronized until it is successfully connected to the App next time. (i.e. if the device is disconnected from the App, the device unit will remain unchanged despite any adjustments in the App). Therefore, it is recommended to set the unit in the App with the device connected.

Power on problems

The scale cannot be powered on.	Power issues: The insulating sheet is not removed, or the rechargeable battery is not activated.	Cancel the shipping mode, press and hold the reset button and ensure that the battery is well-functioning.
	Pressure sensing problem: The loaded pressure fails to reach the minimum weighing range (>5kg).	The scale will only be turned on when the loaded pressure reaches the minimum range.

Charging problems

The scale cannot be charged.	Mismatched charger is used.	It is recommended to use a charger within a specification of 5 V, 1 A or 5 V and a current higher than 1 A, such as 5 V, 2 A. The charger is generally backward compatible. That is, the 5 V, 2 A charger will be compatible with the 5 V, 1 A charger. However, if the charger's input is not 5 V, for example, 10 V, this charger cannot be used.
	The provided charging cable is not used.	Cables with poor quality or cables that are too long can affect the charging or cause charging failure. It is recommended to use the provided charging cable.

FCC STATEMENT

FCC Warning:

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

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Prolinx GmbH
Prolinx Global LTD