DR-HO'S

circulation promoter

Reduce pain and become more active.

We hope you enjoy this total body pain relief system for many years to come!

User manual
P. 2

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Bedienungsanleitung
P. 50

Manuel d'instructions
P. 34

Manual de instrucciones
P. 67

Thank you

Thank you and congratulations on the purchase of your new DR-HO'S Circulation Promoter! You have made a smart and important decision with regard to self-treatment and care. Many people have been pleasantly surprised with the degree of improvement they have experienced with this state-of-the-art device, and you too will discover how it helps to improve your quality of life and reduce tension.

DR-HO'S Circulation Promoter is an elaborate but simple-to-use device for the whole family. From chronic to acute muscle tension and aches, this unique device provides quick and effective relief from many of your conditions. Its microchip program simulates deep pain therapy and muscle movement exercises for various common muscle strains and soft tissue sprains.

Our aim is to develop unique and effective products that help people improve their quality of life. I would like to take this opportunity to express our sincere thanks to you, our customer, for your support.

Please read the information in this manual carefully before using the device for the first time and keep it for future reference.

If you have any questions or require assistance with your DR-HO'S Circulation Promoter, please contact us via info@drhonow.eu or  .

Yours in health & happiness,

Dr. Michael Ho
Pain & Headache Specialist

TABLE OF CONTENTS

| | | | |
|--------------------------------|----|------------------------------------|----|
| Safety precautions..... | 3 | Cleaning & maintenance | 12 |
| Intended use..... | 8 | Problem solving..... | 13 |
| General description..... | 8 | Disposal..... | 15 |
| Specifications | 9 | Symbols..... | 15 |
| Treatment programs..... | 10 | Electromagnetic compatibility..... | 16 |
| Batteries | 11 | Warranty & customer service | 17 |
| Electrode pad guidelines | 12 | | |

SAFETY PRECAUTIONS

Warnings

The long-term effects of intense electrical stimulation are unknown.

PLEASE CONSULT YOUR DOCTOR BEFORE USING THIS DEVICE. THE DEVICE MAY BE DANGEROUS IF USED BY CERTAIN SUSCEPTIBLE INDIVIDUALS.

DO NOT USE THIS DEVICE

- If you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Using the device could then cause electric shock, burns, electrical interference, or death.
- Near magnetic resonance imaging (MRI) equipment. The device is MR Unsafe.

- If you have
 - abnormal high blood pressure;
 - an undiagnosed pain syndrome;
 - a tendency for internal bleeding;
 - if hemorrhage is likely.
- Following recent surgical procedures if stimulation may disrupt the healing process.
- Near your eyes, to your head or face.
- Together with medical electronic devices such as an artificial heart or lung or respirator.
- In the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), because it may cause them to malfunction.
- On open wounds, rashes, swollen, red, infected or inflamed areas or skin eruptions (e.g.,

SAFETY PRECAUTIONS

phlebitis, thrombophlebitis, varicose veins) or on top of or close to cancerous lesions.

DO NOT APPLY THE ELECTRODE BODY PADS

- Over insensitive areas of skin.
- On the opposite sides of your head, because the effects of brain stimulation are unknown.
- To your head, face, the front and sides of your neck, your throat, or on both sides of your thorax.
- To stimulate the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- To stimulate your front neck or mouth, because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Across your chest because the effect of electrical current on the heart may cause cardiac arrhythmias.
- On open wounds, rashes, swollen, red, infected or inflamed areas or skin eruptions or on top of or close to cancerous lesions.

- If you experience any skin irritation or redness after using, stop stimulation in that specific area until the skin returns to normal.

DO NOT USE THE CIRCULATION PROMOTER SYSTEM ON

- Children or infants, because the device has not been evaluated for pediatric use.
- People who have lost their mental capacity.
- Pregnant women, because it is not yet known whether it is safe to do so.

Seek advice from your doctor if you have suspected or diagnosed heart disease or epilepsy.

IN GENERAL: IF YOU FEEL YOU WOULD BENEFIT FROM USING DR-HO'S CIRCULATION PROMOTER AND WOULD LIKE TO GIVE IT A TRY, PLEASE CONSULT YOUR DOCTOR.

DO NOT USE THE DEVICE WITH THE FOOT BASE OR BODY PADS

- when in bath or shower;
- while sleeping;
- while driving, operating machinery, or any activity which might cause involuntary muscle contractions that may put the user at undue risk of injury.

SAFETY PRECAUTIONS

PRECAUTIONS - PAIN THERAPY DEVICE & FOOT BASE

- If you feel discomfort or abnormal symptoms during treatment, stop using the device and consult a doctor before resuming treatment.
- Do not use this device if you are allergic to plastic.
- If you have any problems setting up, maintaining or using the device, contact customer service and report anything unusual.
- Keep the device out of the reach of infants, children or pets. Any inhalation or swallowing of small parts is dangerous and can be potentially fatal.
- Only use parts and accessories that have been specified or authorized by the manufacturer. The use of unauthorized accessories may cause damage to the device or danger to the user.
- Do not use the device in areas of high humidity, such as the bathroom or in the shower, as this may cause damage.
- If using rechargeable 1.5 V batteries, make sure you use a CE approved battery charger. Never connect the device

directly to a battery charger or to other powered equipment.

- Do not use the device in the presence of a flammable anesthetic gas mixture, compressed oxygen or nitrous oxide.
- If the Pain Therapy Device is not functioning properly, you should stop using it.
- Only use the device for its intended purpose.
- Do not use the device while wearing other electronic devices, as this may cause damage.
- Do not use the Pain Therapy Device near a mobile phone, as this may cause it to malfunction.
- Do not throw the batteries into a fire because they may explode.
- Dispose of the device, its batteries and components in accordance with applicable legal regulations. Unlawful disposal may cause environmental pollution.

PRECAUTIONS - PAIN THERAPY DEVICE & ELECTRODE BODY PADS

- Apply pads to normal, healthy, dry, clean skin only.

SAFETY PRECAUTIONS

- If you wish to move the pads to another location of your body during treatment, turn the device off first, then remove the pads and resume treatment.
 - Ensure that metallic items with conductive capabilities, such as necklaces, belts and rings, are kept away from the electrodes during treatment.
 - Do not bend or fold the pads, as this may cause the pads to malfunction.
 - Always store the pads on the transparent protective film when not in use.
 - Do not apply any solvent or any other liquids (except conductive gel or water) to the pads or to your skin because these will cause the pads not to stick and to malfunction. The pads are gelled and will adhere to your skin when moist.
 - Always apply a few sprays of water (or conductive gel) onto the sticky surface of the electrodes whenever these feel dry, to ensure the pads will provide a comfortable, powerful and effective treatment.
 - The simultaneous connection of a patient to high frequency surgical equipment may result in burns where the electrode pads have been attached and possibly damage the Pain Therapy Device.
 - Do not leave pads attached to the skin after treatment.
 - Placement of the electrode pads and stimulation settings should be based on the instructions of the prescribing practitioner.
 - Do not share electrode pads with anyone else because this may cause skin irritation or infection. People should always have their own set of pads for personal use.
- IMPORTANT:**
- Both of the pads in a set must be placed on the body at the same time in order to work. Make sure the pads are not touching each other because otherwise no impulses will be felt.
- Stop using the pads if the individual packaging has been open for more than three months, or if the pads lose their adhesiveness. You could buy new sets of pads via your supplier.

SAFETY PRECAUTIONS

Please refer to the Quick Guide for user instructions and pad placement.

PRECAUTIONS - WIRES

- Do not bend or pull the end of the wires.
- Replace wires when broken or damaged.
- Long wires have the potential to strangulate. Be careful when using these in the neck and shoulder areas.
- Make sure the wires are properly plugged into the Pain Therapy Device and, when using the pads, that the pads are stuck tightly onto the skin area you wish to treat to ensure effective therapy.

NOTE:

Do not treat one area for extended periods of time (more than 40 minutes per session, up to 3 times per day) because your muscles in that region may become exhausted and sore.

Only use original electrode pads, gel and wires supplied by the manufacturer because other electrode pads and wires may not be electrically compatible with your device.

ADVERSE REACTIONS

You may experience:

1. skin irritation or hypersensitivity due to the electrical stimulation (or conductive gel);
2. skin irritation and burns underneath the electrode pads when applied to the skin.

YOU SHOULD STOP USING THE DEVICE AND CONSULT YOUR DOCTOR IF YOU EXPERIENCE ANY ADVERSE REACTIONS.

SERVICE

Do not attempt to repair or service this system, but contact Customer Service at info@drhonow.eu or your local supplier in the case you have problems or questions.

CUSTOMER CARE

If any of the instructions or warnings are not understood please contact Customer Service for assistance at info@drhonow.eu or contact your local supplier.

SAVE THESE INSTRUCTIONS!

INTENDED USE

DR-HO'S Circulation Promoter has been developed as a self-care household device for the whole family and is based on physics, modern microelectronics and clinical practices this system to provide head to toe pain therapy.

It can stimulate healthy muscles and increases local blood circulation, thereby improving and facilitating muscle performance and temporarily relieving pain associated with sore and aching muscles. It uses electrical current to stimulate nerves and relieve chronic, intractable pain and pain associated with arthritis relating to tension at the back of the neck.

IMPORTANT!

The Pain Therapy Device, model CP-1, is intended for the relief of pain associated with sore or aching muscles in the following anatomical locations: the upper and lower back,

back of the neck, upper extremities (shoulders and arms) and lower extremities (legs and feet) due to strain from exercise or activities of daily living principally regular household and work activities.

Note! This product is certified as medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, exclusively for the indication(s) of relief of pain associated with sore or aching muscles. Other non-medical uses escribed to this device are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

SPECIFICATIONS

Basic specifications

| | |
|----------------------------------|--|
| Power source | DC 3.0V, 2 AAA batteries |
| Method of line current isolation | Type BF applied part |
| Number of output channels | 1 |
| Number of output programs | 3 |
| Output intensity level | 5 steps |
| Timer range | 10, 20, 40 min |
| Weight | Pain Therapy System: 66g, Foot base: 1500g, Small body pads: 10g, Foot relieve pads: 200g, FlexTone pads: 37g |
| Dimensions | Pain Therapy System: 90 x 78 x 16 mm, Foot base: 39 cm Ø, Small body pads: 46 x 46 mm, Foot relieve pads: 260 x 110 mm, FlexTone pads: 80 x 120 mm |

Output specifications

| | |
|-------------------------------|---|
| Waveform | Pulsed, symmetric, biphasic, rectangular |
| Maximum output voltage | 40Vp @ 500Ω 80Vp @ 2kΩ 95Vp @ 10kΩ |
| Maximum output current | 80mA @ 500Ω 40mA @ 2kΩ 9.5mA @ 10kΩ |
| Pulse duration | 200μs |
| Pulse frequency | 13.7~48.5 Hz |
| Net charge (per pulse) | 0μC @ 500Ω (Method: Balanced waveform) |
| Maximum phase charge | 19.2μC @ 500Ω |
| Maximum average current | 1.53mA @ 500Ω |
| Maximum current density (rms) | 0.073mA/cm ² @ 500Ω |
| Maximum average power density | 0.056mW/cm ² @ 500Ω |

SPECIFICATIONS

| Additional specifications | |
|---------------------------|--|
| Operating environment | Temperature 5~40°C, Humidity ≤80%RH, Atmospheric pressure 86~106kPa |
| Storage environment | Temperature: Pain Therapy System -20 ~ 55°C, Electrode 5~40°C, Humidity ≤95% RH, Atmospheric pressure 50~106 kPa |
| Shelf-life | Electrode pads: 1 year, Foot relieve pads: 3 years, Pain Therapy System and Foot base: 3 - 10 years (depending on how well it is maintained) |

TREATMENT PROGRAMS

| Therapies: | A | B | C |
|---|--|--|---|
| | To temporarily increase local blood circulation and relaxation. | Symptomatic relief and management of chronic, intractable pain and of pain associated with arthritis. Stimulation of healthy muscles to improve and facilitate muscle performance. | |
| How long does the therapy last? | This depends on the timer setting chosen. The therapy ends automatically. | | |
| How do I personalize the therapy? | You set the program yourself and select the intensity (level 1 to level 5). | | |
| How long should I use the device and how often? | Each session should last 30 to 40 minutes. No more than 3 sessions per day. | | |
| When should I stop using it? | When the pain has subsided or stopped, if treatment is uncomfortable, in the event of an adverse reaction, and if the condition worsens. | | |

Note! Many people experience immediate relief while others require regular use to feel the benefits. The results vary and will depend on your underlying conditions and how often

you use the device. If your pain does not subside, you can try increasing the intensity or you can change the program.

BATTERIES

Insert 2 AAA batteries before use:

1. Remove the battery cover on the back.
2. Insert batteries correctly so that the markings correspond.
3. Refit the battery cover.

Important Notes

Low battery power will cause the device to lose its effectiveness. Please change to new batteries whenever stimulation weakens. Alkaline batteries are recommended as a stable power source.

If the indicator lights fade, or if the stimulation feels weak on your skin, this may be an indication of low battery power.

WARNING!

Replace only with the same or equivalent type of batteries, as recommended by the manufacturer: 2 x AAA (1.5 V) alkaline batteries.

Do not dispose of the batteries with normal household waste, but hand them in at an official collection point in accordance with local regulations governing the disposal of such materials.

Remove the batteries to prevent them from running down whenever you are planning to store the device for more than 3 months.

If the device is used once a day for about 30 minutes per treatment, one pair of AAA Alkaline batteries should last 3 to 6 months.

We recommend the use of Alkaline batteries because they provide a steady power level until it runs out. Non-Alkaline and rechargeable batteries have a tendency to decrease in power level over time, and you might mistakenly think the device is not working right because the treatment is not as strong.

Please refer to the Quick Guide for pad placement and user instructions.

ELECTRODE PAD GUIDELINES

- Wash skin thoroughly and then dry it before applying the electrode body pads.
- Apply the whole surface of the electrode pads firmly to the skin.
- Do not use electrode pads that do not stick properly to the skin, or only partially stick to the skin.
- If the skin looks red after treatment, you should not treat that area until the redness has disappeared.
- The electrode pads can be re-used 30 to 50 times depending on your skin type and other conditions.
- Whenever you notice that the pads are not sticking properly to your skin, or if the pads feel dry, you should moisten the pads with a

small amount of DR-HO'S Electro Therapy Conductive Gel or water using the spray bottle (do not wash them with water directly). This will help restore the adhesiveness and conductivity of the electrode pads. If moistening the pads fails to make them sticky again, they may need to be replaced.

To order new sets of electrode body pads and DR-HO'S Electro Therapy Conductive Gel, please visit the website of your local supplier.

CLEANING & MAINTENANCE

DR-HO'S Foot base and Pain Therapy Device:

- Use a wet cloth or neutral cleanser to wipe the surface.
- Do not use volatile liquids such as alcohol, thinner or gasoline.
- Keep well away from moisture, high temperatures, direct sunlight or water droplets.
- Please do not disassemble, repair or rebuild the device yourself because this will invalidate the warranty.
- Remove the batteries from the device if it is not used for more than 3 months.

Body and Foot relieve pads:

- Keep the electrode pads clean and never expose them to high temperature and direct sunshine.
- Keep the electrode pads at temperatures of 5~40°C.
- Gently clean the electrode pads with a moist cloth or use a brush to remove dust.
- Do not clean the pads with hot water or chemicals.
- Do not touch the gel of the pads with your fingers, as this could reduce their adhesiveness.

PROBLEM SOLVING

| Problem | Possible causes | Try this solution |
|--|--|--|
| The device does not appear to be working properly. | The two pads are not attached properly. | Attach both pads to the body and make sure the pads are not touching each other. |
| | Your bare feet are not making proper contact with the Foot pads. | Make sure you remove all socks and shoes. Your feet must be bare. Ensure both bare feet are on both Foot pads. |
| | The wires from the electrode pads or the Foot pads are not properly connected to the wires of the Pain Therapy Device. | Make sure the wire from the electrodes or Foot pads and main wire to the device are connected to each other. |
| | The wire plug is not firmly connected to the output port at the top of the device. | Push the wire plug firmly into the output port. |
| | Weakening electrical pulses. | Turn the ON/OFF intensity dial in a clockwise direction. |
| | Power is running out. | Replace batteries, even if you bought them recently. |
| | The electrode pads are not properly in contact with the skin. | Press the electrodes firmly onto the skin. |
| | The gel surface of the electrode pads is dry or contaminated. | Moisten the gel pads with a little water or replace the electrode pads. |
| | Foot pads and bare feet are not sprayed with water. | Spray water onto the soles of your feet and onto the Foot pads prior to use. |
| The device stops working before treatment has ended. | Your skin is covered with oil or lotion that is preventing the pads from sticking. | Wash your skin with water and try putting the pads on again. |
| | The batteries have run out. | Change the batteries. |
| | The wire plug has not been fully inserted into the output port. | Reconnect the plug to the output port. |
| | The wires between the electrodes or Foot pads and the main wire are not connected. | Reconnect the wires from the electrodes or Foot pads to the main wire leading to the device. |

| Problem | Possible causes | Try this solution |
|--|--|---|
| The skin has turned red and/or itches or prickles. | You are allergic to the electrodes. | Stop using the machine until you ascertain the cause of the skin irritation. |
| | You have a small cut or acne. | Do not put electrodes on recently shaved skin, areas that have suffered insect bites, scrapes, acne or any open wounds. |
| | The adhesive surface of the pads is dirty or too dry. | Wash the adhesive surface of the pads gently with your fingertips for about 3 seconds under slow running water. |
| | The therapy has gone on too long. | Use the device for a maximum of 40 minutes per treatment and no more than 3 times a day. |
| | The 2 pads are not attached firmly and fully to the skin. | Attach the 2 pads firmly to the skin to ensure full contact. |
| | The surface of the pads is worn out. | Replace the set of pads. |
| | Three yellow lights on the Pain Therapy Device are flashing and you cannot feel any stimulation. | One of the pads has become disconnected from your skin or a wire is unplugged. |
| Three green lights are flashing and you cannot feel any stimulation. | The auto shutdown timer has activated and the device has powered of. | Turn the power dial until it clicks to the OFF position first. Then turn the power dial back to the ON position. |
| Batteries have been installed but the light does not turn on. | The battery compartment is dirty and the connectors are blocked. | Examine and clean. |
| | The batteries are not at full power and have not been installed correctly. | Replace the batteries or install them correctly in the compartment. |
| | There is a problem with the connectors in the battery compartment. | Try pulling the connectors outward slightly using a small pair of pliers. |
| | The power dial has been turned OFF and then ON. | Turn the power dial OFF and then ON again. |

DISPOSAL

1. The Pain Therapy Device, Foot base and accessories should not be thrown away with household waste, but should be recycled by you or the supplier.
2. Dispose of packing materials appropriately in accordance with the rules and regulations in your area to prevent harm to the environment.

SYMBOLS

| Symbol | Meaning |
|---|---|
|  | SERIAL NUMBER. This symbol will be accompanied by the manufacturer's serial number. |
|  | CONSULT INSTRUCTIONS FOR USE. |
|  | Unique Device Identifier (UDI); the Application Identifier (01) represents the device identifier (DI) ("069707422XXXX"), the Application Identifier (11) the production date/date of manufacture ("YYMMDD", with "YY" the last two digits of the production year, "MM" the production month and "DD" the production day), the Application Identifier (21) the serial number of the device ("WWWWW-XXXXXX"). |
|  | Medical device. |
|  | Indication of conformity with the essential health and safety requirements set out in European Directives 0120 and 1639 – Notified body identification (SGS). |
|  | DATE OF MANUFACTURE. This symbol will be accompanied by the date of manufacture. |
|  | MANUFACTURER. This symbol will be accompanied by the name and the address of the manufacturer. |
|  | AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY. This symbol will be accompanied by the name and the address of the authorized representative in the European Community, adjacent to the symbol. |
|  | Dispose of separately from other household. |
|  | Type BF applied part (patient isolation from electrical shock). |

| Symbol | Meaning |
|---|--|
| IP22 | IP classification. |
|  | ATTENTION! Caution! Follow operating instructions! Failure to do so could place the patient or operator at risk. |
|  | Batch code. |
|  | MR UNSAFE. This device cannot be taken into the controlled access area of an MRI system. |

ELECTROMAGNETIC COMPATIBILITY

DR-HO'S Pain Therapy Device has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

CAUTION! Do not use this device simultaneously with devices that have high EMI levels.

“Reduce pain and become more active. We hope you enjoy this total body pain relief system for many years to come!”

Dr. Michael Ho

WARRANTY & CUSTOMER SERVICE

We hereby guarantee that if this product malfunctions or fails in any way within one year from the date of purchase, the purchaser will not be charged for any of the parts needed and labor costs incurred to repair the defective device, if it is used under normal conditions and in accordance with this instruction manual. This warranty is only valid if the purchaser can provide proof of purchase.

Please contact your local supplier in case you have questions.

Imported and distributed by
Dr Ho Europe B.V.
1261 XR - 5, the Netherlands
info@drhonow.eu | www.drhonow.eu

For customer service:
E: info@drhonow.eu
or your local supplier

Are you happy with your DR-HO'S Circulation Promoter®? Share your posts, tips & tricks! Like & follow us @drhonow & tag your posts

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#MyDrHo #DrHoLovers #DrHoNow

Pain relief from head to toe.

WARNING: DR-HO'S equipment is not intended to diagnose, treat or cure any disease or condition. If you experience pain or health issues, always consult your doctor for the right diagnosis and treatment.

PRECAUTIONARY MEASURES: Do NOT use DR-HO'S equipment or other electrostimulators if you wear a pacemaker or are pregnant. If you have any questions, please consult your doctor before using the device.

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Guang Zhou Xinbo Electronic Co Ltd.
MedNet EC-REP GmbH