




# Biohackn's Beacon

Users Manual





The **BIOHACKN BEACON** is one of our most potent portable products, with high output multi-wave LED technology. The Beacon provides expert treatment which comes in a tiny, mobile, cost-effective form.

It works by releasing a mix of Red Light (660nm) absorbed by your skin, leading to enhanced collagen formation and ideal skin health. Near-Infrared Light (850nm) penetrates deeper into your body, improves muscle recovery, and reduces joint discomfort. The beacon gives optimal adaptability for general health by emitting an equal quantity of red Light and Near-Infrared Light.

The Beacon is the world's best-targeted pain relief device. We designed it to work for common problematic areas such as the specific back, neck, shoulder, knee pain, and muscular tiredness or soreness.



## **How to use:**

1. Hold the beacon over a painful area about 1 inch from the skin.
2. We recommend not using the beacon for more than 5 minutes over a targeted area. Please turn off the device for at least a minute after operating it for 5 minutes.

## **Warning:**

1. Do not use the beacon around water or wet areas.
2. Do not stare at the Light.
3. Do not shine the Light near the eyes.

Thank you for choosing  BIOHACKN<sup>®</sup>

[www.Biohackn.com](http://www.Biohackn.com)