



# The Biohackn Mask

Users Manual





# The Biohackn Mask


Biohackns Anti-aging Mask uses our same anti-aging light technology to deliver Low-level light therapy (LLLT), delivered through best-in-class light-emitting diodes for optimal efficacy. In a flash, you'll have brighter, firmer, younger-looking skin. We designed our mask keeping simplicity in mind, so we chose 3 of the most effective light spectrums for a younger-looking you.



# **The world's most effective Anti-aging mask**

The Biohackn Mask delivers red, blue, and orange light sources using the biological phototherapy principle. Causing enhanced skin metabolism causes the skin to be transparent and luster effectively reduces acne, stimulates collagen production, improves cell microcirculation, and combines the essence of facial mask to achieve the effect of skin beautification.

# Product Features





## Usage Suggestion:

- After washing your face, wrap the strap over the head of the device and drape it over the face.
- Press and hold the power button for two seconds to switch on the device, then pick a light setting by pressing M (red, blue, red, or blue) and use
- After treatment, the device will shut down automatically.

## Methods of application:

1. Power on/off with a short press of the power key
2. There are four settings available: "automatic mode," "Red light," "Blue light," and "Orange light."
3. Red light: A lot of red light and graphene heat up, with some blue light thrown in for good measure
4. Blue light: Lots of blue light and graphene heat up, with some orange light thrown in for good measure

5. Orange light: Lots of orange light and graphene heat up, with some red light thrown in for good measure
6. Automatic Mode: Cycle through the three modes in succession. Every minute switches automatically.
7. Working Time: Each mode has a working time of 10 minutes. After 10 minutes, the machine shuts down, restarts, and returns to automated mode.
8. You may use the gadget while it is charging. When the device is in use, the charger is attached, which impacts the brightness of the light and the heating temperature.
9. When the battery power is low, the lighting brightness and heating temperature may be lowered; this is a typical occurrence that may be restored after the battery is fully charged.
10. Working Hours: 10 minutes per mode. After 10 minutes, the machine will shut down automatically and restart in "automatic mode."

## Maintenance and storage:

1. After use, wipe the mask with a dry cloth, do not wash it with water directly.
2. If this mask is damaged, please do not continue to use
3. Disassembling this mask is strictly prohibited.
4. Store in a cool, dry place, away from moisture and intense sunlight.
5. Please do not put the mask into the liquid.

## People Specific:

1. We do not recommend the mask to be used by pregnant women, and the Nano light source will not affect the baby, but it is still not recommended
2. People with epilepsy should not use it.
3. We do not recommend the mask be used with people with thyroid diseases and thyroid drug users
4. Not For people with light allergies

Thank you for choosing  BIOHACKN<sup>®</sup>

[www.Biohackn.com](http://www.Biohackn.com)