

USER MANUAL


Dura-Glass Personal Scale

Model : PS 115




Sensing Health Everywhere

Care and Maintenance



- Always keep the scale in horizontal position, and not put it in wet or hot condition.
- Avoid hitting or dropping the scale otherwise damage may occur.
- Do not put heavy things on the scale when not use.
- Do not overload otherwise the sensor may damage.



- Treat the scale with care as it is a highly precise electronic instrument.
- Do not immerse this product in water.
- Clean the scale with a soft damp cloth, do not use chemicals or harsh abrasives.



- Keep the scale surface dry and do not use it on slippery floor.
- Do not weigh with wet feet to avoid slipping.
- The product is not for trade use.

Warranty & Support

Entitle yourself for the 1 year HealthSense warranty by registering your purchase and easily reach us to submit your product related issues & queries, when you follow these hassle free ways.

At HealthSense our top priority has always been our valued customers. We very well know your HealthSense product is a part and parcel of your everyday life, which is why we are more easily accessible than earlier.

SCAN - For Warranty Registration



(Must register within 30 days from date of purchase)

SCAN - To File Complaint



(Fill & submit form if you are facing any product related issues)

* For detailed warranty terms please visit the link : www.healthsense.in/warranty-terms

* For additional support : Mail us at support@healthsense.in or Call us on : 080-41262836



www.healthsense.in



Healthsense is a proud Indian brand creating quality products idealised and conceptualised for India and manufactured by exclusive overseas partners.


Product Specification

| | |
|--------------------------|------------------------|
| Product size | 300x260x18mm |
| Power | 2 x 1.5V AAA batteries |
| Display | LCD readout display |
| Max Capacity | 180kg/396lb |
| Division | 100g/0.2lb |
| Minimum effective weight | 5kg |

Features

- Step-ON technology with Auto-Off feature
- Auto-Calibration (Reset to zero)
- Measurement units: kg & lb
- Temperature feature
- Low battery, Error & Overload Indicators

Product Structure



Battery Installation



1. On the base of the scale, gently remove the battery cover as per the arrow direction.



2. Insert 2 x AAA batteries into the compartment observing the polarity markings in the base.



3. After inserting the batteries, replace the battery cover and push the cover backward.

Kind notice:
Please be aware that the incorrect polarity will cause the permanent impact on the electronic parts.

How to Use



1. Place the scale on a solid, level surface.



2. Gently step on to the scale until it automatically turns ON. Once a steady reading is taken the display will lock the weight and flash to indicate the final reading



3. The scale will automatically switch off approximately in 8 seconds when no operation.

Kind notice:
1. Do not place the scale on carpet or similar soft surfaces otherwise accuracy will be compromised.
2. Please keep balance in weighing. Do not let the body wobble.

Unit Conversion

If the scale is with unit systems of Metric(kilograms) and Imperial (pounds), users can select units from kg(kilograms) to lb(pounds).



1. Open the battery cover.



2. Find the unit switch key in the battery compartment. Under weighing mode, press the unit switch key to select the desired unit.

Error Indication



The "Lo" symbol on the display indicates that the battery power is running low. Please replace new battery.



The "Err" symbol on the display indicates that the scale is over-loaded. Please step off to avoid possible damage to the product.



It indicates that scale has been calibrated. Please wait until it shows '0' and step back to weigh again

Battery Advice



- If the scale is not to be used for long periods, it's recommended to remove the battery to avoid damage to the scale due to possible battery leakage.
- Keep the batteries out of reach of children.



- Do not mix old and new batteries, with different compositions or of different brands in order to prevent possible leakage, explosion.
- Do not heat or deform the batteries or explore to fire.



- Waste batteries should not disposed of with household waste.
- Please check with your local authority for battery recycling advice.