

Warranty & Support >>>

Entitle yourself for the 1 year HealthSense warranty by registering your purchase and easily reach us to submit your product related issues & queries, when you follow these hassle free ways.

At HealthSense our top priority has always been our valued customers. We very well know your HealthSense product is a part and parcel of your everyday life, which is why we are more easily accessible than earlier.

SCAN - For Warranty Registration



(Must register within 30 days from date of purchase)

SCAN - To File Complaint



(Fill & submit form if you are facing any product related issues)

* For detailed warranty terms please visit the link : www.healthsense.in/warranty-terms
* For additional support : Mail us at support@healthsense.in or Call us on : 080-41262836

Caution >>>

- Do not use the device if you are pregnant or suffering from heart problems or any physical injuries, and consult your physician if you are unsure.
- Only use the device one hour after meals.
- Keep away from children.
- Do not use the device on the same area of body for more than 15 minutes at a time.
- Do not use the device in an environment that exceeds 37 degrees C or 100 degrees F.
- Do not use the device in bathroom or any humid places.
- If the device fails to work, please disconnect the same immediately and contact HealthSense support team.
- Do not use the device when you feel unwell.
- Confirm the device stop working when exchange the massage heads
- If the power wire is damaged, please ask the after-sales maintenance technician to change it.
- If the power cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.

Maintenance & Storage >>>

1. To clean, wipe the machine with a soft cloth and neutral detergent only.
2. Store it in dry place after use, avoid high temperatures, direct sunlight and corrosive gases.
3. Always unplug the unit if not in use



www.healthsense.in



Healthsense is a proud Indian brand creating quality products idealised and conceptualised for India and manufactured by exclusive overseas partners.

Connect with us @healthsenseindia

HEALTHsense®

USER MANUAL

Toner-Pro Handheld Body Massager

Model: HM 210



Sensing Health Everywhere



Introduction >>>

Dear Customer,
Thank you for choosing Toner-Pro. Our brand stands renowned for excellence in quality & craftsmanship, backed with well-grounded testing for applications related to weight, blood pressure, massage & fitness bands. You may please glance through the forthcoming instructions carefully & retain them for future use. Also ensure that this manual is easily accessible for fellow users' convenience.
Winning regards,
Your HealthSense Team

Toner-Pro helps relax your aching muscles while targeting those excess fat pockets to effectively reduce weight and improve stamina. Suitable for usage on waist, neck, shoulder, arms, hips, thighs, buttocks, calf and feet.



Parts and Functions >>>



Operation >>>

1. Connect the plug to a wall electrical outlet.
2. Switch power dial clockwise to turn "ON" the device, and continue in the same direction to speed up. To slow down, turn dial counter-clockwise.
3. After usage, switch the power dial counter-clockwise until it can not turn anymore to turn "OFF" the device, then remove plug from the electrical outlet.

Multiple massage heads with usage >>>

- Flat Massage Head** : Decompression & weight-loss massage for face and abdomen
- Wavy Massage Head** : Acupuncture point & deep tissue massage accelerates blood circulation
- Ball Massage Head** : Combine the essential oil of your choice for deeper absorption & toning of muscles
- Mesh Cover** : To avoid hair entanglement & added comfort during massage
- Micro Filer Massage Head** : Helpful in exfoliating dead skin & effectively remove callus from dry feet
(Notice : Do not use mesh cover when micro filer massage head is attached.)

Product Specification >>>

Model Name : Toner-Pro Handheld Body Massager
 Model No : HM 210
 Rated Voltage : AC 220V, ~50 Hz
 Rated Power : 28W
 Product Material : ABS
 Product Weight : N.W. 0.80kg, G.W. 1kg
 Product Dimension : L 226 x W 110 x H 135 mm

Directions for use >>>

1. **DISCONNECT FROM POWER SOURCE** : Holding the product in one hand, gently press the edge of the mesh cover and turn it anti-clockwise to pull it out.
2. **REMOVE MESH COVER** : Holding the product in one hand, gently press the edge of the mesh cover and turn it anti-clockwise to pull it out.
3. **REMOVE MESSAGE HEAD** : Find the two small gaps between the head and fuselage and gently pull the head upwards.
4. **PLACE DESIRED HEAD** : Replace the desired massage head by pressing it between the gaps on the fuselage
5. **REPLACE MESH COVER** : Holding the product in one hand, gently press the edge of the mesh cover and turn it clockwise until tightly seated.
6. **PLUG IN** : Reconnect the plug to the power source to start using the massager.