




## BATTERY & CHARGING



## REPLACING BRUSH HEAD



## STAY CONNECTED



Access your specialty SPM-1 enamel protector brush heads at [getevo.net/subscription](https://www.getevo.net/subscription)

**In case of any issues that you may have...** our customer support team would love to help. Reach out by emailing [support@getevo.net](mailto:support@getevo.net).

Life happens. Your toothbrush is covered with a 1 year warranty\*. Register for a subscription plan and receive a lifetime warranty\*.

\*See Evo-Instruction Manual for warranty details.

# SPM-1

## RECHARGEABLE SONIC TOOTHBRUSH

COLLECTOR'S EDITION

evo.  
GETEVO.NET

evo.



RECHARGEABLE SONIC TOOTHBRUSH  
USER GUIDE - COLLECTOR'S EDITION

## ABOUT EVO

### EVO® stands for evolution.

Evo was born from the dream of creating futuristic products today. Conceptual products, such as concept cars, are inspiring, amazing to look at and usually something we are unable to purchase for ourselves. So we simply asked "Why?"

This was the start of our journey. Working with some of the most creative and intelligent designers and engineers around the world, we examine every aspect of our products to make them a reality. We are dreamers, creators, and designers infusing innovative technology with great design to produce "dream products" that one may have never imagined possible.


### INTRODUCING EVO-SENSE TECHNOLOGY

Change the way you control your toothbrush by using EVO-SENSE. This feature allows you to control the intensity of the vibration by simply swiping up or down with your finger.

### Why Marvel?

Marvel has always inspired us with the thought of futuristic technology that could be possible in our world. Looking to the future, we partnered with Marvel to create SPM-I as the next part of our collaboration.

## INTRODUCING...




## EVO-SENSE TECHNOLOGY


**Let us introduce to you our latest feature: EVO-SENSE**  
EVO-SENSE allows you to control the intensity level of the vibration by swiping up and/or down with your finger.

**How does EVO-SENSE work?**  
When selecting a mode, you have 8 seconds to select the preferred intensity level for the current mode. Swiping up increases intensity. Swiping down decreases intensity. We created a locking system on EVO-SENSE, that automatically activates 8 seconds after selecting a mode. 8 seconds after selecting a mode, EVO-SENSE will lock. Reset if desired intensity was not chosen. To reset: turn the brush off then back on. You will have 8 seconds to select the desired level of intensity. We created this locking mechanism to avoid accidentally changing the intensity while brushing.

For more information go to [getevo.net/subscription](http://getevo.net/subscription)



## GETTING STARTED



## HOW TO USE

