Tesmed+  *Cellulite*

Patent pending

Manuale d'uso

ITALIANO

User Manual

ENGLISH

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DEUTSCH

Mode d'emploi

FRANÇAIS

Manual de instrucciones

ESPAÑOL

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PORTUGUÊS

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POLSKIE

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NEDERLANDS

Instruktionsbok

SVENSKA

TESMED CELLULITE

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TESMED CELLULITE

*Dear Customer, thank you for choosing Tesmed Cellulite.
Please read the following information carefully.*

Clinically proven to be effective.

Tesmed Cellulite, patent pending, is the only massager against cellulite blemishes that treats the entire thigh. It takes just a few minutes per day to obtain a rapid and significant improvement of micro-circulation, extensibility of the skin and overall elasticity, as well as an evident reduction in the appearance of cellulite. An effective reduction in the circumference of the thigh** and a notable increase in smoothness can be seen after just two weeks.

After just 4 weeks, 87% of users reported an improvement in circulation, 70% of users efficiency in improving firmness, and 67% of users efficiency in reducing the appearance of cellulite.*

* instrumental and subjective assessment of efficiency on 30 volunteers with cellulite on the thighs and buttocks who used the product on one thigh for 5/10 minutes per day for a period of 4 weeks - instrumental assessment after 2 and 4 weeks.

Ref. ISPE: 49/21/01 – 66/21

** -1.7 cm. The best result encountered among 30 volunteers

CONTRAINDICATIONS:

Massage with TESMED CELLULITE is not indicated for:

- patients with malignant tumour pathologies
- people with current cancer or those having had cancer less than 5 years ago
- patients undergoing anticoagulant therapies
- pregnant women
- patients suffering from major cardiovascular diseases
- patients with haematomas, phlebitis, angiomas or varicose veins.

This massager must be intended only for the use for which it was expressly designed. Any other use is to be considered improper and therefore dangerous. The manufacturer cannot be held responsible for any damage resulting from improper, erroneous and unreasonable use.

Massages using the cellulite massager are not a substitute for treatments administered by doctors .

PRECAUTIONS:

- A - Tesmed Cellulite is not recommended in the presence of fragile capillaries.
- B - On sensitive skin, easy bruising skin, and or in the presence of fragile capillaries, Tesmed Cellulite may cause bruising, and therefore should be used with extreme caution, avoiding pressure on the skin.
- C - It is not necessary to press the rollers too hard on the skin. Do not press too hard on the skin to avoid bruising, redness and peeling, which may appear at the beginning of the treatments and then disappear in a few days.

- D - Temporarily suspend treatments in case of widespread bruising. The hematomas can generally appear at the beginning of the treatments, they disappear spontaneously within a few days.
- E - Light reddening of the skin during and after the massage may indicate that your skin microcirculation has been stimulated. The reddening will gradually disappear.
- F - In case of significant and / or persistent complaints, stop the treatments.
- G - In order to allow the rollers to perform optimally, carry out the treatment on perfectly clean and dry skin and do not use any creams, oils, etc.
In the case of skin that is easy to bruise or discretion, it may be advisable to use in conjunction with a body or anti-cellulite cream, so the wringing action will be weakened. In that case, first spread the cream on the area to be treated and then use the device as directed.
- H - For best results, follow a balanced diet and exercise regularly.

Cellulite

Cellulite affects 80% of women.

It manifests as visible blemishes, the cause of which may be a deficiency of microcirculation, an alteration of the subcutaneous tissue, which is rich in adipose cells that can increase in volume after puberty, with a consequent accumulation of excess liquids causing the so-called orange peel skin effect.

Blemishes can be more or less evident and affect certain areas more than others, such as the thighs, buttocks and abdomen.

TESMED CELLULITE TECHNOLOGY

International patent pending

Tesmed Cellulite is based on the best professional massage technology.

The device has multiple rollers arranged in **converging pairs** followed immediately by **diverging pairs**, set out in such a way that as they slide over the skin in a certain direction, **the converging rollers lift and tighten the part being treated with a “pinching” and “kneading” effect**, while the pairs of **divergent rollers simultaneously stretch the skin**.

When moving in the opposite direction, the action of the rollers is reversed.

The repeated action of **the converging rollers** provides significant benefits in combating cellulite blemishes as it promotes the expulsion of excess liquids while applying a potent toning effect.

The repeated action of **the diverging rollers** stretches the skin, consequently providing a smoothing action.

The particular shape of Tesmed Cellulite combined with the quantity of its rollers allows between 80 and 95% of the area requiring treatment to be covered in just one pass, meaning that effective treatment takes just a few minutes.

Tesmed Cellulite provides a kneading massage effect, while at the same time stretching the skin. Its articulated mechanical action can stimulate skin microcirculation and promotes the dispersion of excess liquids.

The shape of Tesmed Cellulite has been designed to suit the majority of thighs and buttocks.

BENEFITS

Using Tesmed Cellulite can help you obtain softer, more elastic skin, helping to significantly reduce the visibility of blemishes.

The treated area can appear more tonic and smooth.

TREATMENTS

THIGHS:

Use for 3 or 5 minutes, 1 or 2 times a day.

Place your foot on a stool and slide the Tesmed Cellulite, exerting light pressure (such as not to cause bruising or hematomas), repeatedly along the length of the thigh, **fig. 1 - pag. 3**, in both directions.

While using Tesmed Cellulite, rotate it a little (between passes) to cover the entire thigh area if necessary.

For skin that is easy to bruise or discretionary, it may be advisable, prior to use, to spread a body or anti-cellulite cream on the area to be treated, so the wringing action will be weakened.

BUTTOCKS:

Use for 3 or 5 minutes, 1 or 2 times a day.

In a prone position, lying down, slide Tesmed Cellulite over the buttocks repeatedly in both directions, trying to keep it perpendicular to the buttocks, **fig. 2 - pag. 3**, exerting light pressure (such as not to cause bruising or hematomas).

If it is easier, treatment on the buttocks can also be performed with the assistance of a second person to apply the massage with Tesmed Cellulite.

MAINTENANCE AND CLEANING

The device is maintenance free, use soap and water for washing. Do not use solvents.

HOW THE WARRANTY OF CONFORMITY WILL BE EXERCISED

The consumer will be able to contact the point of sale where he purchased the product. The point of sale where the purchase was made may not refuse to take delivery of the non-conforming product for repair or replacement.

If the purchase is made via websites or distance selling, the user will have to keep the shipping confirmation received or any documents that can prove the order placed and received and contact directly to the seller for repair or replacement under warranty .

For warranty conditions, please visit the **tesmed.com** site in the section entitled Legal Notices, Sales Terms, and Warranty.

For assistance write to **info@tesmed.com**

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Feldi S.r.l.