MOMAX®

Health Tracker

IoT Body Scale

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(EW1S)

Instruction Manual

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
N

Introduction p.2 - 27

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Introduction

Health Tracker IoT Body Scale (EW1S)

- In addition to weight and body BMI calculation, you can check 17 different body data analysis and automatically track daily data trends.
- With original sense to patent technology, EW1S will automatically switch on as you step on the platform barefooted. When Wi-Fi is connecting, the analysis results will automatically identify the User ID which you had set. The results are also sent to the phone's Momax Smart APP.
- Please open the App and turn on your Wi-Fi, so that the EW1S can successfully connect to your smartphone. Keep your smartphone and EW1S in transmission range. Step on the EW1S and take a measurement. Once the measurement is finished the data is automatically transferred to the app.
- Applies BIA (Bio-impedance Analysis) technology. A small amount of weak current flows through the human body so as to detect the Bio-impedance and estimate BMI, body fat, body water, muscle mass, bone mass and BMR. The electrical current is small and may not be felt. This BIA technology is cheap, safe, non-invasive, toxic-free and harmless. It also possesses the characteristics of simple operation and abundant information.
- The current mentioned above is less than 0.5mA. However, please be aware that anyone with any wearable or implantable medical electronic instrument, such as a pacemaker, must avoid using this device. The intended use of this device is for healthy children 10-17 years old and healthy adults.

Before Using

To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the EW1S on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Step onto the platform with bare feet. Stand still and keep full contact with the electrodes until the measurement is complete.
- Start measurement at least two hours after getting up or dining.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dining.
- Always start measurement at the same time slot and on the same scale located on the same flat, hard surface.
- For maximum accuracy and repeatability, it is recommended that you should use the device at the same time of the day and on the same location.
- The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take the reading with clean, slightly damp feet for best accuracy. If you are having a problem on operating this scale please contact customer service.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to factors such as alcohol consumption, menstruation, illness, intense exercise, etc.
- Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.

Quick Installation & user guide

1. Open the battery cover on the back of the EW1S and insert the battery (3x 1.5V AAA battery), pay attention to the polarity symbols (+ and-) in the battery compartment. Close the battery cover.
2. Place the EW1S on a flat, hard surface. Carpets or uneven floors may affect accuracy.
3. Wait until the "0.0kg" number is displayed on the LCD.
4. Choose the unit of measurement in kilograms, pounds and stones.
5. Step on the body fat scale and stay still, the EW1S will start to work.
6. EW1S will show your weight value.
7. After a few seconds, the EW1S will automatically turn off.

Quick Start Guide

1. Download 'Momax Smart' APP

Search 'Momax Smart' app on iOS or Android or any platform.

MOMAX smart


2. Register

Register an account and log in using your mobile number or email address.

 | The Home page of all connected devices.

 | To set scenarios or conditions and respective actions of the new smart device.

 | Profile management, FAQ, setting, etc.


3. Turn On the Device Wi-Fi

The app is ready to add a device, turn on Wi-Fi on the EW1S.

Preparing mode:

Press and hold ON/OFF for 3 secs until LCD display "  " flickers and tap "Confirm light blinks rapidly".



4. Connect to the Internet Wi-Fi

Choose your Wi-Fi network and enter the password. The smart device will start connecting to your device via Wi-Fi network.


Tips:

- Stable Wi-Fi network is required
- Link to the Wi-Fi network again when the device is disconnected

5. Successfully connected

Once the connection is completed. You can edit the name of the EW1S. Make sure your Wi-Fi is connecting, the results will send to the phone's APP.

6. Input personal information

- Enter personal information including sex, birthday, height
- Up to 10 users can be added
- * Please make sure all blanks have been filled

<cancel Add personal information Save	
Head portrait	 >
Full name	MOMAX Smart >
Sex	Female >
Birthday	1994-12-28 >
Height	160 >
<input type="button" value="delete user"/>	


7. Place the EW1S on a flat, hard surface. Carpets or uneven floors may affect accuracy.
8. Wait until the "0.0kg" number is displayed on the LCD.
9. Choose the unit of measurement in kilograms, pounds and stones.
10. Step on the EW1S and stay still, the EW1S will start to work.
11. EW1S will show your weight value.
12. Sync data from EW1S to Smart App.

	Successful Wi-Fi Connection
	Data ready to transmit to the App [10 sec approx.] <ul style="list-style-type: none">• If SUCCEED, the symbol disappears• If FAIL, the symbol remains

13. View and claim data

1. Open SMART MOMAX APP and click into the device.
2. Click "View" as shown in the picture below.
3. Find the data that synced to the APP and click "claim". The data will be recorded into your user account.


The screenshot shows a 'Data claim' table with two rows of data. Each row includes a name, a date, a weight value, a unit, a BMI value, and a 'Claim' button.



< Back		Data claim				
Jong	04/16	88,9	weight(kg)	30	BMI	Claim
Lam	04/14	45,6	weight(kg)	18	BMI	Claim

Smart Home Ecology System

Adding Multiple Devices in interactive smart mode.

Product Overview

LCD Display


	Wi-Fi Connection
	Data pending to transmit to smartphone
	Low Battery
st	Stone
lb	Pound
kg	Kilogram

Maintenance

When carrying out usual maintenance, please ensure the practice of the following Dos and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use a propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.

Specification




Model No.	: EW1S
Power supply	: 4.5V (AAA batteries X 3)
Water resistant level	: IP21
Unit display	: kg/lb/st
Weighing range	: 5 - 180 kg [max. weighing 183kg] 11 - 396 lbs [400lbs max.]
Accuracy standard	: 50kg (±0.3kg) 100kg (±0.4kg) 150kg (±0.5kg) 180kg (±0.7kg)
Wireless frequency	: 2412MHz - 2472MHz
Wireless distance	: 10m
Size	: 300mm X 300mm X 28mm
Weigh	: 1.8kg [not including batteries]

Contains

1. Health Tracker IoT Body Scale
2. AAA batteries (3 pcs)
3. Instruction Manual And Warranty Card

Trouble Shooting

Problem	Root Cause	Solution
No display on LCD when the device powers on.	Batteries not yet installed.	Install the batteries. (Please refer to Insert the Batteries)
	Worn batteries.	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.
CANNOT proceed to analyze BMI, body fat, total body water, muscle mass and bone mass.	Step onto the platform wearing socks or shoes.	Please keep barefooted during the measurement, and keep full contact with the electrodes as well.
	The user fails to select the User ID from what the system found.	Please assign a User ID.
The device powers off automatically.	Low battery	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.

LCD Display	Description	Solution
	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off in four seconds.	Replace all three batteries at the same time. Please purchase the authorized batteries for replacement.
	Wi-Fi connection fail.	Rejoin Wi-Fi network.

Appendix - Body Composition Data

1. Health Tips - About Body Fat

Human body is made up of a percentage of fat. Body fat is vital for a healthy, functioning body, protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat will damage your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror.

This is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measure of fitness than weight alone—the composition of your weight loss could mean you are losing muscle mass rather than fat—you could still have a high percentage of fat even when a scale indicates 'normal weight'.

The body fat percentage (%) : 5%-60%/0.1%

Table of body fat rating					
Standard for Men	Age				
Rating	20-29	30-39	40-49	50-59	60+
Low	<13.0	<14.0	<16.0	<17.0	<18.0
Normal	13.0-20.0	14.0-21.0	16.0-23.0	17.0-24.0	18.0-25.0
Moderately High	20.1-23.0	21.1-24.0	23.1-26.0	24.1-27.0	25.1-28.0
High	>23.0	>24.0	>26.0	>27.0	>28.0
Standard for Women	Age				
Rating	20-29	30-39	40-49	50-59	60+
Low	<19.0	<20.0	<21.0	<22.0	<23.0
Normal	19.0-28.0	20.0-29.0	21.0-30.0	22.0-31.0	23.0-32.0
Moderately High	28.1-31.0	29.1-32.0	30.1-33.0	31.1-34.0	32.1-35.0
High	>31.0	>32.0	>33.0	>34.0	>35.0

2. Health Tips - About Body Water

Body water is the single most important component of body weight. It represents over half of your total weight and almost two third of your lean body mass (predominantly muscle). Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body's temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be 'normal' but your body hydration level could be insufficient for healthy living.

The body water percentage (%) : 43% -73% /0.1%

Table of body water level (Unit: %)		
Rating	BF % RANGE	OPTIMAL TBW % RANGE
Men	4% to 14%	70% to 63%
	15% to 21%	63% to 57%
	22% to 24%	57% to 55%
	25% and over	55% to 37%
Women	4% to 20%	70% to 58%
	21% to 29%	58% to 52%
	30% to 32%	52% to 49%
	33% and over	49% to 37%

3. Health Tips - About Muscle Mass

According to the American College of Sports Medicine (ACSM), lean muscle mass may decrease by nearly 50 percent between the age of 20 and 90. If you do nothing with it you're losing muscle and increasing fat. It is also important to know your muscle mass % during weight reduction. At rest, the body burns approximately 110 additional calories for each kilo of muscle gained. Benefits of gaining muscle mass include:

- Reversing the decline in strength, bone density and muscle mass with aging
- Maintenance of flexible joints
- Guide weight reduction when combined with a healthy diet.

Muscle Mass Percentage	
Men	Approximately 40% of total body weight
Women	Approximately 30% of total body weight

4. Health Tips - About Bone Mass

Regular exercise and a balanced diet can help maintain healthy bones. Like muscle, bone is a living tissue that can respond to exercise by becoming stronger. For most people, bone mass peaks in their thirties. Then people begin to lose bone mass. Regular exercise can help prevent that loss.

Calcium and vitamin D, good sources of which are dairy products, green leafy vegetables and fish, contribute to healthy bones.

The bone mass readings given by this product are an estimation of the amount of bone in your body. Those with osteoporosis or low bone densities may not get accurate estimations. If you have any concern regarding your bones please consult your doctor timely.

Bone Mass Percentage

- The average bone mass percentage for both men and women is between 4 to 5%.

Table of bone mass between gender and age					
Sex	Age	Low (-)	Normal (0)	High (+)	Very High (++)
Men	20-39	<8.0	8.0 - 19.9	20.0 - 24.9	≥ 25.0
	40-59	<11.0	11.0 - 21.9	22.0 - 27.9	≥ 28.0
	60-79	<13.0	13.0 - 24.9	25.0 - 29.9	≥ 30.0
Female	20-39	<21.0	21.0 - 32.9	33.0 - 38.9	≥ 39.0
	40-59	<23.0	23.0 - 33.9	34.0 - 39.9	≥ 40.0
	60-79	<24.0	24.0 - 35.9	36.0 - 41.9	≥ 42.0

Source: NIH/WHO guidelines for BMI
Source: Gallagher et al., American Journal of Clinical Nutrition, Vol. 72, Sept. 2000

Warnings

Before using, please read this user manual carefully and keep it for future reference.

5. Health Tips - About BMI

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. The Analyzer will calculate your BMI on your height and weight.

6. Health Tips - About Calorie

The device uses your age, height, gender and an activity level setting to calculate the estimated number of calories you can consume a day to maintain your present weight. This estimation may gradually help you with a weight loss, gain or maintenance plan.

EMC Guidance

1. MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment.

(Note. As indicated in Table 6 of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m)

1. MOMAX Smart Body Scale EW1S offers you a seamless way to manage your health.
2. Please note that this device is designed for self-measurement and self-monitoring of physical health. Any information provided by this device is in no way meant to treat, cure or prevent any disease. If in doubt, please contact your doctor.
3. This device is prohibited to be used by any female user who is suspected pregnant or being pregnant. Otherwise, the effect of the device on the fetus is unknown.
4. This device is prohibited to be used by any user who is wearing or implanted with an electronic device or instrument (such as a pacemaker or defibrillator).
5. This device should not be used by anyone who is acutely or chronically ill because of suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.
6. WARNING: No modifications of this equipment are allowed. This may result in increased EMISSIONS or decreased immunity of scale.
7. Please use the device according to the user manual. Any misuse can cause electric shock, burns, fire and other unexpected hazards.
8. Please use and store the device under the environment which was provided in the user manual.
9. Don't expose the device to extreme temperatures, direct sunlight, moist or corrosive environment.
10. Do not step on the scale when your body or feet are wet, especially after bathing or showering to prevent slipping.

11. Please keep the device out of reach of infants, children or pets, since inhalation or swallowing of small parts is dangerous or even fatal.
12. Please do not use this device, for example, due to the performance of the sensor and the electrode is degraded or the electrode is loose, it will reduce the performance or cause other problems.
13. If you have an allergy to stainless steel, please avoid contacting the electrodes of the device.
14. Please do maintenance according to the instructions in the user manual.

Service Agreement of Smart Mobile Application Platform

Please reference Mobile Apps content: Before accepting the Agreement, you are supposed to thoroughly read all the contents of the Agreement, and fully understand its terms, especially restrictive clauses or exceptions and Privacy Notice. Restrictive clauses or exceptions are bold or highlighted in other forms to catch your attention. In case of any doubt about the terms of the Agreement, please contact the relevant business department of Tuya. You are not allowed to use the service before you have read and accepted all the terms of the Agreement, relevant agreements, and rules, etc. Once you select "agree and submit the Agreement" (see registration page for detailed wordings) and complete the registration procedure, or you use the service in any form, it will be deemed that you have read and agreed with the restriction of the Agreement and rules above. In case of any breach of the Agreement, Tuya has the right to unilaterally restrict, suspend or terminate the service for you, and has the right to investigate your relevant responsibilities.

Disposal of Electronic Products

1. Do not place the product as household waste;
2. Please dispose of discarded electronic products properly according to local regulations;
3. MOMAX encourages you to participate in the electronic recycling program.

Warranty Conditions

1. For any defect, in the judgment of our technician, caused under normal use, we are responsible for repairing free of charge 12 months commencing from the date of purchase at the discretion of our company;
2. A valid warranty card and original purchase invoice must be presented upon request of warranty service;
3. This warranty card shall not apply to the product if:
 - (a) It is consumable part(s), such as Filter, or;
 - (b) it has explicitly modified or repaired in any way by persons other than technicians of this company, or;
 - (c) it has been damaged through misuse, malicious acts, negligence, the inflow of any liquid, or;
 - (d) losses or damages are caused by force majeure or natural disasters.
4. No alteration to this warranty card will be allowed;
5. The warranty terms and conditions are subject to change without prior notice;
6. Please retain your proof of purchase for warranty service;
7. Please consult the dealer about the product compatibility with your cell phone/ PDA before purchase. Product specifications are subject to change without prior notice;
8. Please visit www.momax.net for the latest product information.

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Copyright © 2020 MOMAX Technology Ltd. All rights reserved.
MOMAX is registered trademark of MOMAX Technology Ltd.
All other trademarks are the property of their respective owners.