**MOMAX** smart


## Lite Tracker

IoT Body Scale

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(EW2S)


Instruction Manual

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Introduction p.2 - 26

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## Introduction

### How does Lite Tracker (EW2S) work?

This model applies BIA (Bioelectrical Impedance Analysis) technology, which makes measurement based on the study of conductive composition of muscle and water contents within the body. A small amount of weak current flow through the human body so as to detect the bio-impedance and estimate BMI, body fat, body water, muscle mass, bone mass, etc. The BIA technology is safe, non-invasive, toxic-free and harmless. It is easy to operate and provide comprehensive information.

### 13 Body Analysis

- Body weight
- Body fat
- Body water
- Visceral fat index (VFI)
- Basic Metabolic Rate (BMR)
- Body age
- Body scoring
- Body mass index (BMI)
- Muscle
- Fat free mass (FFM)
- Bone volume
- Protein
- Body type

### Before Using

To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- This scale not be used to diagnose or treat any medical condition. You should consult your primary care doctor or physician before undertaking your diet or exercise program.
- Not recommended for any person using a pacemaker or other internal medical device.
- Not recommended for women who are pregnant, the measurements will be incorrect.
- Place the EW2S on a flat, hard surface.
- Step onto the scale with bare feet. Do not jump onto the scale or stand on the edge. Stand still and keep full contact with the electrodes until the measurement is complete.
- Start measurement at least 2 hours after getting up or dining.
- Avoid measurement immediately after strenuous workouts, sauna or bathing, drinking or dining.
- Take measurement at the same or similar time of the day.
- Do not overload the scale (max. 180 kg / 396 lb).
- Install the batteries according to the correct polarities indicated.

## Quick Start Guide

1. Open the battery cover on the back of the EW2S and insert the batteries (3 x AAA battery) according to the polarity indicated in the battery compartment. Close the battery cover.
2. After battery installation/replacement and before measuring body weight, place the scale on the hard floor and the scale will show its weight. Wait for about 10 secs until it is OFF, step on the scale with the bared foot for 3 secs. The scale display will start up from 0.0 then show your body weight.

### 3. Download 'Momax Smart' APP

Download 'MOMAX smart' from iOS App Store, Android Google Play Store or other platform. Enable Bluetooth function in your device.

### 4. Register

Register an account and log in using your mobile number or email address.



The Home page of all connected devices.

To set scenarios or conditions and respective actions of the new smart device.

Profile management, FAQ, setting, etc.


## 5. Bluetooth Pairing

Click "+", then tap "MOMAX Smart Lite Tracker" to add device.

Step on the scale barefoot to start pairing the EW2S to the App. You will see a Bluetooth icon flashes quickly on the scale display. Tap "Confirm" the indicator is blinking rapidly. Then, tap 'Next' to proceed.

\* This product cannot be paired directly through the Bluetooth list of the mobile phone.

\* The pairing process may take 2 minutes, please do not perform other operation.


## 6. Connection successful

If successful you will be prompted with a scale icon. You can edit the name of EW2S. Please do not disable Bluetooth on your device because EW2S will sync your data to 'MOMAX Smart' app.

## 7. User Information


Fill up all user information such as Name, Gender, DOB, Height and Target Weight. Information such as Gender, Age and Height is crucial for the App to calculate your body composition.

\* Please make sure you filled in all information.


## 8. Add Members

You can add up to 20 members.


## Regular Measurements

After pairing is done, you can perform measurement anytime from the home page of the scale. Tap the profile pic located at the upper right side of the page, select a member, and then step on the scale. Stay on the scale until you see the LED displays your body weight, and follows by 3 blinking dash lines '----'. At this point the body analysis has been completed, you can check the analysis and description of each body composition by scrolling up the page.


## Claiming User Data

If there are multiple users, please manually match the body weight data to relevant users. Otherwise, you will be prompted with "You have xxx pieces of data unclaimed!"


## Settings

- Conversion of weight unit (kg/lb)
- Check the battery level
- Share data to other users


## Smart Home Ecology System

Adding Multiple Devices in interactive smart mode.

## Product Overview

## LED Display


|                                                                                   |           |
|-----------------------------------------------------------------------------------|-----------|
|  | Bluetooth |
| <b>lb</b>                                                                         | Pound     |
| <b>kg</b>                                                                         | Kilogram  |

| LED Display                                                                      | Description                                                | Solution                                                                                                   |
|----------------------------------------------------------------------------------|------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
|  | Overload.<br>The device will power off.                    | Stop using this scale for measurement.                                                                     |
|  | Low Battery.<br>The device will power off in four seconds. | Replace all three batteries at the same time.<br>Please purchase the authorized batteries for replacement. |

## Maintenance

When carrying out usual maintenance, please ensure the practice of the following Dos and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use a propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.

## Specification

|                              |                                                                                                                                                   |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Model No.                    | : EW2S                                                                                                                                            |
| Power supply                 | : 4.5V (AAA batteries X 3)                                                                                                                        |
| Unit display                 | : kg / lb                                                                                                                                         |
| Weighing range               | : 5 - 180 kg (max. weighing 180kg)<br>11 - 396 lbs (396lbs max. )                                                                                 |
| Accuracy standard            | : 2 - 5kg ( $\pm 0.2$ kg)<br>5 - 50kg ( $\pm 0.2$ kg)<br>50 - 100kg ( $\pm 0.3$ kg)<br>100 - 150kg ( $\pm 0.4$ kg)<br>150 - 180kg ( $\pm 0.6$ kg) |
| Bluetooth version            | : 4.2                                                                                                                                             |
| Bluetooth effective distance | : 10m                                                                                                                                             |
| Size                         | : 280mm X 280mm X 25mm                                                                                                                            |
| Weigh                        | : 1.25kg (not including batteries)                                                                                                                |

### Contains

1. Lite Tracker IoT Body Scale
2. AAA batteries (3 pcs)
3. Instruction Manual And Warranty Card

## Trouble Shooting

| Problem                                                                          | Possible cause and solution                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The scale does not work.<br>The data on the scale display disappears in a flash. | <ul style="list-style-type: none"><li>• Please check that if the batteries are installed.</li><li>• Make sure batteries are placed in the correct polarities.</li><li>• The power of the batteries are too low and need replacing.</li></ul>                                                                                                                                                                              |
| The scale cannot connect to the App.                                             | <ul style="list-style-type: none"><li>• Make sure to enable the Bluetooth function on your smartphone.</li><li>• The scale can only be paired within the 'MOMAX Smart' App.</li><li>• Please reset the Bluetooth and restart the 'MOMAX Smart' App on your smartphone.</li><li>• There should not be more than one Bluetooth device trying to connect to the scale. Only one at a time for the pairing to work.</li></ul> |
| Body fat and most data are missing in the App after a measurement.               | <ul style="list-style-type: none"><li>• You are wearing shoes or socks, you need to be barefooted.</li><li>• You are overloaded the maximum capacity of 180 kg / 396 lb.</li><li>• You left the scale too soon, you need to stand on it for at least 20 secs until you see the heart rate measurement is done.</li></ul>                                                                                                  |

## Trouble Shooting

| Problem                                   | Possible cause and solution                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The weight results varied drastically.    | <ul style="list-style-type: none"> <li>The scale is placed in a compromised position (leaned against a wall/turned over/ placed on an uneven or soft surface/being stack under other items, etc.)</li> <li>0KG Reset<br/>After battery installation/replacement and before measuring body weight, place the scale on the hard floor and the scale will show its weight. Wait for about 10 secs until it is OFF, step on the scale with the bared foot for 3 secs. The scale display will start up from 0.0 then show your body weight.</li> </ul> |
| How to change the profile picture?        | <ul style="list-style-type: none"> <li>Open 'MOMAX Smart' app → Tap </li> </ul>                                                                                                                                                                                                                                                                                                                                                                                  |
| How to pair with a new smartphone device? | <ul style="list-style-type: none"> <li>Remove EW2S from 'MOMAX Smart' app from the previous device. OR</li> <li>Remove the battery from battery compartment and insert again after 10mins. Follow the instruction of the app to pair with a new device.</li> </ul>                                                                                                                                                                                                                                                                                |

## Appendix - Body Composition Data

### 1. Health Tips - About Body Fat

Human body is made up of a percentage of fat. Body fat is vital for a healthy, functioning body, protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat will damage your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror.

This is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measure of fitness than weight alone—the composition of your weight loss could mean you are losing muscle mass rather than fat—you could still have a high percentage of fat even when a scale indicates 'normal weight'.

The body fat percentage (%) : 5%-60%/0.1%

| Table of body fat rating |           |           |           |           |           |
|--------------------------|-----------|-----------|-----------|-----------|-----------|
| Standard for Men         | Age       |           |           |           |           |
| Rating                   | 20-29     | 30-39     | 40-49     | 50-59     | 60+       |
| Low                      | <13.0     | <14.0     | <16.0     | <17.0     | <18.0     |
| Normal                   | 13.0-20.0 | 14.0-21.0 | 16.0-23.0 | 17.0-24.0 | 18.0-25.0 |
| Moderately High          | 20.1-23.0 | 21.1-24.0 | 23.1-26.0 | 24.1-27.0 | 25.1-28.0 |
| High                     | >23.0     | >24.0     | >26.0     | >27.0     | >28.0     |
| Standard for Women       | Age       |           |           |           |           |
| Rating                   | 20-29     | 30-39     | 40-49     | 50-59     | 60+       |
| Low                      | <19.0     | <20.0     | <21.0     | <22.0     | <23.0     |
| Normal                   | 19.0-28.0 | 20.0-29.0 | 21.0-30.0 | 22.0-31.0 | 23.0-32.0 |
| Moderately High          | 28.1-31.0 | 29.1-32.0 | 30.1-33.0 | 31.1-34.0 | 32.1-35.0 |
| High                     | >31.0     | >32.0     | >33.0     | >34.0     | >35.0     |

## 2. Health Tips - About Body Water

Body water is the single most important component of body weight. It represents over half of your total weight and almost two-third of your lean body mass (predominantly muscle). Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body's temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be 'normal' but your body hydration level could be insufficient for healthy living.

The body water percentage (%) : 43% -73% /0.1%

Table of body water level (Unit: %)

| Rating | BF % RANGE   | OPTIMAL TBW % RANGE |
|--------|--------------|---------------------|
| Men    | 4% to 14%    | 70% to 63%          |
|        | 15% to 21%   | 63% to 57%          |
|        | 22% to 24%   | 57% to 55%          |
|        | 25% and over | 55% to 37%          |
| Women  | 4% to 20%    | 70% to 58%          |
|        | 21% to 29%   | 58% to 52%          |
|        | 30% to 32%   | 52% to 49%          |
|        | 33% and over | 49% to 37%          |

## 3. Health Tips - About Muscle Mass

According to the American College of Sports Medicine (ACSM), lean muscle mass may decrease by nearly 50 percent between the age of 20 and 90. If you do nothing with it you're losing muscle and increasing fat. It is also important to know your muscle mass % during weight reduction. At rest, the body burns approximately 110 additional calories for each kilo of muscle gained. Benefits of gaining muscle mass include:

- Reversing the decline in strength, bone density and muscle mass with aging
- Maintenance of flexible joints
- Guide weight reduction when combined with a healthy diet.

Muscle Mass Percentage

|       |                                        |
|-------|----------------------------------------|
| Men   | Approximately 40% of total body weight |
| Women | Approximately 30% of total body weight |

#### 4. Health Tips - About Bone Mass

Regular exercise and a balanced diet can help maintain healthy bones. Like muscle, bone is a living tissue that can respond to exercise by becoming stronger. For most people, bone mass peaks in their thirties. Then people begin to lose bone mass. Regular exercise can help prevent that loss.

Calcium and vitamin D, good sources of which are dairy products, green leafy vegetables and fish, contribute to healthy bones.

The bone mass readings given by this product are an estimation of the amount of bone in your body. Those with osteoporosis or low bone densities may not get accurate estimations. If you have any concern regarding your bones please consult your doctor timely.

Bone Mass Percentage

- The average bone mass percentage for both men and women is between 4 to 5%.

Table of bone mass between gender and age

| Sex    | Age   | Low [-] | Normal [0]  | High [+]    | Very High [++] |
|--------|-------|---------|-------------|-------------|----------------|
| Men    | 20-39 | <8.0    | 8.0 - 19.9  | 20.0 - 24.9 | ≥ 25.0         |
|        | 40-59 | <11.0   | 11.0 - 21.9 | 22.0 - 27.9 | ≥ 28.0         |
|        | 60-79 | <13.0   | 13.0 - 24.9 | 25.0 - 29.9 | ≥ 30.0         |
| Female | 20-39 | <21.0   | 21.0 - 32.9 | 33.0 - 38.9 | ≥ 39.0         |
|        | 40-59 | <23.0   | 23.0 - 33.9 | 34.0 - 39.9 | ≥ 40.0         |
|        | 60-79 | <24.0   | 24.0 - 35.9 | 36.0 - 41.9 | ≥ 42.0         |

Source: NIH/WHO guidelines for BMI

Source: Gallagher et al., American Journal of Clinical Nutrition, Vol. 72, Sept. 2000

#### 5. Health Tips - About BMI

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. The Analyzer will calculate your BMI on your height and weight.

#### 6. Health Tips - About Calorie

The device uses your age, height, gender and an activity level setting to calculate the estimated number of calories you can consume a day to maintain your present weight. This estimation may gradually help you with a weight loss, gain or maintenance plan.

#### EMC Guidance

1. MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance  $d = 3,3$  m away from the equipment.

(Note. As indicated in Table 6 of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields  $d = 3,3$  m at an IMMUNITY LEVEL of 3 V/m)

## Warnings

Before using, please read this user manual carefully and keep it for future reference.

1. MOMAX Smart Body Scale EW2S offers you a seamless way to manage your health.
2. Please note that this device is designed for self-measurement and self-monitoring of physical health. Any information provided by this device is in no way meant to treat, cure or prevent any disease. If in doubt, please contact your doctor.
3. This device is prohibited to be used by any female user who is suspected pregnant or being pregnant. Otherwise, the effect of the device on the fetus is unknown.
4. This device is prohibited to be used by any user who is wearing or implanted with an electronic device or instrument (such as a pacemaker or defibrillator).
5. This device should not be used by anyone who is acutely or chronically ill because of suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.
6. **WARNING:** No modifications of this equipment are allowed. This may result in increased EMISSIONS or decreased immunity of scale.
7. Please use the device according to the user manual. Any misuse can cause electric shock, burns, fire and other unexpected hazards.
8. Please use and store the device under the environment which was provided in the user manual.
9. Don't expose the device to extreme temperatures, direct sunlight, moist or corrosive environment.
10. Do not step on the scale when your body or feet are wet, especially after bathing or showering to prevent slipping.

11. Please keep the device out of reach of infants, children or pets, since inhalation or swallowing of small parts is dangerous or even fatal.
12. Please do not use this device, for example, due to the performance of the sensor and the electrode is degraded or the electrode is loose, it will reduce the performance or cause other problems.
13. If you have an allergy to stainless steel, please avoid contacting the electrodes of the device.
14. Please do maintenance according to the instructions in the user manual.

## Service Agreement of Smart Mobile Application Platform

Please reference Mobile Apps content: Before accepting the Agreement, you are supposed to thoroughly read all the contents of the Agreement, and fully understand its terms, especially restrictive clauses or exceptions and Privacy Notice. Restrictive clauses or exceptions are bold or highlighted in other forms to catch your attention. In case of any doubt about the terms of the Agreement, please contact the relevant business department of Tuya. You are not allowed to use the service before you have read and accepted all the terms of the Agreement, relevant agreements, and rules, etc. Once you select "agree and submit the Agreement" (see registration page for detailed wordings) and complete the registration procedure, or you use the service in any form, it will be deemed that you have read and agreed with the restriction of the Agreement and rules above. In case of any breach of the Agreement, Tuya has the right to unilaterally restrict, suspend or terminate the service for you, and has the right to investigate your relevant responsibilities.

## Disposal of Electronic Products

1. Do not place the product as household waste;
2. Please dispose of discarded electronic products properly according to local regulations;
3. MOMAX encourages you to participate in the electronic recycling program.

## Warranty Conditions

1. For any defect, in the judgment of our technician, caused under normal use, we are responsible for repairing free of charge 12 months commencing from the date of purchase at the discretion of our company;
2. A valid warranty card and original purchase invoice must be presented upon request of warranty service;
3. This warranty card shall not apply to the product if:
  - (a) It is consumable part(s), such as Filter, or;
  - (b) it has explicitly modified or repaired in any way by persons other than technicians of this company, or;
  - (c) it has been damaged through misuse, malicious acts, negligence, the inflow of any liquid, or;
  - (d) losses or damages are caused by force majeure or natural disasters.
4. No alteration to this warranty card will be allowed;
5. The warranty terms and conditions are subject to change without prior notice;
6. Please retain your proof of purchase for warranty service;
7. Please consult the dealer about the product compatibility with your cell phone/ PDA before purchase. Product specifications are subject to change without prior notice;
8. Please visit [www.momax.net](http://www.momax.net) for the latest product information.

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Copyright © 2020 MOMAX Technology Ltd. All rights reserved.
MOMAX is registered trademark of MOMAX Technology Ltd.
All other trademarks are the property of their respective owners.