

Charge port Power button

Audio output



Receive auscultation sounds from Stemoscope PRO (Do NOT support Stemoscope or Stemoscope II)

To Computer Microphone Jack

+ : increase volume
- : decrease volume

Steps to use the Bluetooth adapter with Stemoscope PRO

Step 1. Power on the Bluetooth adapter by pressing  and hold on until the red and blue LEDs flash alternately. The Bluetooth adapter is waiting for pairing.

Step 2. Slide the switch of the Stemoscope PRO to the left and the Stemoscope PRO will work in Earphone Mode. Stemoscope PRO will search for the Bluetooth adapter and pair automatically.

Step 3. Stemoscope PRO will transmit sounds to the Bluetooth adapter after pairing.

Mobile phone charger/computer can be used to charge the Bluetooth adapter through the microUSB port.