

FACIAL INFLAMMATORY DISEASE

The woman (age 27) first presented with an extensive inflammatory face infection after 6 months of Anti-biotic treatment.

Starting state before treatment:

- very dry skin with scales
- pigment spots after purulent manifestations
- very widespread pores
- newly inflamed pimples
- black dots
- painful bulky skin depth
- tingling sensation in the face
- at the pigment site, the skin's sensitivity is heightened

First 10 treatment days, BIOPTRON light was applied twice daily (10 minutes), and starting from the 11th treatment day, treatment was once per day for 12 min

After 14 treatment days, the skin is not as dry. There is a lack of tension in the face, visible loss of inflamed bearings. The skin is more flexible and "lighter".

After 26 treatment days, dry scales cease to form. The pores are no longer open. Less inclination to create new pimples, black dots cease to form. The affected pigment spots gradually return to normal skin colour.



BIOPTRON 
HYPERLIGHT THERAPY SYSTEM by Zipler Group

Before treatment




BIOPTRON 
HYPERLIGHT THERAPY SYSTEM by Zipler Group

After 10 treatments
11th April 2018



BIOPTRON 
HYPERLIGHT THERAPY SYSTEM by Zipler Group

After 14 treatments
27th April 2018



BIOPTRON 
HYPERLIGHT THERAPY SYSTEM by Zipler Group

After 26 treatments
11th May 2018