

USER MANUAL



Our advanced teeth whitening kit features a special LED mouthpiece that uses cold light technology to help whiten teeth. Blue light interacts with the professional whitening gel, which can effectively remove stains on the tooth's surface.

The soft mouthpiece is designed based on principles of sound oral engineering. When you bite down on the mouthpiece, it will mold completely to your teeth. This will allow you to feel more comfortable over the course of the whitening process.

Before starting treatment, it is very important to know how this kit works and what needs to be prepared. Please read through these instructions carefully.

WARRANTY

10-Year-Warranty

The Perfect Smile Whitening Kit is designed to last so that we can do our part in cutting down annual environmental waste. Is your Perfect Smile Whitening Kit not working? We offer replacements for these technical issues. Please reach out to the CRI Naturals customer service department. We'd be more than happy to assist you.



100% Money Back 90 Day Guarantee

• **Refund Policy** - Here at CRI Naturals we believe that you should be 100% happy with your purchase or you don't pay a penny. If for whatever reason you're not happy with your Perfect Smile Whitening Kit, simply return your order within 90 days of receiving your order(s) and we'll give you a full refund. No hassles. No headaches. No extra hoops to jump through. It's that easy.



• **Your Information is Secure** - All online orders are 100% secure. We do not share your personal information without your explicit consent.

• **Customer Satisfaction is Our #1 Priority** - The Perfect Smile Whitening Kit can help you clean your teeth while also reducing your impact on the environment. Don't miss out!

USER MANUAL



START TO WHITEN



Rinse mouth with water

Apply the gel on tooth surface



Put mouthpiece into your mouth and bite it, then press power button for 1.5 seconds to start the treatment

16 minute treatment



Rinse mouth with water

Rinse mouthpiece

HOW TO USE

1. Rinse your mouth with water and dry your teeth.
2. Rinse the mouthpiece with water.
3. Look in the mirror to record your whitening stage with the shade guide paper.
4. Twist whitening pen until the gel flows. Apply gel evenly on both upper and lower teeth. Use about 0.5mL-0.7mL each time.
5. Place the mouthpiece in your mouth and bite it, then press the power button for 1.5 seconds to start the treatment. (Tips: start with blue color for five seconds, then press the power button to turn red; press the power button twice after five seconds to turn red and blue.)
6. The LED light will turn off automatically after 16 minutes.
7. Rinse the mouthpiece and your mouth with water.
8. Compare your teeth with the shade guide paper again to see how many stages you have improved.
9. Don't eat or drink anything except water for 30 minutes after treatment.

INGREDIENTS

Carbamide Peroxide, Glycerol, Propylene Glycol, Deionized Water, Carbomer, Carboxymethyl, Triethanolamine, Poly Vinylpyrrolidone, Menthol.

RANGE OF APPLICATION

Teeth stains from coffee, tobacco, genetic, tetracycline, food coloring, aging and other factors that may cause discoloration.

STORAGE

Avoid sunlight. Store in a cool, dry place.

SHELF LIFE

12 Months

NOTE

1. Charging: Takes two hours to fully charge.
2. Charging Status: Light flash
3. Charged Full: Light off
4. Press the power button 1.5 seconds to turn on/off.

REMINDERS

1. Dry teeth and mouthpiece thoroughly; water will dilute the whitening gel.
2. Avoid applying the gel on gums or lips. Rinse immediately with water if this happens. Any burning sensation will cease within 24 hours.
3. It's normal to feel slight sensitivity during treatment.
4. Keep out of reach of children.
5. Seek doctor's advice prior to use if pregnant, breast-feeding, experiencing problems with gums or having recently undergone oral surgery.

CAUTION

1. Not suitable for ceramic or false teeth.
2. Not suitable for tooth discoloration caused by lesions or medication.
3. Not suitable for severe tetracycline and decayed teeth.
4. Not suitable for defective enamel, exerted dentin and damaged teeth.
5. Not suitable for children under 16 years old.