ENGLISH

Invisalign® cheek retractor: Instructions for use

Hold the handles, squeeze and bend inward.

Ease the curved edges into your mouth and slowly release.

Close your jaw once the cheek retractor is in place.

- To remove, hold the handles again, squeeze and bend inward, and then pull out.

Intended use

Invisalign cheek retractor is intended for pulling the cheeks back from the surfaces of the teeth during activities and procedures related to dentistry.

Cleaning (for household or clinical setting)

- This device should be cleaned with soap and warm water before first use.
- After use rinse the device with soap and warm water.
- Visually inspect the device for cleanliness and repeat the process above, if necessary.
- Allow the device to dry thoroughly.
- Store in a clean and dry location.
- Inspect the device for defects such as cracks and breaks.
- Dispose device if defects are observed. And contact your doctor on how to obtain a replacement.

Warnings and precautions

Keep the device protected from dust and contaminants. Device intended for single patient use, do not share to prevent contamination.

Symbols

Symbols	Title	Explanatory text	Standard reference
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	EN ISO 15223-1*; 5.4.3
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC.	EN ISO 15223-1*; 5.1.1
	Authorized representative in the European Community	Indicates the authorized representative in the European Community.	EN ISO 15223-1*; 5.1.2
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	EN ISO 15223-1*; 5.1.6
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	EN ISO 15223-1*; 5.1.5
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	EN ISO 15223-1*; 5.4.4
	Do not re-use	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.	EN ISO 15223-1* 5.4.2

*Medical devices-Symbols to be used with medical device labels, labelling and information to be supplied-Part 1: General requirements.

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© 2021 Align Technology, Inc. Align, Invisalign, SmartTrack, among others, are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries. 204783 Rev B
Align Technology, Inc.
Align Technology B.V.
Invisalign Australia Pty Limited