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English

MICRO LINE, TWIN LINE, GENTLE LINE

DIRECTIONS FOR FLOSSING Guide floss between your teeth. Gently slide floss into the space between the gums and the tooth, rubbing the floss gently up and down, keeping it pressed against the tooth.

DIRECTIONS FOR FLIPPICK™ Unfold pick and guide into space between teeth. Gently use an in and out motion to clean teeth and stimulate gums.

ORTHOPICK

DIRECTIONS FOR FLOSSING Slide thin arm between teeth and brace wire and guide floss between your teeth. Gently slide floss into the space between the gum and the tooth, rubbing the floss gently up and down, keeping it pressed against the tooth.

DIRECTIONS FOR FLIPPICK™ Unfold pick and guide into space between teeth. Gently use in and out motion to clean teeth and stimulate gums.

KIDS

DIRECTIONS FOR FLOSSING Hold flosser between fingers. Gently guide floss between your teeth. Floss in gentle back and forth, up and down motion alongside each tooth. Discard after use.

Caution: Adult supervision recommended under age 10.

Plackers™

Manufactured for Ranir, LLC, A Perrigo Company
4701 East Paris Ave SE, Grand Rapids, MI 49512 USA
Made in China

@plackersdental
Plackers.com

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Fabriqué pour Ranir, LLC, A Perrigo Company