

plackers

happy mouth. happy life.™

micro line™



*Plackers Floss Picks with PE floss as compared to the Oral-B® Glide Floss Picks +Scope® Outlast® with PTFE floss. Oral-B® is a registered trademark of The Gillette Company LLC. Scope® and Outlast® are registered trademarks of The Procter & Gamble Company

Use Instructions

DIRECTIONS FOR FLOSSING

While holding the handle, gently guide the floss between your teeth using a zig-zag motion. Move floss away from the gum line in a gentle back and forth, up and down motion alongside each tooth to remove food particles and plaque. Rinse floss as needed and repeat for each tooth.

DIRECTIONS FOR FLIPPICK™

Unfold pick and guide into space between teeth. Gently use an in and out motion to clean teeth and stimulate gums.

Caution: Adult supervision recommended under age 10.

#HappyMouthHappyLife
@plackersdental

click here for all
the hot flossip!

