

# plackers

happy mouth. happy life.™

twin line®



## Use Instructions

### DIRECTIONS FOR FLOSSING

While holding the handle, gently guide the floss between your teeth using a zig-zag motion. Move floss away from the gum line in a gentle back and forth, up and down motion alongside each tooth to remove food particles and plaque. Rinse floss as needed and repeat for each tooth.

### DIRECTIONS FOR FLIPPICK™

Unfold pick and guide into space between teeth. Gently use an in and out motion to clean teeth and stimulate gums.

**Caution: Adult supervision recommended under age 10.**

#HappyMouthHappyLife  
@plackersdental

click here for all  
the hot flossip!

