**FIT**

**MANUAL**

**NED • ENG • DEU • SPA**

**FRA • POL • ITA**

**- KYNETI**

# KYNETT

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**ENG** • Mount and use the Kynett FIT as described in the manual.

The Kynett FIT can be attached to a pole or tree but can also be screwed to a wall. When screwing the Kynett FIT to a wall, make sure it is a solid (concrete) wall.

Always check before use whether all parts are fixed properly and replace worn out or broken parts as quickly as possible.

Always keep the Kynett FIT doors closed before use. Use maximum one resistance disc per side.

**WARNING!** Always carefully read the instructions of the Kynett FIT before you start training. KYNETT B.V. is not liable for personal injury or material damage caused by improper use. However, please report any serious incident to us and to the competent authority of your Member State. Users with injuries or heart patients are expected to train under the supervision of a qualified (physio) therapist.

For maintenance, only use the tooling described in the manual. Kynett B.V. is not responsible for damage to the machine and its parts when defective tooling has been used.

**Message to the user and seller:** Every serious event must be reported to the manufacturer and the authority of the member state where the event took place.

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**ENG** · Flywheel training is a unique development in sports and physiotherapy. Training with a flywheel is often referred to as kinetic training. Instead of traditional weights, a rotating disc is used during training. In the concentric phase, the flywheel is set in motion and accelerated by pulling a webbing. When the webbing is completely unwound, so much kinetic energy has been built up that the webbing rewinds itself, just like with the well-known yo-yo. The goal is to brake the flywheel and then start it again in a concentric motion. The resistance you experience is therefore not due to gravity, but due to inertia, also known as (mass) inertia or (rotating) kinetic energy. The harder you pull, the faster the flywheel will turn, and the more inertia will be created. So, you will need a lot of power to stop the flywheel afterwards!

### **Eccentric training and overload**

Eccentric training improves strength, speed, power and performance-related skills such as jump height and running speed to a greater extent than concentric or traditional training. What do we mean with eccentric training and eccentric overload? Lifting a weight with your biceps is a great example of a concentric movement. Your muscles shorten and therefore lift the weight toward the ceiling. When you extend your muscles again, the weight returns to the starting position. Lengthening the

muscles and thus lowering the weight is called an eccentric movement. It is possible to make an eccentric movement without much force. For example, by lowering a lightweight dumbbell. However, when we talk about eccentric overload, we are not talking about a lightweight movement. A movement contains eccentric overload when you produce more force in the eccentric movement than you can produce concentrically. In practice, this often means that you lower more weight than you can lift (> 1RM). Therefore, slowly lowering a weight does not have an eccentric overload if you could also lift the weight. It should come as no surprise that slow eccentric movements do not yield the same positive results as true eccentric overloads such as flywheel training. Exercises in flywheel training are also more multi-functional because at the beginning of the eccentric phase you have to slow down considerably and remain in balance and therefore you train more muscles (especially your abdominal muscles) and improve your balance.

### **For all target groups**

Training with Kynett is suitable for many target groups such as top sports, rehabilitation, physiotherapy, personal training, fitness and sport-specific training. You can train lightly or intensively, and for speed and strength.

Hoe werkt het?

How does it work?

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Comment ça marche?

Jak to działa?

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### Effect

With flywheel training you mainly train the explosive white muscle fibers (type 2). You train optimally in both directions where there is an "eccentric overload" at the end of the movement, this is greater than with traditional strength training and ensures that you also become stronger in the concentric phase. The "eccentric overload" is in particular a good stimulus for the connective tissue of the muscle and tendon (myofascial). It becomes stronger and thus you increase the load capacity. Flywheel training is proven effective eccentric strength training for the whole body and a very suitable treatment method in physiotherapy and rehabilitation. Common injuries such as hamstring injuries, Achilles tendon injury, patellar tendinopathy, shoulder instability, rotator cuff tendinopathy can be treated well with flywheel training. The flywheel training is also very suitable for fall prevention, postoperative, the prevention of sarcopenia or people who had a CVA. Several studies have shown that flywheel training has significantly more effect on muscle strength, muscle size, jump power, running speed, and strengthening connective tissue and tendons.

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**ENG** • Place the bracket against a **rigid** pole or tree and wrap the strap around it!

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**ENG** • Open the ratchet and guide the strap through the strap slit

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**ENG** • Pull the strap completely through the strap slit. Pull until it is tight around the pole

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**ENG** • The bracket set is supplied with a mountingset enabling to fix the bracket to a solid wall.

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**ENG** - Unfastening of the lashing strap starts with unbolting the ratchet. Pull the ratchet handle toward you until it has turned 180 degrees.

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**ENG** - Pull the clip on the inside of the ratchet toward you until the tension is released from the strap. Pull the strap from the strap slit.

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**ENG** • With every repetition of each exercise, the webbing must remain tight and must be unrolled as much as possible until the end of the webbing is felt. This, of course, within a safe range of motion, so no overstretching of joints, muscles and tendons. If someone is too close to the device, the webbing cannot remain tight until it is completely unwound and it takes too long to rewind so the effect of flywheel training is lost. So, to determine the correct distance where you must stand, find a place where the exercise has been performed until the end of the movement and unwind and tighten the webbing at that position.

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**ENG** · With one hand you hold the disc and with the other hand you loosen the quick release button.

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**ENG** · The quick release button can be removed

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**ENG** · Grab the disc with both hands and slide it off the center axle. Replace the disc with another disc. A thicker disc makes the exercises more difficult and slower; a thinner disc makes the exercises easier and more explosive.

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**ENG** • Reattach the replaced disc with the quick release button and close the door.

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**ENG** · For some exercises, or in some rooms, it is necessary to train with a shorter webbing. Use the strap adjust for this.

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**ENG** · Attach the strap adjust to the carabiner · Secure the strap adjust with the clips facing the Kynett

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**ENG** · Thread the loop of the Kynett's webbing through both clips from the bottom

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**ENG** · Pull the loop back between both clips and pull it tight

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
**ENG** • Take the new webbing and push it through the **narrow** opening of the shaft. Place the key back into the loop. **Take note!** With the **narrow** side of the key downward!

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**ENG** • Tighten the webbing in the shaft and screw the housing parts together

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**ENG** · During the development and improvement of the Kynett HOME, FIT, PRO and ULTIMATE and all other training accessories, we constantly choose strong and durable materials. For the mobile version; the Kynett HOME has been chosen for lightweight materials on the one hand and the highest quality in strength and weather resistance on the other.

The Kynett HOME is equipped with high-quality materials such as the engineering plastics polycarbonate and polyamide. The crucial parts that are subject to great forces are made of special aluminum alloys. These have subsequently undergone a high-quality post-treatment. The specific parts have been tested under extreme conditions to guarantee the highest quality.

The materials of the Kynett FIT, PRO, ULTIMATE and Squatplate have been adapted for heavier loads.

To ensure our quality vision, we have ensured that the production and assembly of the Kynett devices are carried out in-house. The most modern production methods are chosen for this and, where possible, production takes place fully automated in the Netherlands.

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