


Model C5 Quick start

For experienced I.C.E.S. users, getting started is easy:

- 1 – Plug in the coils, up to 4 in any combination.
- 2 – Switch it ON
- 3 – Plug in the power with the supplied USB charger cable
- 4 – Let it RUN. The C5 will run automatically in 10 seconds
- 5 – The screen will look like this:
just let it run on these settings
- 6 – Use your C5 immediately without changing anything. Then over time you can learn how to adjust the intensity, and select a protocol from the list.
- 7 – C5 User Manual is available on the C5 product page:



<https://www.micro-pulse.com/collections/frontpage/products/ices-digiceutical-c5-model-system>

NOT FDA APPROVED
Intended for experimental purposes only

Control Panel Overview:

Powering ON display

Press RED button to start previous program immediately

Adjust Power Level display

Press <YELLOW> or <BLUE> to adjust power UP or DOWN 15 (high) to 1 (low), press <WHITE> to go to next screen, or wait 4 seconds, or Press <RED> to begin immediately

Select Protocol display

Press <YELLOW> or <BLUE> to adjust protocol UP / DOWN
See all protocols in the list to the right →→
<WHITE> or <RED> to begin pulsing immediately

System Running display

Displays Protocol and Power Level, buttons are disabled
You must switch the power OFF then ON to change power or protocol settings

LEDs:   Green flashing = pulses being sent
Red steady = ERROR or malfunction

For general use:
Experiment on pain and inflammation

- C5 -- (5, 100+, 100-, 10, 5/100+, 5/100-); +/- = unipolar
- A9 -- (classic, before 5/2016) (5, 100+, 100-); +/- = unipolar
- P2 -- (5, 100+, 100-, rest 20 min) = SomaPulse, AllewaWave...
- Omni 8 -- (5, 100+, 100-, 3.9, 7.1, 10.4, 13.7, 16.9) = A9 now

Experiment with Schumann frequencies

- Schumann 1 -- (7.83 pps)
- Schumann 2 -- (7.83, 14.3 pps)
- Schumann 3 -- (7.83, 14.3, 20.8 pps)
- Schumann 4 -- (7.83, 14.3, 20.8, 27.3 pps)
- Schumann 5 -- (7.83, 14.3, 20.8, 27.3, 33.8 pps)

Experiment with fixed constant frequencies

- 1 pps -- continuous bipolar pulses at 1 pulse per second
- 2 pps -- continuous bipolar pulses at 2 pulses per second
- 3 pps -- continuous bipolar pulses at 3 pulses per second
- 4 pps -- continuous bipolar pulses at 4 pulses per second
- 5 pps -- continuous bipolar pulses at 5 pulses per second
- 10 pps -- continuous bipolar pulses at 10 pulses per second

Experiment with low-power TMS

- scTMS 10pps 30 minutes -- Low power TMS simulator
- scTMS 10pps 60 minutes -- Low power TMS simulator

Experiment with brainwave entrainment

- alpha wave -- 10 Hz to 13 Hz (10 minute cycle time)
- beta1 wave -- **LOW:** 12.5 Hz to 16 Hz (4 minute cycle)
- beta2 wave -- **MED:** 16.5 Hz to 20 Hz (4 minute cycle)
- beta3 wave -- **HIGH:** 20 Hz to 28 Hz (4 minute cycle)
- delta wave -- 1.5 Hz to 3 Hz (20 minute cycle time)
- theta wave -- 5 Hz to 6 Hz (20 minute cycle time)
- mu wave -- 8.5 Hz to 11 Hz (10 minute cycle time)
- SMA wave -- 13 Hz to 15 Hz (10 minute cycle time)
- gamma wave -- 32 Hz to 47 Hz (10 minute cycle time)

NOTE: "Hz" means pulses per second (pps); +/- = unipolar pulses