INSTALLATION
BEDIENUNG
WARTUNG

ARTIKEL 1817

DEA
DEAKTIVIERUNGS-
SCHALTER

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INSTALLATION
OPERATION
MAINTENANCE

ARTIKEL 1817

DEA
DEACTIVATION
SWITCH

Dear users,

you have chosen the IQcare DEA deactivation switch to manually control the IQcare system. The sensor mattress can be deactivated at any time. The system is automatically deactivated again after 15 minutes at the latest.

If handled properly, you can rely on this quality product at all times. Therefore, please observe the information in these installation and operating instructions.

Your IQcare team

FOR QUESTIONS:
JUST CALL OR EMAIL!

FON 0 23 54 / 944 99 69

E-MAIL INFO@IQFY.DE

NO ELEKTROSMOG

All IQcare sensors and switches transmit only very briefly when they are activated. This means that there is no electrosmog!

AREAS OF APPLICATION

In combination with an existing nursing mattress DEKU3 and the sensor system SESY or the existing alarm system PORTA (outward walking tendency).

SCOPE OF DELIVERY

DEA
Deactivation
switch

Adhesive pad

Sticker to
identify the
individual
systems

(Original may differ
from picture!)

PREPARATION

Please decide now which situation you have and where you would like to use this switch:

1 Preparation

Continue on page 17

2 Deactivation switch (working current)

Connecting/joining to existing bed escape system SESY

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Connecting/joining to existing bed escape system SESY

Continue on page 20

4 Function check	21
5 Error handling	22
6 Manufacturer's warranty	23
7 Maintenance and care	23
8 Disinfection	23
9 Technical data	23


1 PREPARATION


Stick the supplied adhesive pad to the back of the switch (or to the DEA switch support plate - not included) ...

... or disassemble the switch (no tools necessary) ...

... and screw the support plate into the desired position:

2 DEA deactivation switch on a call system (working current)

- **Important: Prerequisite for deactivation switch is that a SESY bed escape system has already been taught-in!**


- Unplug the radio receiver NAVI from the socket and wait until the LRN-button goes out.
- Plug the radio receiver NAVI radio receiver back into the socket. The LRN-button lights up permanently red.
- Press the red LRN-button permanently for approx. 10 seconds until the light goes out and the CLR-button flashes green.
- Now the NAVI is ready to learn and can be connected to the deactivation switch.

10 sec

Teach-in with existing IQcare bed escape system SESY with three or five sensors ...

... or the alarm system PORTA

- Now press 3 x quickly three times in succession the 0 of the deactivation-switch.

3 X QUICKLY

- The CLR-button lights up briefly (2 sec.) permanently green. The renewed flashing again indicates that the the switch has been recognised.

- In addition, the NAVI confirms with a signal tone that the learning operation completed.

- To save the setting briefly press the CLR- button. The LRN-key lights up permanently red.

BRIEFLY


3 DEA deactivation switch on a call system (closed-circuit current)

- **Important: Prerequisite for deactivation switch is that a SESY bed escape system has already been taught-in!**
- Unplug the radio receiver NAVI from the socket and wait until the LRN button goes out.
- Plug the radio receiver NAVI radio receiver back into the socket.
- The LRN-button lights up permanently red.
- Press the red LRN-button permanently for approx. 10 seconds until the light goes out and the CLR-button flashes red.
- Now the NAVI is ready to learn and can be connected to the deactivation switch connected to the deactivation switch.

**Teach-in with existing IQcare bed escape system
SESY with three or five sensors ...**

... or the alarm system PORTA


- Now press 3 x quickly in succession the 0 of the deactivation switch. The CLR-button lights up briefly (2 sec.) permanently red. The renewed flashing again indicates that the the switch has been recognised.
- In addition, the NAVI confirms with a signal tone that the learning process completed.
- To save the setting briefly press the CLR-button. The LRN-key lights up permanently red.

4 FUNCTION CHECK

Now press the 0 on the deactivation switch and deactivate the IQCare system. LRN and CLR must flash red simultaneously. The system is immediately reactivated by pressing 1 on the switch. If you do not actively activate the system via the switch, the system switches on again automatically after 15 minutes again.

TIP: If you have several IQcare wireless solutions in your home, number these units using our helpful identification stickers.

5 ERROR HANDLING

- In rare cases, the installation may not work as described.
- **Then it is recommended to completely reset the radio receiver to factory settings. To do this, please proceed as follows:**
 1. The setting can be deleted within 30 seconds by unplugging the radio receiver and plugging it back into the socket.
 2. Press the LRN-button for approx. 10 seconds, until the CLR-button flashes.
 3. Press the CLR-button for approx. 3 seconds until the light goes out briefly before it flashes green again. The radio receiver is now back in the delivery state and can be reset according to the operating instructions.

Attention: Please note that you must then reprogram all the previously connected components must then be reprogrammed!

ATTENTION:

If you have tuned in several accessory components to a Navi / Navi Go radio receiver and would like to remove some of them, please call our technical support and you will receive instructions on how to proceed by telephone!

6 MANUFACTURER'S WARRANTY

We provide a warranty within the scope of the statutory provisions. In the event of a complaint, please return the product to us with a meaningful description of the fault. Please return the goods carriage paid!

7 MAINTENANCE AND CARE

- IQCare products can be cleaned carefully with a soft cloth and normal household cleaning agents.

8 DISINFECTION

- Use commercially available disinfectants. A list with recommended disinfectants can be found on the homepage under the downloads ...
- ... or you can download it directly via the QR code on the right

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IQfy GmbH