INSTALLATION
BEDIENUNG
WARTUNG

ARTIKEL 8006/8007

NOBU GO

NOTRUFUHR/-KETTE

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INSTALLATION
OPERATION
MAINTENANCE

ARTICLE 8006/8007

NOBU GO

EMERGENCY CALL
WATCH/NECKLACE

Dear users,

You have chosen the IQcare NOBU GO emergency call watch/chain so that you can be informed via our Navi radio receiver when people requiring care in care facilities, shared flats or at home are in an emergency situation and need your help.

If handled properly, you can rely on this quality product at all times. Therefore, please observe the information in these installation and operating instructions.

Your IQcare-Team

FOR QUESTIONS:
JUST CALL OR EMAIL!

FON +49 23 54 / 944 99 69

E-MAIL INFO@IQFY.DE

NO ELEKTROSMOG

All IQcare sensors and switches transmit only very briefly when they are activated. This means that there is no electrosmog!

AREAS OF APPLICATION

Indoors and outdoors, e.g. room, apartment, shared flat, fenced-off outdoor area, balcony. As a stand-alone alarm system, wherever the person to be cared for is currently present.

SCOPE OF DELIVERY (depending on order)

NOBU

**NOBU GO
Watch**

**NOBU GO
Necklace**

**NAVI
Radio receiver**

(in connection with
your call system)

**NAVI GO
Mobile Radio receiver**

(e.g. in private rooms or shared
flats without a technical
prerequisite)

START

Please decide now which situation you have and which products you want to install in addition:

1 Preparation

Continue on page 17

2 NOBU GO as an independent emergency call on a call system (working current)

Continue on page 18

3 NOBU GO as stand-alone emergency call on a call system (closed-circuit current)

Continue on page 19

4 NOBU GO to existing radio receiver NAVI

Continue on page 20

5 NOBU GO on mobile receiver NAVI GO

Continue on page 21

6 Error handling

22

7 Maintenance

23

8 Desinfection

23

9 Technical Data

23

1 PREPARATION

If there is an adapter cable in the scope of delivery that has not yet been connected, please do so.

Important note: This is a **plug-in connection**, please do not twist the plug on to the connector. You will find a white line on both sides: the radio receiver and the adapter cable - these have to fit over each other.

TIP

If you have several IQcare wireless solutions in your home, number these units using our helpful identification stickers.

2 NOBU GO as an independent emergency call on a call system (working current)

- Find a socket that is close to the connection of your alarm system so that you can plug in the free end of the adapter-cable **at the end of the teach-in process.**
- Plug the radio receiver into this socket.
- The CLR-button now flashes green, the radio receiver is in the teach-in mode.

- Press 3 x quick and briefly the NOBU GO emergency button.

- The LRN-button now lights up permanently red. NOBU GO and NAVI now form one unit. The connection is confirmed by an acoustic signal.
- Connect the adapter cable.

3 NOBU GO as an independent emergency call on a call system (closed-circuit current)

- Find a socket that is close to the connection of your alarm system so that you can plug in the free end of the adapter cable **at the end of the teach-in process.**
- Plug the radio receiver into this socket.
- The CLR-button now flashes green, the radio receiver is in the teach-in mode.

- Press once on the LRN-button, the receiver is now in standby mode. The CLR-button flashes red.

- Press 3 x quick and briefly the NOBU GO emergency button. 3X

- The LRN-button now lights up permanently red. NOBU GO and NAVI now form one unit. The connection is confirmed by an acoustic signal.

4 NOBU GO on existing radio receiver

- You already use an IQcare system such as the bed escape system SESY- accordingly the radio receiver NAVI is already available.

Then proceed as follows:

- Unplug the Navi from the power socket, wait 3 seconds (until the light on the LRN button goes out) and plug it back in.
- Press the LRN-button for approx. 10 sec. until the CLR-button flashes:
 - green for working current
 - red for closed circuit

- Press the NOBU GO button briefly 3 times.

- The receiver signals the connection with an acoustic signal.
- Finally, press the CLR-button once.

5 NOBU GO with mobile radio receiver NAVI GO

- Plug the radio receiver NAVI GO into a socket.

- The CLR-button flashes green.

- Press 3 x briefly in succession on the emergency call button.

- The CLR-button goes off, the NAVI GO lights up red and is ready for use!

- As soon as the NOBU is pressed an acoustic signal sounds and the NAVI GO flashes red. The signal must be confirmed by pressing the LRN / Off-button to confirm the signal.

ERROR HANDLING

- In rare cases, the installation may not work as described.
 - **In this case, it is recommended to completely reset the radio receiver to factory settings. To do this, please proceed as follows:**
 1. by unplugging the radio receiver and plugging it back into the socket, the setting can be deleted within 30 seconds.
 2. Press the LRN-button for approx. 10 seconds, until the CLR-button flashes.
- 3a. **NAVI**
Press the CLR-button for approx. 3 seconds until the light briefly goes out before it flashes green again. The radio receiver is now back in the delivery state and can be reset according to the operating instructions.
- 3b. **NAVI GO**
Press the CLR-button for approx. 3 seconds. The Navi Go lights up red again briefly before it is back in the delivery state. The CLR-button flashes green again!

Attention: Please note that you must then reprogram all the previously connected components must then be reprogrammed!

ATTENTION:

If you have tuned in several accessory components to a Navi / Navi Go radio receiver and would like to remove some of them, please call our technical support and you will receive instructions on how to proceed by telephone!

MAINTENANCE AND CARE

- IQcare products can be cleaned carefully with a soft cloth and normal household cleaning agents.

DESINFECTATION

- Use commercially available disinfectants. A list with recommended disinfectants can be found on the homepage under the downloads ...
- ... or you can download it directly via the QR code on the right.

TECHNICAL DATA

	RADIO RECEIVER NAVI/NAVI GO	EMERGENCY CALL WATCH/NECKLACE NOBU GO
POWER SUPPLY	230VAC/50HZ	3V BATTERY
TRANSMITTING POWER:	10 MW	10 MW
AMBIENT TEMPERATURE:	0° BIS +40°C	-25 °C TO +65 °C
RELATIVE HUMIDITY:	0 % TO 93 %	0 % TO 93 %
APPROVALS:	CE	CE
DIMENSIONS (D X H):	65 X 90 X 70 MM	48 X 63 MM
SWITCHING TRAVEL:	1.5 MM	2 MM
PROTECTION CLASS:	IP 20	IP 20
RADIO SYSTEM	868.3 MHZ ENOCEAN	868.3 MHZ ENOCEAN
TRANSMIT TELEGRAM:		EEPF6-01-01, BASED ON RPS

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IQfy GmbH