manual de instrucciones

ES

instructions manual

EN

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PT

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## IMPORTANT WARNING

Read the instructions carefully before using the product. Keep them in a safe place for future reference.

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Only specialized personnel should perform repairs to this electrical appliance. Any incorrectly performed repair can lead to considerable hazards for the user. Always contact the after-sales service authorized by B&B Trends for repair. Failure to observe these instructions will lead to the total loss of the product warranty.

In no way should this device be used as a substitute for medical treatment. The measurements obtained are informative. Provide this information to your doctor. Do not take a medical decision on your own (for example modify the medication)

Children must be supervised by an adult to prevent them using the appliance as a toy.

Only use the device for the purposes indicated in this instruction manual.

Never point towards the eyes. The infrared light can damage your optical system.

This device is designed for domestic use, not professional use.

Do not use the device near other equipment that might cause interference

Do not allow the product to come into contact with any type of liquid.

Do not expose the device to extreme electrical charges or temperatures (>50°C or <0°C).

This device should be used at an ambient temperature of between 10°C and 40°C.

The device must not be subjected to heavy falls or knocks. The shelf life of the product depends on how the device is handled (beating, cleaning, etc.)

Do not use the device if any internal or external damage is observed. Check, that the lens is undamaged before each use. In the event of any damage contact an authorized technical service.

Check that the low battery indicator does not flash when turning on the product to be used. If so, change the batteries.

In order to prevent any risk of danger, do not handle the product. Only qualified technical personnel from the brand's Official Technical support service may carry out repairs or procedures on the device.

Any incorrect use or improper handling of the electric blanket shall render the warranty null and void. Only an authorized technical support service centre may carry out repairs on this product.

**B&B TRENDS SL** accepts no liability for any damage that may be caused to people, animals or objects owing to failure to observe this warning correctly.

## GENERAL INFORMATION:

The IT-122 thermometer is designed to read the body temperature on the forehead as well as of objects.

First check that the appliance is complete and undamaged. Contact the store where you purchased the product in the event of any doubt.

With your purchase you have acquired:

- 1 IT-122 Infrared thermometer
- 1 Instruction manual.

### Appliance description:

1. Trigger to turn on and take reading
2. Cover for the battery compartment.
3. Reading mode. Surface (objects) or body.
4. Results.
5. Storage position.
6. Up key
7. Lighting key.
8. Reading mode key.
9. Down key
10. Value stored in the memory.
11. Level of charge in the batteries.
12. Units (°C Celsius).
13. Sound signals.

## INSTRUCTIONS FOR USE

Before using the product for the first time, make sure that there is no apparent damage to any of the parts. If the device is damaged, please contact our aftersales service or one of the authorized distributors.

Remove the product from its original packaging taking special care not to damage any of its parts. Do not use excessive force to remove the product from the box.

Before taking any measurements, the thermometer should be in place where the measurement is to be taken for at least 30 minutes. This prevents temperature changes in the components of the product from affecting the result.

This product performs temperature measurement in direct mode.

### **Regular body temperature values:**

The IT-122-contactless thermometer allows the taking of temperature of different types of objects such as the body temperature of a person. Body temperature varies in each individual due to various factors such as sex, age, skin type, time of day, regular level of physical activity or diet.

### **Take a reading:**

Follow the instructions that are presented below in order to take a reading:

1. Place on the surface of whatever you wish to know the temperature and take a reading by pulling the trigger of the gun (1). If it is a body temperature it is recommended to take the temperature by directly applying it to the forehead from a distance of between 5 and 8cm. Make sure that there are no elements that can influence the incorrect reading of the temperature (sweat, hair, make-up, piercings or other objects).
2. Read the value of the reading given by the thermometer. Make sure that the reading mode indicated is correct (3). If the temperature of an body has been taken Body temp, should appear and if on the contrary an object temperature has been taken Body temp. should not appear. In the case that the correct mode is not indicated, it should be changed by pressing the Mode key (8) and the reading taken again.

Repeat the reading if you do not trust the validity of the reading values. In the event that you believe that the results are not correct, take the appliance to a Technical Support Service. This product cannot be calibrated by the user.

The product will automatically turn off after a few seconds without pressing the keys.

### Replacing the batteries:

Open the cover of the compartment located on the handle of the thermometer and insert two AAA batteries. Close the cover once the battery has been introduced. Carry out this task when the battery has run out or when the low battery symbol appears on the LCD screen (11).

### Error messages:

When you take a temperature reading using this thermometer one of two error messages shown below may appear:

**Hi** “Hi” (High) indicates that the average value surpasses the maximum value of the range of readings. If the Body Temp mode is selected it indicates that the temperature is higher than 42.9°C. Otherwise, if Surface/object mode is selected the device is indicating that the temperature is higher than 100°C.

**Lu** “Lu” (Low) indicates that the average value is lower than the maximum value of the range of readings. If the Body Temp mode is selected it indicates that the temperature is lower than 32°C.

Otherwise, if Surface/object mode is selected the device is indicating that the temperature is higher than 0°C.

Occasionally the surface temperature that is being measured is not really found outside the range of temperatures but the thermometer interprets it as such and one of the above error messages is shown. Below some of the reasons for this and how to solve them are shown:

CAUSE	SOLUTION
<b>Too much distance between the thermometer and the surface.</b>	Take the reading again respecting the reading distance (5-8cm).
<b>There has been a sudden change in temperature on the thermometer.</b>	Leave the thermometer in the room where the reading is to be taken for around 30 minutes so it can acclimatise and take the temperature again.

<p><b>The temperature has been taken on a bright surface.</b></p>	<p>The reading of polished or shiny surfaces cannot be taken with the infrared thermometer. Take the temperature using other means.</p>
<p><b>When taking the body temperature sweat, hairs, piercings or any other element that can distort the thermometer reading is detected.</b></p>	<p>Remove any possible element from the area that can cause an incorrect reading of thermometer and take the reading again.</p>
<p><b>The lens has become dirty</b></p>	<p>Carefully clean the lens</p>

In case of doubt in the use, maintenance, configuration or device malfunction, contact an authorized service centre or contact us via the web [www.ufesa.es](http://www.ufesa.es)

### Memory management:

This device relies on space to store a total of 32 temperature values. The values are automatically saved as the readings are taken and in the case that the memory is full the newest reading will replace the oldest. Each time a reading is taken this appears on the lower right part (10) of the LCD screen and its storage position appears on the lower left hand corner (5).

Press on the keys (6) and (9) when the thermometer is on to access the memory. You can move through all stored readings using both keys.

### Settings:

You can access the settings menu of the device to change the temperature units, to configure an alarm, to calibrate the device or to activate/deactivate the sound. To access this menu you must press the MODE key (8) for a couple of seconds. "F1" will appear on the screen after a beep. It is the first element that can be adjusted. By pressing the Mode (8) key again you will be able to move through the settings menu and go through the rest of the configurable elements.

**F1. Setting the alarm:** once inside the settings menu, when "F1" appears on the screen the device will be ready for you to select the body temperature from which you wish the acoustic alarm sound when a reading is taken. By pressing the up (6) and down (9) keys you can modify the warning temperature by intervals of 0.1°C. The predetermined value is 38°C.

**F2. Calibrate the device:** Contact an authorized technical service in case the product shows abnormal temperatures.

**F3. Activate/deactivate sound:** Once inside the settings menu, when “F3” appears on the screen the device will be ready for you to activate or deactivate the sound from the device. By pressing the up key (6) or down key (9) the sound will be activated or deactivated.

**F4. Changing modes:** Once inside the settings menu, when “F4” appears on the screen the device will be ready for you to calibrate. By pressing the up key (6) or down key (9) the device will change between body temperature or surface/object (SurF).

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Range of reading	Body mode: 32.0°C - 42.9°C (109.2°F) Surface mode: 0.0°C - 100.0°C (212.0°F)
Precision	Body mode: ±0.2°C Surface mode: ±1°C
Screen resolution	0,1°C
Memories	32
Operating Conditions	From 10°C to 40°C Humidity less than 80%
Storage and transport conditions	From -25°C to 55°C Humidity less than 90%
Power	2 AAA batteries (LR03)
Auto-power off	30 s
Time of reading	0,5 s

This product complies with the European Directive of medical product 93/42/EEC.

# ENGLISH

## STORAGE AND DISPOSAL OF THE PRODUCT

### PRODUCT STORAGE

Store the product in a clean, cool, dry environment.

Remove the batteries if you do not plan to use the appliance for an extended period of time.

### PRODUCT CLEANING

This product can only be cleaned by hand, using a soft damp cloth and mild soap. The use of abrasive cleaning products or brushes, solvents, benzene, thinner or alcohol is completely prohibited. Never immerse the appliance or any of its parts in water or excessively moisten it in order to clean it. Do not re-use the appliance until it is completely dry.

Make sure to dry the product well after cleaning it, ready for use or storage.

Keep the lens clean and free of dust.

### PRODUCT DISPOSAL

Do not dispose of the product in the rubbish at the end of its shelf-life. Contact the organic collection centre nearest your home in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU to contribute to the protection of the environment.

Request detailed information regarding the current methods for disposing of packaging and used devices from your distributor.

We hope that you will be satisfied with this product.

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**Carefully read the instructions for use.**

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**The product may contain batteries. Remove them before disposing of the product and dispose in special containers approved for this purpose.**

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**The EC marking certifies that complies with directive 93/42/EEC.**

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---

Representante autorizado en la UE.

**Authorized Representative in the EU.**

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**Elimination of waste according to the European directive 2012/19/UE on waste electrical and electronic equipment.**

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**CE** 2460

**Shenzhen Brav Electronic Technologies Co., Ltd**

ADD: 4th Floor, 11block, Tongfuyu Industrial Park, Lezhujiao, Baoan, Shenzhen, China.

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## WARRANTY REPORT

**B&B TRENDS, SL.** guarantees compliance of this product for the purpose for which it is intended for a period of two years for the entire Spanish territory. In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by **B&B TRENDS, SL.** This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

## USE OF WARRANTY

Consumers should contact a technical service centre authorised by **B&B TRENDS, SL.** for repair of the product as any tampering thereof by persons not authorised by **B&B TRENDS, SL.**, as well as misuse or inappropriate use, will lead to the cancellation of this warranty. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims must take place within three years from delivery of the product according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting **B&B TRENDS, SL.**

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**JULEN S.A.T.**