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## Court appui :

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## Long appui :

- Arrêter

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## Long press:

- Shutdown

 Silica gel head: suitable for stubborn tooth stains ;

 Flat head brush: tooth stain suitable for plane ;

 Tip brush: suitable for tooth stains between teeth ;

 Polishing cup: increase tooth finish ;

## Cleaning:

Please avoid water flowing into this product. If you need cleaning, please remove the head then wash it.

## Power supply:

This product is powered by USB by default; In case of emergency, you can select 2 batteries of AA battery to supply power; Battery can not be charged in this product;

## Products include:

One controller ,one USB cable, four working heads, one Instruction manual 、one Qualification certificate.

## Precautions :

### Warm Tips:

This product is used with toothpaste by rotating the working head. When it is used for the first time, it may cause itching or slight bleeding of gums, because gums are not used to this efficient cleaning method. It is recommended that minimum speed should be used for first use and then switch over for a period of time.

1. Before using, please ensure the plug and socket contact well
2. Tighten the brush head when installing.
3. If you want to change the Mode when using, please keep the brush head away from teeth.
4. If you find the brush head serious wear, please replace the brush head replacement.

The working head of this product is a loss product, and the loss occurred during use is a normal phenomenon. Please replace the new working head in time.

## Prohibition :

Do not immerse this product in water. This product is not waterproof; Do not immerse this product in water. This product is not waterproof; Please do not put this product into the fire to heat or use this product in a high temperature environment; Do not wrap foreign matter around the brush head. ; Do not share the brush head with others.

## Product parameter

Product Name	Multifunctional tooth cleaner
Product Model	WN1806
Rated Voltage	5V
Rated Power	≤2W
Working Temperature	0-40°C
Executive Standard	GB4706.1-2005

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